

BeGaze 2.4 Manual

Version 2.4

February 2010



Table of Contents

Part I	Introduction	2
Part II	How to Read this Document	4
Part III	Important Notice	6
Part IV	Overview	8
	1 Features and Benefits.....	8
	2 General Product Information.....	11
	BeGaze 2.4 Product Variants	11
	Dongle Protection and License Update	12
	3 How to Operate the Program.....	15
	Basic Operation	15
	Use Cases	18
	4 Getting Started.....	21
Part V	Experiment Setup	27
	1 Create Experiment Wizard.....	27
	Overview	27
	Load Experiment from Folder	28
	Experiment Name Tab	29
	Gaze Data Tab	30
	Stimulus Images Tab	34
	Stimulus Association Tab	36
	Event Detection Tab	37
	2 Measurement Scenario.....	39

3	Signal.....	39
4	Manage Experiments.....	41
	Modify Experiment	41
	Save Experiment	42
	Open Experiment	42
	Close Experiment	43
	Experiment Backup	43
	Experiment Restore	45
	Delete Experiment	46
5	Annotations.....	46
Part VI	Experiment Analysis	50
1	Overview of Analysis Data View.....	50
2	Data View Selection.....	51
3	Data Views.....	53
	Overview	53
	Operating the data views	54
	Stimulus Selection	55
	Subjects	56
	Subjects Selection.....	56
	Subject-Trial Details.....	60
	Events	61
	Events Selection.....	61
	Event Details.....	63
	Player	66
	Player Control.....	66
	Playback Control.....	67
	Zoom Control.....	69
	Thumbnail Control.....	70
	Thumbnail Control Context Menu.....	72
	Chart Display Modes	74

4	AOI Editor	75
	Overview	75
	Toolbar	77
	Open AOI Editor and Select Stimulus	78
	Create AOIs	78
	Edit AOIs	80
	Edit AOI Properties	86
	Change AOI's Visibility	90
	Navigate through Key Frames	92
	Delete AOIs	93
	Save and Load AOIs	94
	AOI Format Description	94
5	Gaze Replay	100
	Overview	100
6	Bee Swarm	102
	Overview	102
	Main Data View	104
	Settings	105
	View Settings Dialog.....	105
	Bee Swarm Tab.....	106
	Cursor Tab	107
7	Scan Path	108
	Overview	108
	Main Data View	110
	Settings	113
	View Settings Dialog.....	113
	Scan Path Tab.....	113
	Cursor Tab.....	116
	Fixations Tab.....	117
8	Focus Map	119
	Overview	119
	Main Data View	121
	Settings	122

9 Heat Map.....	125
Overview	125
Main Data View	126
Settings	128
10 Key Performance Indicators.....	130
Overview	130
Main Data View	131
Settings	133
11 Gridded AOIs.....	137
Overview	137
Main Data View	139
Scan Path Strings	141
Settings	143
12 AOI Sequence Chart.....	146
Overview	146
Main Data Tab	147
13 Binning Chart.....	150
Overview	150
Main Data Tab	151
14 Event and Reading Statistics.....	152
Overview	152
Selection Trees	153
Template List	155
Time Interval	157
Results Grid	158
Export Statistics	158
Event Statistics - Definitions and Examples	160
Reading Statistics - Definitions and Examples	184
Reading Statistics - References	199
15 Line Graph.....	203
Overview	203

Events List	204
Graph Area	206
Diagram Cursors	208
Data Table	208
Miniview	209
Settings	210

Part VII Event Detection 215

1 Built-In Event Detector	215
2 Adjust Event Detection.....	216
3 Low Speed Event Detection.....	219
4 High Speed Event Detection.....	221


Part VIII Export and Conversions 225

1 Overview.....	225
2 Export Events.....	225
Export Events	225
Export File Format	228
Export File Format.....	228
Header.....	229
Trial Section.....	229
3 Export Raw Data.....	232
Export Raw Data	232
Export Raw File Format	236
Export Raw File Format.....	236
Header.....	236
Trial Section.....	237
4 Video Files.....	239
Video Export	239
Optimizing AVI Videos	241

Converting Videos with SMI Video Optimizer	242
Background Information	244
Part IX Workspace Reference	246
1 Menu Commands.....	246
2 The Toolbar.....	249
3 Hotkeys Overview.....	251
Part X Appendix	256
1 About Box.....	256
2 Dongle - Installation and Troubleshooting.....	257
3 Experiment Types.....	258
4 Database.....	259
5 System Requirements.....	259
6 Limits.....	261
7 Program Installation.....	263
Part XI Copyright and Trademarks	266
Part XII License Agreement and Warranty	268
Part XIII About SMI	275
Index	277

Introduction

Chapter



1 Introduction


Congratulations on your purchase of SMI BeGaze™ 2.4 behavioral and gaze analysis software for eye tracking data. SMI BeGaze™ 2.4 is designed particularly for researchers working in the fields of reading research, psychology, neurology, cognitive neuroscience, marketing research and usability testing.



Document number: 091222-P-1400-001-000-A

How to Read this Document

Chapter



2 How to Read this Document

This manual is designed to serve both as online help and as printed documentation of BeGaze 2.4.

Latest software versions covered in this document: BeGaze – Version 2.4

You can use this manual in one of these ways:

- Read through the chapters pertaining to particular functions to get background information before using the program.
- Consult the manual as a reference document to find out particular information. You can find a topic either by consulting the table of contents (at the front of the manual), or the index (at the end).

All the information in this manual can also be accessed through the program. Press F1 to get help on the menu-item or the element that has been currently selected.




If you cannot find what you are looking for try searching the index.

Last updated: <February 2010>

Important Notice

Chapter



3 Important Notice

Experiment Responsibility

Make sure the presented visual stimuli do not harm or injure your subjects.

SensoMotoric Instruments GmbH is in no way responsible for the experiments you develop, execute and analyze.

Do not offend against your subject's cultural background, age, psychological condition, or similar.

Photosensitive Epilepsy

Some people may have epileptic seizures triggered by light flashes or patterns.

This may occur while presented successive pictures or video material, even if they have never had a seizure before.

Supervise your subjects during experiments.

Stop immediately and consult a doctor if a subject has the following or similar symptoms:

- Involuntary movements
- Disorientation
- Convulsions
- Loss of awareness
- Altered vision

Overview

Chapter



IV

4 Overview

4.1 Features and Benefits

Meaningful results

The Behavioral and Gaze Analysis (SMI BeGaze™ 2.4) software simplifies monocular and binocular tracking data analysis by structuring the information on experiments and subjects, as well as displaying the results as meaningful graphs – all in one advanced application.

Simultaneous analysis

- Designed to support gaze sampling rates from 50Hz up to 1250Hz
- Processes both eye and head tracking data
- Stores all movement data, subject demographics and graphics in its internal database
- Analyzes several subjects or trials simultaneously
- Changes easily the parameters for reanalyzing previous data

Various Stimuli

SMI BeGaze™ 2.4 displays, analyses and visualizes various kind of stimuli - whether

- text and graphics
- still images
- video clips and screen recordings
- websites

SMI BeGaze™ 2.4 analysis does not limit the choice of stimulus for experiments.

Multiple Subjects

- Designed to handle multiple subjects
- Integrated filter functions allow analyzing subgroups of subjects within trials based on user assigned parameters (e.g. gender, age, etc.)

Smart Visualizations

SMI BeGaze™ 2.4 provides the full spectrum of visualizations

- Gaze plots (scan path, bee swarm, gaze replay)
- Attention maps (focus map, heat map)
- Real time statistics (key performance indicators, gridded AOIs)
- Visualization parameters can be modified "on-the-fly"
- Visualizations can be exported as video (AVI) or bitmap for documentation, presentation etc.

Exploit Optimized Workflow and Interaction

SMI BeGaze™ 2.4 is not only the tool for visualization of gaze interaction with stimuli. It is also a tool to optimize workflow when it comes to quantitative analysis.

- Drill into fixation and saccade event data from scanpath or linegraph
- Find point of regard by time interval of events
- Click on data plot to view detailed information and statistics of selected events
- Customize diagrams and statistical data tables before exporting to file,
- Define your personal visualization standards and apply them across analyses or experiments etc.

AREAS OF INTEREST (AOI) – static and dynamic

- The integrated AOI editor allows definition of zones of interest
- Various geometries can be fitted to the element of interest
- Automatic Move&Morph™ function for dynamic stimuli e.g. video clips ensures the AOI being “on target” even in position and form changing elements of interest
- AOI statistics can be visualized as AOI sequence per subject, or AOI Binning Chart for groups of subjects
- The AOIs can be displayed together with gaze plot or attention map visualization
- Geometric definition of AOIs can be saved to, and loaded from file – e.g. for recurring experiments with same stimuli

Statistical Data – Your way to quantitative Analysis

- Powerful statistics module allows configuration and export of statistical data tables of more than 100 statistical variables (e.g. first fixation duration, number of glances, pupils size, blink frequencies etc.)
- Export AOI transition matrix for single or multiple subject analysis
- Export fixation and saccade parameters to file
- Measure saccade latencies and reaction times in Linegraph diagram
- Adjust event detection parameters as needed

Intelligent integration

- SMI BeGaze™ 2.4 fully integrates with SMI Experiment Center™ 2.4 - the software to make gaze tracking experiments and visual stimuli creation a snap
- Load all experiment data into SMI BeGaze™ 2.4 by 1-click: Fail-safe, fast, convenient

- SMI BeGaze™ 2.4 offers an experiment creation wizard to load manually the experiment data, allow to assign attributes to the subjects for later grouping and filtering
- Assignment of stimulus and subject's gaze data is done manually or automatically

4.2 General Product Information

4.2.1 BeGaze 2.4 Product Variants

BeGaze 2.4 is distributed in various variants that are customized to the variety of research needs.

- The BeGaze 2.4 **Light** version is delivered with the iView X™ system together with the SMI Experiment Center™ 2 Light software. BeGaze 2.4 allows to analyze experiments with two subjects and five still image stimuli and predefined video examples.
- The BeGaze 2.4 **Professional** version offers the full range of program features to analyze and export eye tracking data for still images stimuli, without any restrictions concerning the number of subjects or stimuli.
- The **Video** extension supports video stimuli in addition to still images stimuli.
- The **Reading Package** extension adds detailed statistics for reading experiments.
- The **Observation Package** extension adds the User video and User audio playback.
- The BeGaze 2.4 **Video** version offers the same range of features as the BeGaze 2.4 **Professional** version especially for video stimuli but without still image stimuli support.


4.2.2 Dongle Protection and License Update

BeGaze 2.4 is dongle-protected and requires a license. If you want to update your BeGaze 2.4 version, please contact the [SMI sales department](#) ^[275] to obtain a new license.

Collect license information

The SMI sales department will need your current license information:

1. From the Windows™ start menu, select **Programs: SMI: Experiment Suite 360° Remote Update Utility**.
2. In the **Collect Key Status Information** tab of the Remote Update Utility, click the **Collect information** button. This will acquire the current license information which is currently stored on the dongle device.



3. You will be prompted to save a file identifying your current BeGaze 2.4 license ("Save key status"). Please save the file under your last name for easy identification.


4. Send this file to sales@smi.de.

You will receive a new license key from SMI.


Update license

After you have purchased your new license key (*.v2c file format), update your license as follows:

1. From the Windows™ start menu, select **Programs: SMI: Experiment Suite 360° Remote Update Utility**.
2. Switch to the **Apply License Update** tab.



Ensure that only the BeGaze 2.4 dongle is plugged. Remove all other dongles from the PC.

3. Locate the update file (*.v2c) by clicking on the browse button  and click **Apply Update**. This will write the updated license information to the dongle device.
4. You will be prompted that a receipt has been produced to confirm the update. Please send this receipt file to sales@smi.de.
5. Close the **Remote Update Utility** and start BeGaze 2.4. You can view


detailed licensing information in the BeGaze 2.4 [About Box](#)^[256].



Type and status of your licenses are stored on the dongle device, not on the PC on which BeGaze 2.4 is installed. With the license update procedure, the dongle is updated. That means, that you can run BeGaze 2.4 on any PC when the dongle is plugged in.

Time Limited Dongles

There are dongles that contain time limited licenses for certain features. In such cases the features with time constrains can be checked in the "About" dialog.



A message will also be displayed when a feature's license expires. After the license expires the feature is no longer available for use.

4.3 How to Operate the Program

4.3.1 Basic Operation

In BeGaze 2.4 you process the measurement data with the following steps:

1. Collect and assemble all data which belong to one experiment.
2. Select an analysis, its data sources (stimulus, subjects, time interval).
3. Modify single or multiple dimensions of the data source to adapt the analysis.
4. Role over a selection of data sources to the next analysis for a different perspective or drill down.
5. Evaluate, export and/or print diagrams or data.

Data collection and experiment structure


In a typical eye tracking *Experiment*, a number of subjects are presented with a certain stimulus. For each subject a data file is recorded which is called a *Run*. In order to synchronize the measurement data with changes in stimulus presentation, the data files contain either a *Trial Number* or a *User Message* at the onset time of the stimulus change. This synchronizing information can be used to separate each run into *Trials*, where each trial is associated with a certain stimulus image. So a typical BeGaze 2.4 experiment has the following structure:

Experiment

- Run 1 (a measurement data file)
 - Trial 1 (associated with a certain stimulus image)
 - Trial 2
 - Trial [...]
- Run 2
 - Trial 1
 - Trial [...]
- Run [...]

A BeGaze 2.4 experiment is a data collection which consists of one or several measurement data files (runs), a number of stimulus images and some additional information you have to provide.

The BeGaze 2.4 experiment is assembled with the [Create Experiment wizard](#)^[27] and is stored into a [database](#)^[25], which may consist of a number of different experiments.




Combine stimulus images with the data

The **Create Experiment** wizard has automatically combined the stimulus images with the data.

Select a diagram and its data sources

After the experiment has been created, you can select the desired [diagram](#)^[51] and choose the trials from the experiment that should be displayed.



Export and print the diagram

Finally you can [export](#)^[225] the data to a text file or print the diagram.

4.3.2 Use Cases

BeGaze 2.4 can be used in a broad range of eye tracking data analyzing contexts but there are typical use cases. To get familiar with the powerful features of the program, it will be helpful to know some standard use cases.

Advertising

This use case includes the evaluation of still images (e.g. print ads) or video material (e.g. television commercials) which are presented to the subjects using the SMI Experiment Center. With this use case, you present the same visual stimuli to a larger group of subjects.

- Prerequisites:
 - min. versions for still images: iView X 2.0.23 and Experiment Center 2.0
 - min. versions for videos: iView X 2.1.16 and Experiment Center 2.1
- Experiment design: Experiment Center is used to create and record the experiment. The experiment includes various stimuli, such as videos, still images, and text.
 - Typical image presentation: Images (BMP, JPG, PNG) up to 1280x1024 pixels
 - Typical video presentation: Videos (AVI) with 30 to 300 seconds in length and a typical video size of 320x200, 640x480, or 720x576 pixels
- Experiment recording:
 - Use a proper gaze tracking device (RED, Hi-Speed, or MRI/MEG).
 - During the experiment, the data set is stored in the experiment's results folder. The data set includes the presented stimuli as well as the IDF files (gaze tracking data and user events), the subject protocols, and the meta data (subject properties, experiment design).
- Typical evaluation: The analysis of this common use case is described step-by-step in the [Getting Started](#) ^[21] topic.

Web Testing

Another use case is to evaluate web page perception and/or user navigation during web browsing sessions. This use case features the presentation of web pages to a group of subjects using the SMI Experiment Center. To evaluate the user navigation, Experiment Center provides screen recording of all actions the subjects perform during the web browsing session.

- Prerequisites: min. version is iView X 2.1.16 and Experiment Center 2.1
- Experiment design: Experiment Center is used to create the experiment and to record the subjects' web site perception and/or navigation within the site.
 - Use **Full Website** mode to store the web page as one large picture with automatic scroll compensation
 - Record keystrokes and mouse clicks
 - Optionally, use the screen recording feature to record the user actions.
- Experiment recording:
 - Use a proper gaze tracking device (RED, Hi-Speed, or MRI/MEG).
 - During the experiment, the data set is stored in the experiment's results folder. The data set includes either as a series of still images representing full web pages, or screen shots of landing pages, and (optional) background screen recordings. In the results folder, the IDF files (gaze tracking data and user events), the subject protocols, and the meta data (subject properties, experiment design) are stored also.
- Typical evaluation: Open the experiment in BeGaze 2.4 by using the [Load Experiment from Folder](#) ^[28] command. Evaluate the experiment together with the recorded mouse clicks and key presses (which BeGaze 2.4 indicates as **User Messages**) with the [Gaze Replay](#) ^[100], [Bee Swarm](#) ^[102], [Scan Path](#) ^[108], [Focus Map](#) ^[119], [Heat Map](#) ^[125] and AOI statistics data views ([Key Performance Indicators](#) ^[130], [Gridded AOIs](#) ^[137], [AOI Sequence Chart](#) ^[146] and [Binning Chart](#) ^[150]).

Software Usability

A third use case is to monitor subjects with the objective to improve software usability. For this, a group of subjects is working with a software program while their gaze tracking data and their user actions are recorded to individual videos.

- Prerequisites: min. version: iView X 2.1.16, Experiment Center 2.1
- Experiment design: Experiment Center is used to create the experiment and to record the subjects' actions (mouse clicks and key presses). For each subject, an individual video is recorded.
 - Typical video length: 60 to 300 seconds
 - Typical video size: 1280x1024 pixels / 1680x1050 pixels
- Experiment recording:
 - Use a proper gaze tracking device (RED, Hi-Speed, or MRI/MEG).
 - During the experiment, the data set is stored in the experiment's results folder. This includes the recorded videos as well as the IDF files (gaze tracking data and user events), the subject protocols, and the meta data (subject properties, experiment design).
- Typical evaluation: Open the experiment in BeGaze 2.4 by using the [Load Experiment from Folder](#) ^[28] command. Analyze the videos together with the recorded user actions, such as mouse clicks and key presses (which BeGaze 2.4 indicates as **User Messages**) with the [Gaze Replay](#) ^[100], [Bee Swarm](#) ^[102], [Scan Path](#) ^[108], [Focus Map](#) ^[119], [Heat Map](#) ^[125], and AOI statistics data view ([Key Performance Indicators](#) ^[130], [Gridded AOIs](#) ^[137], [AOI Sequence Chart](#) ^[146] and [Binning Chart](#) ^[150]).

HED Videos

Another use case is to record individual in-the-field videos while monitoring the subjects gaze position. A single subject is monitored, for example while visiting a supermarket, doing sports, or driving a car.

- Prerequisites: min. iView X 2.1
- Experiment design: For each subject, an individual real-world video is

recorded.


- Experiment recording:
 - Use the SMI Head mounted eye tracking device for real-world eye tracking studies.
 - Typical video length: 10 to 60 minutes
 - Typical video size: 752x480 pixels
- Typical evaluation: Use the BeGaze 2.4 analysis data view ([Scan Path](#)^[108] and [Attention Map](#)^[119]) and AOI statistics data view ([Key Performance Indicators](#)^[130], [AOI Sequence Chart](#)^[146] and [Binning Chart](#)^[150]) to analyze the recorded video data.

4.4 Getting Started


The following steps describe how to analyze a typical **Advertising** experiment (see [Use Cases](#)^[18]) recorded using SMI Experiment Center. If you start BeGaze 2.4 for the first time, you may proceed as described below. Alternatively, you can open one of the provided sample experiments (see [Open Experiment](#)^[42]).

1. Create a BeGaze 2.4 experiment directly from the Experiment Center's results folder (see [Load Experiment from Folder](#)^[28]).
2. Open the **Scan Path** plug-in (see [Scan Path Overview](#)^[108]).
 - Select a stimulus (see [Stimulus Selection](#)^[55]).
 - Select subjects, either manual or based on a subject property filter (see [Subjects Selection](#)^[56]).
 - Modify the **Scan Path** settings (see [View Settings Dialog](#)^[113]). For video stimuli, you may configure the "bee swarm" mode. Therefore, change the **Display** setting to **Raw Data** with the **Trailer** switched to **Constant Length** and the length slider set to zero (left image). For still image stimuli, you may change the **Display** setting to **Fixations** with the **Trailer** switched to **From Beginning**. When displaying **Fixations**, you should open the **Fixations** tab and change the **Size of**

fixation circles (right picture).



- Use the [Player Control](#) ^[66] to play the scan path presentation. To move to a specific event, use the **Events** view (see [Events Selection](#) ^[61]).
 - Export the data – either to a picture or to a video (see [Export Overview](#) ^[225]).
3. Now open the **Focus Map** data view (see [Focus Map Overview](#) ^[119]).
 - The **Focus Map** data view inherits the settings of the previously opened **Scan Path** data view. If appropriate, change the stimulus selection and the subjects selection (see above).
 - Modify the **View Settings** (see [Focus Map Settings](#) ^[122]). Change the visible area with the **Kernel Width** slider. Change the **Trailer** setting to **From Start** to see how the AOIs have evolved over time.



– Use the [Player Control](#) ^[66] to play the attention map presentation. To move to a specific event, use the **Events** view (see [Events Selection](#) ^[61]).

– Export the data – either to a picture or to a video (see [Export Overview](#) ^[225]).

4. Open the **AOI Editor** data view (see [AOI Editor Overview](#) ^[75]). This data view allows you to define **Areas Of Interest (AOIs)**. An AOI defines an image area you are interested in. AOIs are painted on top of an object in a video or image. If the subjects gaze position hits the defined area, this is evaluated as an "AOI hit". You need to define AOIs in order to use the subsequent data views (**AOI Sequence Chart** or **Binning Chart**).

– Select a stimulus (see [Stimulus Selection](#) ^[55]).

– If you have selected a video stimulus, move forward to the position in the video where you want to start with an AOI (see [Player Control](#) ^[66]).

– Select an AOI type: rectangle, polygon, or circle and paint it on the object (see [AOI Editor Toolbar](#) ^[77]). To toggle the visibility of an AOI, press the [∇] key. For a video stimulus, use the left and right arrow keys to move within the video. Use the mouse to change the position of the AOI. Note, that AOI key frames are generated when size,

position or visibility changes, while the interpolation between key frames is done automatically (tweening). For still image stimuli, AOIs are always fixed and valid for the whole selected time period.

- Rename the AOI if necessary (see [Rename AOI](#)^[81]).
 - Add more AOIs as required.
5. Open the **Key Performance Indicators** data view (see [Key Performance Indicators Overview](#)^[130]). This data view shows relevant statistical indicators for the defined AOIs.
- Modify the **View Settings** (see [Key Performance Indicators Settings](#)^[133]) to select the desired indicators and the font size used for the display.
 - Select the desired subjects, either manual or based on a subject property filter (see [Subjects Selection](#)^[56]).
 - Select the **Save Image to File...** command from the **Export** menu to export the current visualization as a picture.
6. Open the **AOI Sequence Chart** data view (see [AOI Sequence Chart Overview](#)^[146]). This data view shows the correlation between subject and AOI hits.
- Modify the settings available in the bottom view (see [Chart Display Modes](#)^[74]). It is recommended to select **Raw data** for video stimuli and **Fixations** for still image stimuli.
 - Select the desired subjects, either manual or based on a subject property filter (see [Subjects Selection](#)^[56]).
 - Select the **Save Image to File...** command from the **Export** menu to export the current visualization as a picture.
7. Open the **Binning Chart** data view (see [Binning Chart Overview](#)^[150]). This data view shows a statistical overview of AOI hits for separated time slices (bins).
- Select a stimulus (see [Stimulus Selection](#)^[55]).
 - Select the desired subjects, either manual or based on a subject property filter (see [Subjects Selection](#)^[56]).
 - Modify the settings available in the bottom view (see [Chart Display](#)

[Modes](#)^[74]). It is recommended to select **Raw data** for video stimuli and **Fixations** for still image stimuli. Modify the **Bins integration time** to your needs.

– Select the **Save Image to File...** command from the **Export** menu to export the current visualization as a picture.

Further steps depend on your requirements. For example, you may

- use other data views (see [Overview of Analysis data views](#)^[50]),
- export data to CSV files (see [Export data to files](#)^[225]),
- print or save images of the currently opened diagram (see [Export menu commands](#)^[246]), or
- backup your experiment (see [Backup](#)^[43]).

Experiment Setup

Chapter



V

5 Experiment Setup

5.1 Create Experiment Wizard

5.1.1 Overview

With the **Create Experiment** wizard you assemble all data to be analyzed to a BeGaze 2.4 experiment (see [Basic Operation](#)^[15]). There are two ways to do so.

Load experiment from folder

You can load a results folder which has been stored with the SMI Experiment Center to BeGaze 2.4 and thus easily create your experiment (see [Load Experiment from Folder](#)^[28]).

Create experiment step-by-step

Alternatively, you can create a new experiment step-by-step.

1. Click on the  icon in the [toolbar](#)^[249] or go to the **File** menu and select **New Experiment**.

The **Create experiment** dialog opens with several tabs.

2. You can proceed through the tabs step by step using the **< Back** and **Next >** buttons. You can also immediately jump to a specific tab by clicking on the tab title.
3. Fill in the experiment data in the following tabs:

[Experiment Name](#)^[29]: Experiment name and additional experiment information can be entered here.

[Gaze Data](#)^[30]: Here you select the eye tracker data files to be analyzed, if needed the plane file is selected in this tab.

[Stimulus Images](#)^[34]: All images for one experiment need to be selected in this tab.

[Stimulus Association](#)^[36]: Based on the experiment type the selected stimuli need to be associated with the trials or planes of the experiment.

[Event Detection](#)^[216]: The parameters for the fixation/saccade detection can be changed in this tab.



Note that the [Create experiment](#) button is enabled only if the experiment contains sufficient data to perform the analysis.

5.1.2 Load Experiment from Folder

You can easily create an experiment based on the data generated with the SMI Experiment Center. The stored gaze tracking data will be processed to BeGaze 2.4. During this process the stored meta data such as subject properties and the properties of the presented stimuli will be parsed and the experiment will then be created automatically in BeGaze 2.4.

Load experiment from folder

1. Select **New Experiment from Folder** from the **File** menu.

A file selection dialog opens where you can browse to the folder containing the experiment you want to load.

2. Select the appropriate folder from the directories list.
3. The Create Experiment dialog opens and the experiment is created automatically.

A progress bar indicates the creation of the experiment. After completion the new experiment is already loaded in the interface.

Load experiment from folder with drag and drop

Another way to achieve the same as the above is to simply drag the experiment folder from any file browser and drop it in the main BeGaze window. Creating the experiment then proceeds as explained above.



To load an experiment from folder, you can alternatively use the Load from Folder command which is located in the Experiment Name tab of the [Create Experiment](#) ^[27] dialog which appears when selecting **New Experiment** from the **File** menu. With this method the experiment will not be created automatically and you will be able to adjust the settings in all tabs (as explained in the following chapters) before pushing the **Create Experiment** button.

5.1.3 Experiment Name Tab

In this tab you can enter general information for the experiment. The experiment will be saved in the [database](#) ^[259] with the chosen name and description.




The **Load from Folder** command allows you to automatically fill the data and to create the experiment (see [Load Experiment from Folder](#) ^[28]).

The screenshot shows the 'Create Experiment' dialog box with the 'Experiment Name' tab selected. The dialog has a title bar 'Create Experiment' and five tabs: 'Experiment Name', 'Gaze Data', 'Stimulus Images', 'Stimulus Association', and 'Event Detection'. The 'Experiment Name' tab is active and contains the following fields and controls:

- 'Experiment name:' field with the value 'Experiment001'.
- 'Created by:' field with the value 'vladi'.
- 'Experiment description:' field with the value '12/18/2009 11:57:51 AM'.
- A 'Load from Folder...' button on the left side.
- Navigation buttons at the bottom: '< Back', 'Next >', 'Create Experiment', and 'Cancel'.

5.1.4 Gaze Data Tab

In this tab you select which eye tracker data files should be analyzed.



Select files

BeGaze 2.4 currently supports the iView X data files (*.idf) .

- If you click on **Add Files...**, a file selection dialog opens. Select one or more files for the experiment.
- To remove a file from the list, select the file and click on **Remove Files**.



Multi-Frequency support: IDF files recorded with different sampling rates are allowed in the same experiment.

Add, delete or modify subject properties

You can define individual subject "group" parameters for the experiment.

These parameters are entered as subject properties and serve as additional information to your experiment. Useful properties may be "Age" and "Gender". The first property is already defined as the subject's **Color** and can be changed at this point or later.



Subject properties are taken automatically from results generated with the SMI Experiment Center (see also [Load Experiment from Folder](#)²⁸). You can modify the properties in BeGaze 2.4 as described below.

To add new subject properties proceed as follows:

1. Click on **Add Property**.

The **Add Subject Properties** dialog opens.

The screenshot shows a dialog box titled "Add Subject Properties". It has a dark background with light text. There are two input fields: "Property name:" and "Default value:". Below the fields are two buttons: "OK" and "Cancel".

2. Enter a property name, e.g. "Gender".
3. Optionally, you can enter a default value.
4. Click **OK** to confirm your entry.


The new property will be inserted in the gaze data table. If you didn't enter a default value for the property, you can now enter a value for a selected table entry.

5. Select an entry and enter a value in the property column. If you want to change the value, simply overwrite it.

To remove an existing subject property proceed as follows:

1. Click on **Delete Property**.

The **Delete Subject Properties** dialog opens.




2. Select a property name, e.g. "Webcam".
3. Click **Delete** to delete the property.

The corresponding property column will be removed in the gaze data table.




Properties can also be directly edited in the [gaze replay](#)^[100], [bee swarm](#)^[102], [scan path](#)^[108], [focus map](#)^[119], [heat map](#)^[125], [key performance indicators](#)^[130], [gridded aoi](#)^[137], [aoi sequence chart](#)^[146] or [binning chart](#)^[150] data view when you click on the property.



Information on file entries in the data files table

- **Status:** In order to be analyzed together, all files must be recorded under the same conditions. The file to be first in the list serves as reference. All other files must fit to the reference file. If a file in the list fits the criteria, its status is ok. If a file is rejected, the status will inform of the reason of rejection and the color of the row will be red.
- **File Name** and **Date:** In these columns the file name and date are displayed.
- **Subject** and **Description:** If the files contain subject and description information they will be listed here. In this tab, they can be edited with a single click of the mouse.



- **Trials:** The number of trials in the file are computed and shown in this column.


- **Calibration Area:** Sample rate and calibration area size are presented in this column.
- **Plane file:** If the data files used require a plane stimulus file, then a **Select Plane File** button will be shown on the tab.



The planes description file comes from the Surveyor. The [measurement scenario](#)^[39] is determined by the number of planes in the selected file.

5.1.5 Stimulus Images Tab

All required stimulus images for an experiment need to be selected in this tab.



Copy movies (only for video files)

The **Copy movies** check box is checked by default. This effects that the video files used in the experiment are copied to the [database](#)^[259].

Note that this may cause a high data volume in the database directory. If you [backup an experiment](#) ⁴³, the video data will be stored in the database even if the **Copy movies** check box has been deactivated during the experiment analysis.



Warning: When the **Copy movies** check box is deactivated, videos are taken from their original location. If video files are deleted or moved, the experiment cannot be loaded any more.



Warning: Movies over 1Gb in size (HED experiments) will be automatically split into parts smaller than 1Gb while creating the experiment. The associated trial from the data file will also be split into corresponding trials (one trial per movie part).

Select files

- a) If you click on **Add Files...**, a file selection dialog will open. Select one or more files for the experiment.
- b) To remove a file from the list, select the file and click on **Remove Files**.

Information on file entries in the image files table


- **Status:** To be analyzed together, each stimulus has to meet the following criteria:
 - The format of an image file must be of type: bmp, jpg, jpeg, png.
 - The format of a video file must be of type avi and optimized with the XMP4 encoder provided in the installer (incompatible videos can be optimized with the Video Optimizer tool provided in the package)
 - The image size must be at least as large as the calibration area of the reference data, which is the first data file in the [gaze data file list](#) ³⁰.

If the stimulus fits the criteria, the status is ok. If the stimulus fails, the status will give a clue about the reason of failure and the color of the row will be red.

- **AOI File:** Images and Videos can be associated with AOI files. The AOI files should have the .xml extension (see also [AOI Format Description](#) ⁹⁴⁾) and be located in the same folder as the images. If an AOI file has the same name as an image file, except for the extension, it will be automatically added to the experiment and listed in the **AOI Files** column next to the respective image file.

5.1.6 Stimulus Association Tab

In this tab you can associate each trial (or plane in the case of a multiple plane [Measurement Scenario](#) ³⁹⁾) with a stimulus image, that will be used as background for the single views. It is recommended to set suitable associations between stimulus images and trials at an early stage of the analysis process, as it will allow an easy handling with the experiment data later on.



It's not required to make the associations. Items that have no stimulus associated will get a default gray image as background.

In the left part of the window all stimulus images of the experiment are displayed in an image pool. In the right part all trials (or planes) are listed

in the **Association** list. If the trials are separated by [trial separator messages](#) ^[39], every trial should already be associated with the appropriate stimulus image. Otherwise, the stimulus images will be sorted and associated with the trials in alphanumerical order.

Associate a stimulus image

1. Click the image you want to associate.
2. Click the trial (or plane) you want to associate.
3. Click the **Associate to selected** button.

You can also associate stimulus images with the following actions:

- a) If a trial is selected then you can simply double-click the image you want associated with it.
- b) To clear an association, select a trial and use the **Clear Association** button.
- c) All actions that can be done on one trial, can be done on multiple trials by selecting multiple trials in the trials list.
- d) With the **Associate alphabetically** button, all associations are redone by associating images to all trials in alphabetical order.

5.1.7 Event Detection Tab

In this tab you can adjust the event detection parameters for the trials loaded within the experiment. You can also adjust these settings during analysis. For information on the event detection parameters, see [Adjust Event Detection](#) ^[216].

Low - Speed data (<200Hz):

Experiment Name | Gaze Data | Stimulus Images | Stimulus Association | **Event Detection**

Low Speed Event Detection

Fixation detection parameters

Min. duration: 80 ms
Max. dispersion: 100 px

Exclude first fixation

Geometry

Use geometry from file

Stimulus screen resolution: 1280x1024

Horizontal: 1280 px Vertical: 1024 px

Phys. stimulus dimensions: Horizontal: 376 mm Vertical: 301 mm

Distance monitor-head: 700 mm

< Back | Next > | Create Experiment | Cancel

Z

High- Speed data (≥ 200 Hz) with selectable event detection algorithms, either low speed or hi-speed algorithm:

Experiment Name | Gaze Data | Stimulus Images | Stimulus Association | **Event Detection**

Event Detection: Low Speed High Speed

Saccade detection parameters

Min. duration: Auto 24 ms
Peak velocity threshold: 75 %/s

Peak velocity

Start: 20 % of saccade length
End: 80 % of saccade length

Exclude first fixation

Geometry

Use geometry from file

Stimulus screen resolution: 1280x1024

Horizontal: 1280 px Vertical: 1024 px

Phys. stimulus dimensions: Horizontal: 378 mm Vertical: 301 mm

Distance monitor-head: 655 mm

< Back | Next > | Modify Experiment | Cancel

5.2 Measurement Scenario

There are three scenarios that BeGaze 2.4 can handle:

Non Head Tracking survey:

No head tracking system was used and the raw data is mapped directly on the selected stimulus.

Single Plane survey:

Only one plane is surveyed. All measurements are performed on one single plane. The raw data is mapped on the surveyed plane. The contents of the plane may change during the experiment. Possible use case: subjects reads a newspaper.

Multiple Plane survey:

Several planes are surveyed. Each plane has a fixed content, that does not change during the experiment. The raw data is mapped to it's associated plane. Possible use case: subject sits in a cockpit and watches the various panels.

5.3 Signal

Data Trial Separator

For a better overview each BeGaze 2.4 experiment run is separated into *Trials* (see [Basic Operation](#) ¹⁵). The separation is performed automatically by "Trial Number" or by "Trial Separator Message", according to the recorded data.

The trial number and/or trial separator message was recorded by the eye tracker together with the data. Note, that iView X allows both trial number and trial separator message recording. If trial separator messages are present, BeGaze 2.4 automatically performs the separation by trial

separator message. Otherwise, the trial number separation is used.

Separation by trial number. If you use a trial number you have to set [associations](#) ^[36] between stimulus image and trials manually.

Separation by trial separator message. If you use an trial separator message it must have a specific format:

```
<Timestamp>MSG# Message: <image name>
```

Example:

```
28437864110MSG# Message: image01.bmp
```

This allows an automatic [association](#) ^[36] between stimulus images and trials. The following image and video formats are supported: bmp, jpg, jpeg, png, avi.

The separator message can be inserted in the IDF file during recording by sending the remote command ET_REM to iViewX. The format has to be:

```
ET_REM "filename.suffix"
```

Example:

```
ET_REM "image01.bmp"
```

Auxiliary Events

You can choose if *Trigger Events* should be created by *Trigger Message*. If so, the trigger message must have a specific format:

```
<Timestamp>MSG# Message: TRG: <trigger message>
```

Example:

```
28437864110MSG# Message: TRG: left Button up
```

The trigger message can be inserted in the IDF file during recording by sending the remote command ET_REM to iViewX. The format has to be:

```
ET_REM "TRG:<trigger message>"
```

Example:

```
ET_REM "TRG: left Button up"
```

5.4 Manage Experiments

5.4.1 Modify Experiment

With the **Modify Experiment** wizard you modify the data to be analyzed in the current experiment.

1. From the **File** menu, select the **Modify Experiment** command.
A dialog opens with several tabs.
2. You can proceed through the tabs step by step using the **< Back** and **Next >** buttons. You can also immediately jump to a specific tab by clicking on the tab title.
3. Fill in the experiment data in the following tabs:

Experiment Name ²⁹: Experiment name and additional experiment information can be entered here.

Gaze Data ³⁰: Here you can select the new eye tracker data files to be analyzed, and also remove from the data base the existing data. The existing data will be removed permanently. You can also add new subject properties or modify the content of existing subject properties.



[Stimulus Images](#)^[34]: Here you can add new stimuli and also remove existing stimuli from the data base. The existing stimuli will be removed permanently.

[Stimulus Association](#)^[36]: Based on the experiment type the selected stimuli need to be (re)associated with the trials or planes of the Experiment.


[Event Detection](#)^[216]: The parameters for the fixation/saccade detection can be changed in this tab.



Note that the **Modify Experiment** button is enabled only if the experiment contains sufficient data to perform the analysis.

5.4.2 Save Experiment


To save an experiment proceed as follows:

1. Click on the  icon in the [toolbar](#)^[249] or go to the File menu and select **Save Experiment**.
2. To save the experiment to a new name, click **Save Experiment As**. Enter a new name and click **Save**.

The experiment will be saved with it's current settings, for example the opened data views, in the [database](#)^[259] directory.

5.4.3 Open Experiment

To open an experiment proceed as follows:

1. Click on the  icon in the [toolbar](#)^[249] or go to the **File** menu and select **Open Experiment**.
2. The **Open Experiment** dialog opens.
3. Select the experiment you want to open.

4. Click **Ok**.

5.4.4 Close Experiment

You can interrupt the creation and analysis of an experiment by closing it. To close an experiment proceed as follows:

1. From the **File** menu, select the **Close Experiment** command.
2. Click **Save** if you want to save the experiment with its current settings, for example the opened data views. Otherwise click **Don't Save**.
3. To continue the experiment, simply [open](#)⁴² it again.

5.4.5 Experiment Backup

You can backup a saved experiment to a file. To backup an experiment proceed as follows:

1. [Close](#)⁴³ all experiments.
2. From the **File** menu, select the **Backup Experiment to File** command.



The **Backup Experiment to File** command can be performed only if all experiments are closed.

The **Select Experiment** dialog opens.

3. Select the experiment you want to backup.

Backup Experiment(s)							
Name	Description	Created By	Created On	Last Saved On	Subjects	Stimuli	Size [MB]
Sample Exp Ads Lite	2 subjects, 2 stimul...		9/28/2009 5:31:4...	10/9/2009 3:12...	2	2	18
Sample Exp Reading Light	2 subjects, 2 pages...		Not available	Not available	2	2	41
Example_Movieclip	SMI Example with...		Not available	Not available	4	2	19
CV_TEST_23			10/29/2009 12:1...	12/7/2009 12:0...	2	8	43
Reading test	- Experiment Center		Not available	Not available	2	2	41
Partybiene	5/14/2009 11:09:2...		5/14/2009 11:09...	5/14/2009 11:0...	10	1	18
Example_Slideshow	SMI example show...		9/24/2009 5:41:1...	9/24/2009 5:41...	3	12	54
CV - Videoclips			11/3/2009 2:39:2...	11/11/2009 4:5...	41	11	242
752x480_HED	12/11/2009 1:57:0... vi		12/11/2009 1:58...	12/11/2009 1:5...	1	2	27


Please select one or more experiments to backup.

Ok Cancel

4. Enter the desired experiment file name. Browse for the folder or create a new folder where the backup will be stored.

The **Experiment Backup** dialog will be presented, showing the following information:

- path of the file
- remaining time
- progress bar




5.4.6 Experiment Restore

To restore an experiment proceed as follows:

1. From the **File** menu, select the **Restore Experiment from File** command. No experiment must be loaded for the option to be available.
2. In the file selection dialog, browse for the file corresponding to the experiment you want to restore.
3. Select the experiment you want to restore.

The **Experiment Restore** dialog will be presented, showing the following information:

- path of the file
- remaining time
- progress bar



4. At the end of the process you'll be asked if you want to open the experiment.

Alternatively you can drag a backed-up experiment from a file browser and drop it in the main BeGaze window. Restoring the experiment starts automatically.



Note that the "BeGaze2\SampleExperiments" folder from the

Installation CD contains sample experiments that can be restored and used in BeGaze 2.4.

5.4.7 Delete Experiment

To delete a [saved](#) ⁴² experiment from the database proceed as follows:

1. From the **File** menu, select the **Delete Experiment from Database** command.

The **Delete Experiment** dialog opens.

2. Select one or more experiments you want to delete.
3. Click **Delete Experiment**.




The experiment will be removed from the database. This process is irreversible.

5.5 Annotations

Annotations are user defined notes associated with a certain moment of time in a data recording. They can either be previously defined during gaze recording in Experiment Center or they can be defined offline during analysis in any of the Data Views that offers a [Player Control](#) ⁶⁶.

Annotation types


All annotations have an associated *Type* to allow various filtering scenarios. Types can be defined beforehand by selecting **Define Annotations...** from the **File** menu.



The type can range from A to Z. Additionally a color and a type definition can be associated to a particular type. Types can be added and deleted from here (except for the "default" type which is always present for annotations that don't need a specific type). Types are also automatically added here when new annotations are created.

Creating and editing annotations

When adding a new annotation or editing an existing one from the context menu of the *Annotations* line in the [Player Control](#) the following window appears:



The image shows a dialog box with a grey background. At the top, there is a label 'Type' above a dropdown menu containing the text '(default)'. Below this is a label 'Text Notes' above a large, empty text area. Underneath is a label 'Operator Name' above a single-line text input field. At the bottom right, there are two buttons: 'OK' and 'Cancel'.

Here one can define the following fields:

- **Type:** any type from A to Z.
- **Text Notes:** note content.
- **Operator Name:** name of person placing the note.

Defined annotations are shown in their separate timeline underneath the [Player Control](#) thumbnails in the color defined for their type.

-

Experiment Analysis

Chapter
















VI

6 Experiment Analysis

6.1 Overview of Analysis Data View

BeGaze 2.4 provides various data views to analyze gaze data. Here is a brief overview of the data views and what they are for:

Toolbar button	Data view description
	In the AOI Editor ^[75] , you define the AOIs (A reas O f I nterest) that should be evaluated for the stimulus.
	The Gaze Replay ^[100] displays a quick overview of all stimuli associated to a subject, with a visualisation similar to the scan path one.
	The Bee Swarm ^[102] displays a raw gaze data overlay over the stimulus image/stimulus video.
	The Scan Path ^[108] displays a gaze data (raw or eye events) overlay over the stimulus image/stimulus video.
	The Focus Map ^[119] shows gaze patterns over the stimulus image visualized as a transparent map.
	The Heat Map ^[125] shows gaze patterns over the stimulus image visualized as a colored map.
	The Key Performance Indicators ^[130] displays relevant statistical data for each defined AOI over the stimulus image
	The Gridded AOIs ^[137] displays relevant statistical data for an automatically defined grid of rectangular AOIs over the stimulus image
	The AOI Sequence Chart ^[146] displays the AOI hit order over time.

	The Binning Chart ^[150] gives a statistical overview of AOI hits per binning frame.
	The Event Statistics ^[152] computes diverse statistics based on events and AOI hits.
	The Reading Statistics ^[184] computes statistics for reading experiments based on automatic generated AOIs
	The Line Graph ^[203] displays x and y directions of gaze data plotted as graphs over time and events displayed in a timeline.

Note on monocular and binocular data: The Line Graph data view shows binocular data. All other data views (except the **AOI Editor**) show monocular data.

6.2 Data View Selection

Select data view

1. Select a data view by clicking on the respective icon of the [toolbar](#)^[249]. Alternatively, you can choose the respective entry from the [Analysis](#)^[246] menu.
The appropriate data view will open in a new tab.
2. If required, you can repeat step 1 to open another data view.

Operating the data views


Each plug-in will open in a separate tab. Note that a plug-in can be opened several times within one experiment, e.g. to examine the scan path for several subjects/trials.




The AOI Editor and Gaze Replay can be opened only once in an experiment.

1. You can switch between the data views by clicking on the tab titles. You can also use the [CTRL] + [Tab] keyboard command to switch between the tabs.

If multiple tabs of a data view are opened, it may be useful to rename them for differentiation.

2. Right click the tab title.
3. In the context menu, click  to expand it.
4. Enter a new name in the **Change name** field.




5. Press [ENTER] to confirm your entry.

6.3 Data Views

6.3.1 Overview

Each visualization consists of several data views. The views contents vary but there is a standard layout:



- **Data selection view:** On the left side of the screen, you find the views to select and restrict the data to evaluate. In the [AOI Editor](#) ^[75], the left view serve to create and edit AOIs.
- **Subject Usercam and Audio:** If user videos (recorded with a webcam in Experiment Center 2.4) are available, the video corresponding to the selected subject is shown here. This view can be minimized to ignore the user video and audio completely. When the view is visible, the recorded audio is played back as well.



Usercam and Audio playback requires the observation package license.

- **Main view:** On the upper right, the main view displays the corresponding diagram, the AOI preview or the statistics.
- **Control view:** On the lower right, a control view offers individual commands for operating the display in the main view. When the webcam view is present and its panel is not minimized the subject video is played in sync with the main stimulus and the subject audio is played instead of any sound the stimulus might have.

6.3.2 Operating the data views

You can adapt the display of the views to your needs.



Resize views

1. To resize a view, position the mouse on it's border.

The mouse cursor changes to .

2. Resize the view by dragging the mouse into the desired direction.


Hide and show views

- a) To hide a view, click on it's  button.
- b) To display the view again, click on it's  button.

Sort and modify order of columns

You can sort the lists displayed in the data selection view (see [Data Views Overview](#)^[53]).

1. To sort columns, click on one of the column titles. An arrow indicates if the order is ascending or descending. To change that, click on the column header again.
2. To modify the order of the columns, click on one of the headers and move the column with the mouse to a new position (Drag & Drop).



Type	Index	Start (ms)	Duration (ms)
User Message	1	0	0
Saccade	1	2	109
Fixation			241
Saccade			22
Fixation	2	375	757
Saccade	3	1132	16
Fixation	3	1149	440
Saccade	4	1589	23

6.3.3 Stimulus Selection

The **Stimulus** selection view allows you to change the stimulus and thus the trials associated with it.

Stimulus: jamba.avi


Change stimulus: 

The stimulus selection is available in the following data views:

- [AOI Editor](#) ^[75]
- [Bee Swarm](#) ^[102]
- [Scan Path](#) ^[108]
- [Focus Map](#) ^[119]
- [Heat Map](#) ^[125]
- [Key Performance Indicators](#) ^[130]
- [Gridded AOIs](#) ^[137]
- [AOI Sequence Chart](#) ^[146]
- [Binning Chart](#) ^[150]

Select stimulus

To select a stimulus proceed as follows:

1. Click on the select stimulus button  to open a view with all available stimuli.

The file name of the currently selected stimulus is highlighted.

2. Double click on the appropriate stimulus thumbnail or click on the select stimulus button again.

The selected stimulus will immediately be displayed in the data view's main view.



You can also use the [CTRL] + [X] keyboard command to open and close the stimulus selection and you can use the left and right arrow keys to move within the stimulus selection.



You can also use the [CTRL] + [T] keyboard command to switch between a list view and a thumbnail view in the stimulus selection.

6.3.4 Subjects

6.3.4.1 Subjects Selection

In the **Subjects** view all subjects together with their associated trials are listed. The list entries are related to the selected stimulus (see [Stimulus Selection](#)^[55]).

The subjects selection is available in the following data views:

- [Gaze Replay](#)^[100]
- [Line Graph](#)^[203]
- [Bee Swarm](#)^[102]
- [Scan Path](#)^[108]
- [Focus Map](#)^[119]
- [Heat Map](#)^[125]


- [Key Performance Indicators](#) ^[130]
- [Gridded AOIs](#) ^[137]
- [AOI Sequence Chart](#) ^[146]
- [Binning Chart](#) ^[150]
- [Event Statistics](#) ^[152]
- [Reading Statistics](#) ^[152]

Select subjects

You can decide whether you want to use all subjects trials gaze data for your analysis or if you want to restrict the analysis to a subset of them by using filters. Filters are based on the subject group properties which have been set with the SMI Experiment Center. They are stored in the experiments IDF files. If no subject properties are given, you can configure them afterwards in BeGaze 2.4 by modifying the experiment (see [Modify Experiment](#) ^[41]) or by double-clicking on the property you would like to change.

You can select one or more subjects/trials with the following procedures:

- a) Click the **Select all** check box to check/uncheck all items presented in the list at once.



Subject - Trial	Color	Gender	Age
<input checked="" type="checkbox"/> cv - Trial002	Green	m	30-40
<input checked="" type="checkbox"/> pp - Trial002	Light Blue	f	20-30
<input checked="" type="checkbox"/> qq - Trial002	Yellow	m	20-30
<input checked="" type="checkbox"/> tt - Trial002	Yellow	f	30-40
<input checked="" type="checkbox"/> zz - Trial002	Orange	m	20-30

- b) To select single items, click the appropriate check box next to an item.




Subject - Trial	Color	Gender	Age
<input checked="" type="checkbox"/> cv - Trial002	Green	m	30-40
<input type="checkbox"/> pp - Trial002	Light Blue	f	20-30
<input checked="" type="checkbox"/> qq - Trial002	Yellow	m	20-30
<input type="checkbox"/> tt - Trial002	Yellow	f	30-40
<input type="checkbox"/> zz - Trial002	Orange	m	20-30

- c) Click the **Filter** check box to enable the filter setting. The subjects list displays the group properties, e.g. age. Click on \pm to open the list of given filters for this property. Select the desired filter(s). The related items will automatically be checked.



Subject - Trial	Color	Gender	Age
<input checked="" type="checkbox"/> cv - Trial002	Green	m	30-40
<input checked="" type="checkbox"/> pp - Trial002	Light Blue	f	20-30
<input checked="" type="checkbox"/> qq - Trial002	Yellow	m	20-30
<input checked="" type="checkbox"/> tt - Trial002	Yellow	f	30-40
<input checked="" type="checkbox"/> zz - Trial002	Orange	m	20-30



The checked items will represent the subjects trials used in the current analysis.

If you select an item (the selected item is highlighted), it becomes the selected trial and will be used to fill:

- and the [Trial Details](#)^[60]
- the [Events List](#)^[61]

Sorting is possible by clicking on the column titles.

Modify properties

While you are operating the [scan path](#)^[108], [attention map](#)^[119], [key performance indicators](#)^[130], [aoi sequence chart](#)^[146] or [binning chart](#)^[150] data view, you can change the properties of a subject if required. To do so:

1. Click on the corresponding property in the **Subjects** view.
2. Overwrite the property value.

Subject - Trial	Color	Gender	Age
<input checked="" type="checkbox"/> cv - Trial002	Green	m	30-40
<input checked="" type="checkbox"/> pp - Trial002	Light Blue	fe	20-30
<input type="checkbox"/> qq - Trial002	Yellow	m	20-30
<input type="checkbox"/> tt - Trial002	Yellow	f	30-40
<input type="checkbox"/> zz - Trial002	Orange	m	20-30



If you have the filter settings dialog open, you can neither select single subjects nor edit properties.




You can edit the **Color** property for several subjects at once by selecting them and clicking any color property of the selected items.

6.3.4.2 Subject-Trial Details

The **Details** view shows detailed information of the currently selected subjects trial.

The trial details view is available in the following data views:

- [Gaze Replay](#)^[100]
- [Line Graph](#)^[203]
- [Bee Swarm](#)^[102]
- [Scan Path](#)^[108]
- [Focus Map](#)^[119]
- [Heat Map](#)^[125]
- [Key Performance Indicators](#)^[130]
- [Gridded AOIs](#)^[137]
- [AOI Sequence Chart](#)^[146]
- [Binning Chart](#)^[150]



The screenshot shows a 'Details' window with a table of trial statistics. The table has two columns: 'Property' and 'Value'. The data is as follows:

Property	Value
Length	00:00:29:815
Sampling rate	50 Hz
Eye(s)	Both
Number of samples	1119
Number of fixations	35
Number of saccades	139
Number of blinks	262

If a subject trial is selected (see [Subjects Selection](#) ^[56]), information will be given about


- duration of the trial,
- sampling rate in [Hz],
- available data channels (left/right/both),
- number of samples,
- number of fixations,
- number of saccades,
- number of blinks.

6.3.5 Events

6.3.5.1 Events Selection

The **Events** views contain the summary of events of the currently selected subjects trial (see [Subjects Selection](#) ^[56]). There are two views available:

- Eye Events



The screenshot shows a table titled "Eye Events for cv - Trial003". The table has four columns: Type, Index, Start (ms), and Duration (ms). The data is as follows:

Type	Index	Start (ms)	Duration (ms)
Fixation	1	14	218
Saccade	1	233	159
Fixation	2	392	178
Saccade	2	571	39
Fixation	3	611	258
Saccade	3	870	19
Fixation	4	889	417
Blink	1	1347	79
Fixation	5	1546	278

- User Events



The screenshot shows a table titled "User Events for cv - Trial003". The table has four columns: Type, Time (ms), Event, and Content. The data is as follows:

Type	Time (ms)	Event	Content
Experiment Event	0	user event	www.smivisio...
Experiment Event	0	user event	fullWebsite...
User Action	56	scroll	scroll 0 0 55
User Action	256	scroll	scroll 0 0 55
Experiment Event	320	URL	URL complet...
Annotation	1467	A	sdsds
User Action	3948	left click	mouseclick le...

The events are listed in chronological order. For detailed information on the various eye events see [Event Details](#)⁶³. For the user events the relevant data is shown directly in the user events view:

- Type: experiment event, user action, annotation
- Event: keyboard presses, mouse clicks, page scrolls, annotation types, etc.
- Content: the relevant content for the specific event

The events views are available in the following data views:

- [Gaze Replay](#) ^[100]
- [Line Graph](#) ^[203]
- [Bee Swarm](#) ^[102]
- [Scan Path](#) ^[108]
- [Focus Map](#) ^[119]
- [Heat Map](#) ^[125]
- [Key Performance Indicators](#) ^[130]
- [Gridded AOIs](#) ^[137]

Select event

1. Mark an item by clicking on it with the left mouse button.

Now more information about the event will be given in the [Event Details](#) ^[63] field.

2. Depending on the selected data view, the main view is being updated as well. For example, when you click on a fixation in the scan path, the corresponding fixation is shown and selected also in the main view.

6.3.5.2 Event Details

In the **Details** view more detailed information of the currently selected event is displayed (see [Events Selection](#) ^[61]).

The events details view is available in the following data views:

- [Gaze Replay](#) ^[100]
- [Line Graph](#) ^[203]
- [Bee Swarm](#) ^[102]
- [Scan Path](#) ^[108]


- [Focus Map](#)^[119]
- [Heat Map](#)^[125]
- [Key Performance Indicators](#)^[130]
- [Gridded AOIs](#)^[137]

Depending on the event type, different parameters will be shown.

Fixation

If you selected a fixation, information will be given about

- start and end time,
- duration of the fixation in [ms],
- the averaged position of the fixation in [pixels],
- the dispersion of the fixation in [pixels].



The screenshot shows a 'Details' window with a table of properties and values for a fixation event. The table has two columns: 'Property' and 'Value'. The data is as follows:

Property	Value
Event type	Fixation
Start	111 ms
End	352 ms
Duration	241 ms
Position x, y	608 px, 589 px
Dispersion x, y	47 px, 23 px


If the [experiment](#)^[258] contains head tracking data in a [multiple plane scenario](#)^[39], additionally image name and plane number are displayed.

Saccade

If you selected a saccade, you will get information about

- start and end time,
- duration of the saccade in [ms],

- the amplitude of the saccade in [°],
- the average and peak velocity of the saccade in [°/sec],
- the average, peak acceleration and deceleration of the saccade in [°/sec²].



The screenshot shows a 'Details' window with a table of saccade parameters. The table has two columns: 'Property' and 'Value'. The data is as follows:

Property	Value
Event type	Saccade
Start	2 ms
End	111 ms
Duration	109 ms
Start pos x, y	630 px, 957 px
End pos x, y	635 px, 584 px
Amplitude	7,87°
Average velocity	71,85°/sec
Peak velocity	182,58°/sec at 39%
Average acceleration	4466,50°/sec ²
Peak acceleration	9129,07°/sec ²
Peak deceleration	-6522,39°/sec ²

Blinks

If you selected a blink, you will get information about


- start and end time,
- duration of the blink in [ms].

Details	
Property	Value
Event type	Blink
Start	1964 ms
End	2065 ms
Duration	101 ms

6.3.6 Player

6.3.6.1 Player Control

The player control contains commands to navigate in a video stimulus displayed in the [AOI Editor](#)^[75] and respectively in a [Gaze Replay](#)^[100], [Line Graph](#)^[203], [Bee Swarm](#)^[102], [Scan Path](#)^[108], [Focus Map](#)^[119], [Heat Map](#)^[125], [Key Performance Indicators](#)^[130] or [Gridded AOIs](#)^[137] stimulus.




Detailed descriptions for the player control elements can be found in the following sections:

- [Playback Control](#)^[67]

- [Zoom Control](#) ^[69]
- [Thumbnail Control](#) ^[70]
- [Thumbnail Control Context Menu](#) ^[72]

6.3.6.2 Playback Control





The playback control allows you to control the presentation of gaze measurement data and videos, both in playback or in single step mode.









In the **AOI Editor**, you can use the toolbar buttons to control the display of a video stimulus in the AOI main view. With the **Scan Path**, **Attention Map** or **Key Performance Indicators** data view, you use the toolbar buttons to control the display of the gaze measurement data.

Playback control buttons and key commands

To control the playback, you can use the following playback control buttons and key commands:

Button	Key command	Description
	[CTRL] + [HOME]	Jumps to the begin of the trial resp. the selected time window (see Thumbnail Control ^[70])
	Right arrow key	Moves presentation one step forward according to the selected step size (see Thumbnail Control Context Menu ^[72])
	[SPACE]	Plays/pauses the presentation
	Left arrow key	Moves presentation one step backward according to the selected step size (see Thumbnail Control Context Menu ^[72])

Button	Key command	Description
	[CTRL] + [END]	Jumps to the end of the trial resp. the selected time window (see Thumbnail Control 70 ^b)
		Repeats the presentation with the chosen playback speed under consideration of the selected start and end time (see Thumbnail Control Context Menu 72 ^b)
		For video stimuli only: activates and deactivates the speaker of the PC on which BeGaze 2.4 is running and plays the audio stream of the video Note that the speaker function only works if the video is played back with 100% playback speed (see Thumbnail Control Context Menu 72 ^b).
		Sets the playback speed.
		Sets the movement step size
		Sets the thumbnail time window size manually or automatically
	Arrow up key	increases the step size (see Thumbnail Control Context Menu 72 ^b)
	Arrow down key	decreases the step size (see Thumbnail Control Context Menu 72 ^b)
	[B]	Sets and resets a bookmark (video stimuli)
	[CTRL] + arrow right	Jumps to the next bookmark

Button	Key command	Description
	[CTRL] + arrow left	Jumps to the previous bookmark
	[ALT] + arrow right	Jumps to the next user event
	[ALT] + arrow left	Jumps to the previous user event
	[SHIFT] + arrow right	Jumps to the next annotation
	[SHIFT] + arrow left	Jumps to the previous annotation
	[CTRL] + [ENTER]	Add/Edit annotation

6.3.6.3 Zoom Control

For large images and videos, you can use the zoom control to adapt the display of the selected stimulus to the size of the data view's main view (e.g. the AOI main view of the **AOI Editor**).



Here is an overview of the buttons and what they are for:



Zooms out



Fits the stimulus display to the size of the main view



Fits the stimulus display to the width of the main view (useful for webpage stimuli)



Displays stimulus in full-scale (= original stimulus size)



Zooms in



Whether the zoom control is active or not, depends on the proportion between the BeGaze 2.4 program window size and the size of the presented stimulus.

You can also navigate in the displayed stimulus using the following procedures if you are using a mouse with a mouse wheel:

- a) Turn the mouse wheel to scroll up and down.
- b) Press the [SHIFT] key, keep it pressed and turn the mouse wheel to zoom in and out.

6.3.6.4 Thumbnail Control

The thumbnail control displays the stimulus presentation over time as a sequence of thumbnails which represent the stimulus' single images at specific timestamps. Using the thumbnail control, you can navigate in the stimulus presentation of the [Gaze Replay](#)^[100], [Line Graph](#)^[203], [Bee Swarm](#)^[102], [Scan Path](#)^[108], [Focus Map](#)^[119], [Heat Map](#)^[125], [Key Performance Indicators](#)^[130] or [Gridded AOIs](#)^[137].

The thumbnail control gives an overview on

- the time window of the trial,
- user defined bookmarks in all stimuli types (video, still image, web),
- user events (mouse clicks, page scrolls, key presses),
- and in case of a video stimulus in the AOI Editor the set [key frames](#)^[92] are shown instead of the user events.

You can adapt the settings of the thumbnail control to your needs. For example, you can restrict the number of displayed thumbnails by increasing the interval in seconds that a single thumbnail represents (see [Thumbnail Control Context Menu](#)^[72]).

Control playback using the mouse

When you grab the navigation slider with the mouse by clicking it the stimulus/video will be played back in the main view of the data view in real-time. The navigation slider moves according to the mouse movement and indicates the current position within the stimulus. You can lock the navigation slider and thus freeze the video with a single click on the appropriate thumbnail.

Add and delete bookmarks

Press B on the keyboard in order to add a bookmark on the current position where the green navigation slider is positioned. A yellow circle is added to show the bookmark positions. You can use the key combination Ctrl + Left/Right to navigate between bookmarks. Press B a second time while you are on a bookmark to delete the bookmark.

Alternatively, position the mouse over the thumbnail or the *Bookmarks* line under the thumbnails and right-click. From the context menu select "Add bookmark". If a bookmark is already present in that position right-clicking shows the option "Delete Bookmark" in the context menu.

Bookmarks are global for all data views within the experiment for the selected stimuli.

Managing annotations

Right-clicking with the mouse over the *Annotations* line under the thumbnails allows **adding** new annotations or managing existing ones from the context menu that appears. See [Annotations](#) ^[46] for more information. Adding a new annotation of a type that is not defined yet automatically adds that type to the list of defined annotation types.

Right-clicking over an existing annotation allows to **delete** it or **edit** its content. The option to filter shown annotations by their type is also available in the context menu.

Defined annotations can also be dragged left or right with the mouse in order to **change** their **position** in time.

Annotations are global for all data views within the experiment for the selected stimuli.

Filtering user events

The *User Events* line under the thumbnails shows the user events read from the recorded trial data. These are read only as they are not user defined in like the bookmarks or the user defined annotations. The context menu shown by right-clicking over the line allows to filter the user events by their type.

User events are global for all data views within the experiment for the selected stimuli.



Hovering with the mouse over a specific bookmark, user event or annotations shows a tool-tip containing relevant information (timestamp, content).

Modify Time Window

It is possible to limit the analysis time and view a smaller time window.

1. Position the mouse cursor at the left border of the first thumbnail.
2. Press the left mouse key and drag the mouse cursor on the timestamp in the thumbnail control which should define the start time.
3. Position the mouse cursor at the right border of the last thumbnail.
4. Press the left mouse key and drag the mouse cursor on the timestamp which should define the end time.
5. Position the mouse cursor on the top or bottom border of the time window.
6. Press the left mouse key and drag the mouse cursor left or right to move the whole selected time window.

Alternatively, you can use the handler to limit the time window:

1. Click on the left handler to activate it.
2. Use the left and right arrow keys to limit the time window.

The selected time window is highlighted. The movement of the navigation slider will now be restricted to this time window. Start and end time of the time window are displayed at the bottom of the thumbnails.

6.3.6.5 Thumbnail Control Context Menu

The context menu of the thumbnail control contains commands to manage the display and the replay of the stimulus.

Right click the thumbnail control. The context menu opens, offering

different commands depending on the area where the click was done:

1. Over the thumbnails and *Bookmarks* line:

- **Playback Speed:** Select one of the entries in the pop-up menu to modify the playback speed (**10%, 25%, 50%, 100%, 200%, 400%, 800%, 1600%**).
- **Thumbnail Time Window:** You can adjust the number of thumbnails which are displayed in the thumbnail control. Select one of the entries in the pop-up menu (**1 second, 2 seconds, 5 seconds, 10 seconds, Fit to Width, Fit to selection**). For example, the **5 seconds** entry will set a thumbnail every 5 seconds, the **Fit to Width** entry will distribute the stimulus' thumbnails according to the available space on the screen whereas the **Fit to Selection** will distribute the thumbnails like **Fit to Width** but only for the selected video area.
- **Step Size:** Video streams are stored as a sequence of single images. The step size determines how many image frames are skipped when you navigate the stimulus presentation with the [Playback Control](#)^[67]. Select one of the entries (**Single Step (Videoframe), Videoframe x 2, Videoframe x 4, Videoframe x 8, Videoframe x 16, Videoframe x 32, Videoframe x 64, Videoframe x 128, Videoframe x 256**).
- **Add/Remove Bookmark:** Allows adding a bookmark, or if one already exist at that timestamp, to remove it.

2. Over the *User Events* line:

- Several checkboxes to enable or disable the display of the following user event types: **Keyboard, Left Click, Right Click, Scroll, URL Loaded**.

3. Over the *Annotations* line:

- **Filter Annotations:** Check-boxes that enable or disable the display of annotations of a certain type. The annotation types are defined manually (see [Annotations](#))^[46] or automatically when defining a new annotation of an inexistent type.
- **Add Annotation:** Add a new annotation if one is not already

present at the given timestamp.


- **Delete Annotation:** Deletes the annotation if one exists at that timestamp.
- **Edit Annotation:** Edits the annotation content (type, text, operator name) if an annotation exists at that timestamp.
- **Move to Cursor Position:** Moves the annotation under the mouse to the navigation slider position.



You can also use the Arrow up and the Arrow down keys to increase/decrease the step size.

6.3.7 Chart Display Modes

In the Chart Display Modes view, you can adapt the settings for the [AOI Sequence Chart](#)^[146] and the [Binning Chart](#)^[150]. If you change a setting, the respective display will update immediately.



The view also displays a thumbnail of the currently selected stimulus to the right. Operate this view with the following steps:

1. **Base data:** Select whether AOI hits percentages are computed using data from calculated **Fixations** or measured **Raw data**.
2. **Data channel:** Select the data channel to be considered for AOI hits. In case of monocular recordings, the channel is selected automatically.
3. **Bins integration time [ms]:** Change the duration for the time slices displayed. You can adjusted the time for single time slices in milliseconds ranging from the sampling interval value up to the trial

duration. Note, that this setting is available with the **Binning Chart** data view only.



You can change, delete or create AOIs with the [AOI Editor](#)^[75].

6.4 AOI Editor

6.4.1 Overview


The following data views in BeGaze 2.4 require the existence of AOIs (**Areas Of Interest**):

- [AOI Sequence Chart](#)^[146]
- [Binning Chart](#)^[150]
- [Event Statistics](#)^[152]
- [Reading Statistics](#)^[152]
- [Key Performance Indicators](#)^[130]

AOIs can be defined for still images stimuli as well as for video stimuli where the AOIs change their position and size during the sequence of single video frames (Move&Morph™ functionality).

If you have already created AOIs for the current stimulus image, they are stored in the database and will be displayed as overlay over the image. Note, that also AOIs that were created with the iView eye tracker will be displayed if they were collected in the [Create Experiment wizard](#)^[27] with the [stimulus images](#)^[34]. If no AOIs are displayed, you have to create them prior to selecting one of the above views.

You can create new AOIs and edit or delete existing ones in the **AOI Editor**. In the following you find a short description of it's interface:



- The **AOI main view** shows all defined AOIs.
- The **AOI list view** lists all AOIs for the selected stimulus image by name. You can create new AOIs and edit existing ones via the **AOI Editor toolbar** on the right of this view. If several stimuli are used within the experiment, you can select another one via the **stimulus selection** area on the top of the AOI list view.
- In the **AOI detailed properties** view, you can view the properties of an AOI selected in the AOI list view and edit it.
- The **AOI player control** view shows the stimulus presentation over time. In case of a video stimulus, this view will show the video's contents image by image.



If the reading package is licensed, reading AOIs for paragraphs, sentences, words and character are automatically generated in Experiment Center and been used in BeGaze. These reading AOIs cannot be self created. For more information, please see [Reading AOI Statistics](#)

6.4.2 Toolbar

The **AOI Editor** toolbar is located on the right of the AOI list view. It gives you short-cuts to create and edit AOIs. Here is an overview of the buttons and what they are for:



Selects an AOI and switches to edit mode



Draws a rectangular AOI



Draws an ellipsoidal AOI



Draws a polygonal AOI



Changes the priority of overlaying AOIs. The selected AOI gets a higher priority.



Changes the priority of overlaying AOIs. The selected AOI gets a lower priority.



Deletes a selected AOI



Duplicates the selected AOI



Undoes the last step



Redoes the last step



Saves AOIs to an XML file



Loads AOIs from an XML file

6.4.3 Open AOI Editor and Select Stimulus

1. Click  in the [toolbar](#) ^[249].

The **AOI Editor** opens, displaying the experiment's stimulus. If several stimuli are used in the experiment, you can now select another one (see [Stimulus Selection](#) ^[55]).



2. Proceed with one of the following steps:
 - [Create AOIs](#) ^[78]
 - [Edit AOIs](#) ^[80]
 - [Delete AOIs](#) ^[93]

6.4.4 Create AOIs


Prerequisite


A stimulus is displayed in the AOI's main view (see also [Stimulus Selection](#) ^[55]).

Create a new AOI

1. Select the shape of the AOI you want to create by clicking on the appropriate button.
 - If you want to create an ellipsoidal AOI, click on the  button. Then left-click in the image to set the start point, keep the mouse button pressed and drag the mouse vertically over the image to define the size of the ellipse. Release the mouse button if the desired size is reached.
 - If you want to create a rectangular AOI, click on the  button. Left-click in the image to set the start point, keep the mouse button pressed and drag the mouse vertically over the image to define the size of the rectangle. Release the mouse button if the desired size is

reached.


– You can also create a polygonal AOI by clicking on the  button. Click in the image to set the starting point of the first straight line. With the second click you set the end point of the first line which is also the starting point of the second line etc. By clicking, moving the mouse, and clicking again you will define the shape of the polygon. When you have completed the AOI except for the last side of the polygon, double click the left mouse button to mark the last corner point. The last corner point of the polygon will automatically be connected with the starting point.

 In case of a video stimulus, BeGaze 2.4 will automatically set a key frame for each new AOI position, a changed AOI shape/size, and a change of the AOI visibility (see also [Navigate through Key Frames](#) ^[92]).

2. Name the AOI. A new AOI is named "AOI" followed by a serial number (e.g. AOI 001). To assign a meaningful name edit it in the box that appears immediately after you draw the AOI. You can double click the AOI afterwards to get the name editing box back.

Alternatively, you can double click the AOI in the AOI list view or click on the desired AOI in the AOI main view and overwrite the given name in the **Name** field of the AOI detailed properties view.

3. You may set another new AOI at a later time position (e.g. with a video stimulus). To do this, position the time cursor in the AOI player control on the appropriate image thumbnail (see [Thumbnail Control](#) ^[70]).
4. To create the new AOI, repeat steps 1 and 2.


 If required, you can change the position, rotation angle or the shape of an AOI. For more information, see the topic entitled [Edit AOIs](#) ^[80].

6.4.5 Edit AOIs

You can edit existing AOIs as follows:


- [rename AOI](#) ^[81],
- [change position and/or shape of a still image stimulus AOI](#) ^[83],
- [change position and/or shape of a video stimulus AOI](#) ^[85],
- [change the AOI priority](#) ^[85],
- change the visibility of a selected AOI, see [Change AOI's Visibility](#) ^[90],
- edit several properties for a selected AOI, see [Edit AOI Properties](#) ^[86].

Prerequisite

If you want to edit an AOI, you have to switch to the edit mode by clicking on the  button.


Enable/Disable AOI

- AOI's are enabled by default and can be disabled if the AOIs shall not be considered in the whole experiment (statistics, ...)
- "Enable all" allows to enable and disable all AOIs in one go or with the filter when clicking on the filter checkbox
- Individual AOIs can be enabled/disabled by clicking on the checkbox left to the AOI name.




Rename AOI

1. Double click the desired AOI in the main view and change the name.




Or you can click the AOI in the AOI list view and overwrite the given

name.



Alternatively, you can click on the desired AOI in the AOI main view and overwrite the given name in the **Name** field of the AOI detailed properties view (after expanding it).

Detailed Properties	
Visible	True
CurrentTimestamp	920000
Name	Eyes
Group	
Enabled	True
Scope	Local
Color	 Coral
Transparency	50
Angle	0
Points	AOI points.
BorderWidth	2
Style	HalfTransparent
Area	13201
Shape	Polygon


Change position and/or shape of a still image AOI

If you want to change the position or the shape of an AOI, proceed as follows:

1. Click on the desired AOI in the AOI main view.

The selected AOI is marked by selection handles (small squares at the corner points of the AOI).

Polygons and group of AOIs are marked in addition with a frame and additional handlers.



2. You can now move the AOI by clicking somewhere in the AOI area and dragging the AOI to the desired position while keeping the left mouse button pressed. To change the shape (e.g. the size) of the AOI, click on the selection handles and drag them in the appropriate directions. The AOI will behave the same as in other graphic programs.
3. AOIs can be rotated by using the round handler on top
4. You can change the size of the selected AOI by pressing the `[Shift]` key and turning the mouse wheel or by using the handlers in the corners.
5. There are two options only available when right-clicking on a polygonal AOI: **Add Point** and **Remove Point**. You can add new points to an existing polygon by hovering over an edge, right-clicking and selecting the **Add Point** option (notice the mouse cursor changing while hovering over an edge). An existing point can be removed by hovering over the point and selecting **Remove Point** from the context menu.

Change position and/or shape of a video stimulus AOI

With a video stimulus, the position and shape of one AOI can change in the course of the video. With the following steps, you adapt the AOI to the changed display detail.

1. Click on the desired AOI in the AOI main view.

The selected AOI is marked by selection handles (small squares at the corner points of the AOI).

2. In the AOI player control view, position the time cursor on the appropriate video frame (see [Thumbnail Control](#)^[70]).

The selected video frame is displayed in the AOI main view. The AOI is located on its former position.

3. Move it to its new position. If necessary, change its shape/size/rotation also (as described in the section [Change position and/or shape of a still image AOI](#)^[83]).

BeGaze 2.4 will automatically set a key frame for the new AOI position (see also [Navigate through Key Frames](#)^[92]).



Tip: It will be efficient to use key commands to navigate in the player control (see [Playback Control](#)^[67]) and to use the mouse for changes on the AOI shape and position.





Removing points from a polygon in a certain key frame affects the shape in all key frames so a warning pops up when using these options on a polygon in a video stimulus.

Change AOI Priority


If you have several AOIs in a stimulus image that overlay upon each other, and the chosen diagram only allows evaluation of one AOI per time (which is the case with the [Binning Chart](#)^[150]), only the one with the highest priority will be validated. The priority of an AOI corresponds to its position in the list view: AOIs that are placed on top of the list have a higher priority than AOIs with a lower position. You can change the priority

of an AOI by proceeding the following two steps:


1. Mark the AOI to be changed in the list view.
2. Click on the  and  buttons to move the AOI to the desired position in the list and, thus, assign it the desired priority.

6.4.6 Edit AOI Properties


You can change the properties of a selected AOI as follows:


1. Click on the  button to switch to the edit mode.
2. Click the desired AOI in the AOI list view. Alternatively, you can click on the desired AOI in the AOI main view. Expand the AOI detailed properties view.



Now you can enter the desired values directly in the AOI detailed properties view.


3. **Visible:** This field is displayed with a video stimulus only. Click on  to open the drop-down menu. Select **True** if the AOI is visible at the current timestamp and select **False** if the AOI gets invisible at this time (this means that AOI of the displayed theme fades out).
4. **Name:** If required, overwrite the given name.
5. **Group:** You can assign a group name to several AOIs and use it to sort of filter the AOI list (useful for reading or web experiments).
6. **Enabled:** This sets whether the AOI is taken into account in the other plugins (KPI, Event Statistics and so on). A disabled AOI is drawn in a dash-dot pattern instead of a full line one. This setting is identical to toggling the checkbox in front of the AOI in the AOI list. The default setting is **True**.
7. **Scope:** Can take the values of **Local** or **Global**. **Local** shows that the AOI is available for the current stimulus only and is the default setting while **Global** means it is available in the whole experiment, maintaining its name and color in all stimuli. When first creating an AOI it is set to **Local** and exists in the current stimulus only and changing it to **Global** replicates it in all the other stimuli in the

experiment. The position and shape can be changed independently in each stimulus afterwards.

- Color:** New AOIs are created with standard colors. It is recommended to change these colors if the AOIs are hardly recognizable on your stimulus image. Click on  to open the color selection drop-down field, offering separate color tabs. Select the desired color.



- Points:** Click on  to display the list of points that define the AOI's position and size. This list is dependent of the type and should contain exactly 2 points for rectangle or ellipse, and at least 3 points for polygon. You can modify the AOI's position and size by entering new values.
- Border Width:** Enter a value between 1 and 10 to define the AOI's border width. The default value is 2.
- Style:** Click on  to open the transparency selection drop-down menu. Select the transparency style.



12. Area is showing the size of the AOI in square-pixel.




The other fields in the AOI detailed properties view, such as Current Timestamp and Shape give further information on the AOI. These properties cannot be edited.

For convenience there are two alternative methods for editing the most commonly used properties rendering the Detailed Properties panel useful for advanced editing only:

1. Edit the Name, Group, Scope, Color and Enabled state (checkbox) directly in the AOI list view.




2. Edit the above and more in the context menu that shows when you right-click on an AOI in the main view. The options that are not available for the specific AOI are grayed out.




6.4.7 Change AOI's Visibility

The visibility of AOIs affects video stimuli only. A video stimulus shows the objects / protagonists / visuals you are interested in, but they may appear or disappear in the course of the video. To reflect this, an AOI can have the visible and invisible status.

1. Click on the  button to switch to the edit mode.
2. Click the desired AOI in the AOI main view.
3. Pressing the [V] key, you can toggle the visibility of the selected AOI.

Alternatively, you can set the visibility of a selected AOI in the AOI property view (see [Edit AOI Properties](#)^[86]).

Invisible AOIs are indicated with a dotted border.



Note, that no AOI hit is counted while the AOI has the invisible status. This is true even if BeGaze 2.4 detects the gaze position meets the AOI area. This means that no AOI hits are emitted in the [AOI Sequence Chart](#)^[146] and the [Binning Chart](#)^[150].


Example: In the course of the video, a new character appears on the screen. At this timestamp you draw the corresponding AOI in the video's fixed-image (the first key frame for this AOI is set). After some seconds, the character disappears. At this timestamp you set the AOI to invisible (the second key frame for this AOI is set). Some seconds later, the character appears again. You set the AOI to visible again (the third key frame for this AOI is set).

BeGaze 2.4 evaluates the AOI in the following manner: The video starts with the AOI invisible until the AOI key frame 1 is reached. Between key frame 1 and key frame 2 and from key frame 3 to the end of the video (the AOI is visible), the hits for this AOI are count. Between the key frames 2 and 3 when the AOI is set to invisible, no hits for this AOI are count even if a subject gazed at the AOI.

6.4.8 Navigate through Key Frames

Move&Morph

With a video stimulus BeGaze 2.4 sets a key frame for each AOI, and also for each changed AOI position, a changed AOI shape/size, and a change of the AOI visibility. Between the successive key frames of an AOI, BeGaze 2.4 automatically calculates the tweening of the AOI's motion and size and adapts it to the single images of the video sequence lying between these key frames. (Move&Morph)

With the help of key frames, you can navigate through a sequence of AOIs, e.g. to change their position, size or shape if necessary. The [Thumbnail Control](#) ⁷⁰ indicates the key frames which are set for a video stimulus with .






Navigate through key frames

The key frames control is located on the bottom of the **AOI Editor**.



1. Position the time cursor in the AOI player control at the beginning of the video or on the appropriate video's single image (see [Thumbnail Control](#) ⁷⁰).
2. If you want to restrict the navigation to one special AOI, now select the appropriate AOI in the AOI list view. If you want to navigate through the complete series of the stimulus' key frames, make sure that no AOI is selected.

3. Navigate through the frames:

- Click  to jump to the next key frame relative to the image currently displayed.
- Click  to move back to the previous key frame.
- Click  to delete the current key frame or press [D]



Navigate through key frames using hotkeys

You can use the following hotkeys for fast navigation through the key frames:

Keys	Description
[HOME]	jumps to first key frame
[END]	jumps to last key frame
[PG Up]	goes to next key frame
[PG Dn]	goes to previous key frame
[D]	deletes the current selected key frame

6.4.9 Delete AOIs

You can delete AOIs as follows:

1. Click on the  button to switch to the edit mode.
2. Mark one or more AOIs that should be deleted either in the stimulus image or in the AOI list view. A selection in the stimulus image will automatically select the appropriate item in the AOI list view and vice versa.
3. Click on the  button.

Alternatively, you can press the [DEL] key or right-click on the AOI and select the **Delete** option in the context menu.



When deleting AOIs that have the **Scope** setting set to **Global** a warning dialog with several options appears informing you that you are about to delete the global AOIs from all the stimuli in the current experiment.

6.4.10 Save and Load AOIs

Save AOIs

AOIs will be automatically saved in the database when you close the **AOI Editor**. You can also save AOIs in an XML file (*.xml), if, for example, you want to reuse a stimulus image with the appropriate AOIs in further experiments.

1. Click on the button and select the name and the storage folder for the XML file.

Load AOIs

1. To load AOIs for the current image click on and select an XML file (*.xml) from the file selection dialog.



To create an XML file using an external tool, follow the AOI Format Description (see [AOI Format Description](#)⁹⁴).

6.4.11 AOI Format Description

The XML file that contains the AOIs has the following structure (except for automatic generated reading AOIs):

```
<?xml version="1.0"?>
```

```

<ArrayOfDynamicAOI
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
  <DynamicAOI[96]>
    <Points[97]>
      <Point>
        <X>1003</X>
        <Y>748</Y>
      </Point>
      <Point>
        <X>1169</X>
        <Y>886</Y>
      </Point>
    </Points>
    <Enabled[96]>true</Enabled>
    <Group[96]>Main Group</Group>
    <Scope[96]>Local</Scope>
    <Angle[97]>0</Angle>
    <BorderWidth[97]>2</BorderWidth>
    <Type[96]>Rectangle</Type>
    <Style[97]>HalfTransparent</Style>
    <Transparency[97]>50</Transparency>
    <Area[97]>22908</Area>
    <Color[97]>NamedColor:Blue</Color>
    <Name[96]>Logo Name</Name>
    <Font[97]>
      <FontName>Tahoma</FontName>
      <FontSize>13</FontSize>
      <FontStyle>Regular</FontStyle>
      <FontUnit>Point</FontUnit>
      <FontGdiCharSet>1</FontGdiCharSet>
      <FontGdiVerticalFont>>false</FontGdiVerticalFont>
    </Font>
    <Visible[97]>true</Visible>
    <CurrentTimestamp[97]>0</CurrentTimestamp>
    <KeyFrames[97]>
      <KeyFrame>
        <Points>
          <Point>
            <X>1</X>
            <Y>37</Y>
          </Point>
        </Points>
      </KeyFrame>
    </KeyFrames>
  </DynamicAOI>
</ArrayOfDynamicAOI>

```

```
        </Point>
        <Point>
            <X>167</X>
            <Y>345</Y>
        </Point>
    </Points>
    <Angle>0</Angle>
    <Area>51128</Area>
    <Visible>true</Visible>
    <Timestamp>0</Timestamp>
</KeyFrame>
...
</KeyFrames>
</DynamicAOI>
...
</ArrayOfDynamicAOI>
```

Description of Elements:

- **ArrayOfDynamicAOI**: the root element, contains one or more [DynamicAOI](#) ⁹⁶ elements.
- **DynamicAOI**: corresponds to one static AOI and has the following child elements:
- **Name**: defines the name of the AOI
- **Type**: defines the shape of the AOI and should have one of the following values:
 - Rectangle
 - Ellipse
 - Polygon
- **Enabled**: defines the state of the AOI. Disabled AOIs are present only in [AOI Editor](#) ⁷⁵. This element is optional and the implicit value is true.
- **Group**: contains the name of the group. This element is optional and the implicit value is empty.
- **Scope**: defines the scope of the AOI. This element is optional and the implicit value is Local. It should have one of the following values:

- Local
- Global
- **Points:** contains the list of points that defines the AOI and it is dependent of the [type](#)^[96]. The list should contain exactly 2 points for Rectangle or Ellipse, and at least 3 points for Polygon.
- **Angle:** defines the rotation angle of each point defining the AOI around the center of gravity of the AOI. It is expressed in degrees.
- **Color:** defines the color of the pen and brush used to draw the AOI. This element is optional and the implicit value is NamedColor:Black.
- **BorderWidth:** defines the width of the pen used to draw the AOI. This element is optional and the implicit value is 2.
- **Font:** defines the font used to draw the name of the AOI. This element is optional and the implicit values for the child elements are FontName = Tahoma and FontSize = 13.
- **Style:** defines the filling style of the brush used to draw the AOI. This element is optional and the implicit value is HalfTransparent. It should have one of the following values:
 - Hatched
 - Transparent
 - HalfTransparent
- **Transparency:** defines the transparency level (0..100) and is taken into account when the [Style](#)^[97] is HalfTransparent. This element is optional and the implicit value is 50.
- **Area:** the size of the AOI expressed in square pixels
- **Visible:** true if the AOI is visible at the [current timestamp](#)^[97].
- **CurrentTimestamp:** defines the current timestamp.
- **KeyFrames:** defines several key frames made up of [Points](#)^[97], [Visible](#)^[97] and [Timestamp](#)^[97]. The Dynamic AOI position is interpolated in time between the defined key frames.

Examples

The minimal structure that describes a static AOI should look like:

```
<DynamicAOI96>
  <Points97>
    <Point>
      <X>1003</X>
      <Y>748</Y>
    </Point>
    <Point>
      <X>1169</X>
      <Y>886</Y>
    </Point>
  </Points>
  <Type96>Rectangle</Type>
  <Name96>Volvic Logo</Name>
  <Visible97>true</Visible>
</DynamicAOI>
```

The minimal structure that describes a dynamic AOI should look like:

```
<DynamicAOI96>
  <Points97>
    <Point>
      <X>1</X>
      <Y>37</Y>
    </Point>
    <Point>
      <X>167</X>
      <Y>345</Y>
    </Point>
  </Points>
  <Type96>Rectangle</Type>
  <Name96>Bee</Name>
  <Visible97>true</Visible>
  <CurrentTimestamp97>0</CurrentTimestamp>
  <KeyFrames97>
    <KeyFrame>
      <Points>
        <Point>
          <X>1</X>
```

```
        <Y>37</Y>
    </Point>
    <Point>
        <X>167</X>
        <Y>345</Y>
    </Point>
</Points>
<Visible>true</Visible>
<Timestamp>0</Timestamp>
</KeyFrame>
<KeyFrame>
    <Points>
        <Point>
            <X>1</X>
            <Y>60</Y>
        </Point>
        <Point>
            <X>221</X>
            <Y>345</Y>
        </Point>
    </Points>
    <Visible>false</Visible>
    <Timestamp>80000</Timestamp>
</KeyFrame>
</KeyFrames>
</DynamicAOI>
```

6.5 Gaze Replay

6.5.1 Overview

The **Gaze Replay** data view shows gaze positions and eye events for the selected subject plotted over all the stimuli included in the experiment. This is useful to get an overview of the subjects general behavior during the recording of the experiment.

The behavior of this data view is identical to the [Scan Path](#) ^[108] data view (except for the fact that the stimuli are concatenated one after the other in a single playback). For more information on the settings see the [Scan Path Settings](#) ^[113].

The screenshot displays the Gaze Replay software interface. On the left, there is a sidebar with several sections:

- All stimuli**: A table with columns for Subject, Color, ID, and Mean.
- Subjects**: A table with columns for Subject, Color, ID, and Mean.
- Eye Events for Slideshow_Webcam Logtech**: A table with columns for Type, Index, Start Time, and Duration.
- User Events for Slideshow_Webcam Logtech**: A table with columns for Type, Index, Start Time, and Duration.
- Details for Slideshow_Webcam Logtech**: A table with columns for Property and Value.

The main area shows a video player with a red dot indicating the gaze position on a blue device. A red line connects the dot to the device. A white arrow points to the red dot with the label "View Scan Path". A white arrow points to the video player with the label "Animate Gaze Replay". A white arrow points to the video player with the label "Stimulus Start". A white arrow points to the video player with the label "Select Trial Data Set".

Subject	Color	ID	Mean
Slideshow_Webcam Logtech	red	108	108
Slideshow_Webcam Logtech	blue	109	109
Slideshow_Webcam Logtech	green	110	110

Type	Index	Start Time	Duration
Fixation	1	6	99
Fixation	1	86	218
Fixation	1	305	59
Saccade	2	304	139
Saccade	2	504	99
Blink	2	543	99
Blink	3	1001	599

Property	Value
Length	00:00:53.586
Sampling rate	50 Hz
Eye(s)	Both
Number of samples	2332
Number of fixations	159
Number of saccades	186
Number of blinks	124

A specific element of the **Gaze Replay** data view is the automatic insertion of hidden bookmarks in the player control at the beginning of each stimulus to ease the navigation. The usual bookmark navigation keyboard shortcuts apply here ([CTRL] + left/right arrow).

Operate the **Gaze Replay** data view with the following steps:

1. Use the [Stimulus Selection](#)^[55] to change to the desired stimulus.

The [Subjects Selection](#)^[56] displays matching subjects together with their trial gaze data sets.

2. In the [Subjects Selection](#)^[56], activate the desired trial or filter combination.

The [Bee Swarm Main Window](#)^[104] is updated and shows the raw data for the activated trial combination.

While selecting trials, the [Events Selection](#)^[61] view and the [Trial Details](#)^[60] view shows information about the currently selected trial or event.

3. If you click on an event in the **Events** selection view, the corresponding event is automatically selected in the main view.
4. Select the bee swarm time position in the [Thumbnail Control](#)^[70]. Use the [Playback Control](#)^[67] to view an animated bee swarm.
5. You can export the animated scan path display to an AVI file. From the **Export** menu, select the **Export Gaze Replay Video** command.

Alternatively, you can export the current view of the bee swarm to an image file. From the **Export** menu, select the **Save Image to File...** command.



All features of this data view are available with gaze tracking data generated with the iView X™ system. Note that the type of stimuli (still images and/or videos) which can be analyzed depends on the acquired [BeGaze 2.4 program version](#)^[11].




Screen recording experiments and HED videos are only compatible with gaze tracking data which have been generated with the iView X™ version 2.1 or higher.

6.6 Bee Swarm

6.6.1 Overview

The **Bee Swarm** data view shows raw data gaze positions of the selected trial data set plotted on the stimulus image or video.



Operate the **Bee Swarm** data view with the following steps:

1. Use the [Stimulus Selection](#) ^[55] to change to the desired stimulus.

The [Subjects Selection](#) ^[56] displays matching subjects together with their trial gaze data sets.

2. In the [Subjects Selection](#) ^[56], activate the desired trial or filter combination.

The [Bee Swarm Main Window](#) ^[104] is updated and shows the raw data for the activated trial combination.

While selecting trials, the [Events Selection](#) ^[61] view and the [Trial Details](#) ^[60] view shows information about the currently selected trial or event.

3. If you click on an event in the **Events** selection view, the corresponding event is automatically selected in the main view.
4. Select the bee swarm time position in the [Thumbnail Control](#)^[70]. Use the [Playback Control](#)^[67] to view an animated bee swarm.
5. You can export the animated scan path display to an AVI file. From the **Export** menu, select the **Export Bee Swarm Video** command.

Alternatively, you can export the current view of the bee swarm to an image file. From the **Export** menu, select the **Save Image to File...** command.




All features of this data view are available with gaze tracking data generated with the iView X™ system. Note that the type of stimuli (still images and/or videos) which can be analyzed depends on the acquired [BeGaze 2.4 program version](#)^[11].



Screen recording experiments and HED videos are only compatible with gaze tracking data which have been generated with the iView X™ version 2.1 or higher.

6.6.2 Main Data View

The **Bee Swarm** main view visualizes the selected trial data set as a 2D plot over the stimulus image or video. The following image shows an example:



The view shows raw gaze data as colored circles (each color corresponds to a subject).

You can change the bee swarm display with the following steps:

1. Right click the bee swarm display to open a context menu.
2. Select the **Settings** command to display the [Bee Swarm Settings](#) ^[105] dialog. Change settings and confirm with **OK**.

The bee swarm display is updated.

3. Select the **Show AOIs** command, to toggle the visibility of AOIs in the bee swarm display.
4. In the **Export** menu, either select the **Save Image to File** ([CTRL] + [S]) or select the **Copy Image to Clipboard** ([CTRL] + [C]) keyboard command to export the current bee

swarm display to a single image. You can also export the bee swarm to a video file using the **Export Bee Swarm Video** command from the **Export** menu.

Select Gaze Cursor

If you click on gaze cursor in the bee swarm, the clicked subject will be highlighted [Subjects Selection](#)^[56].


Modify subject properties

If required, you can edit the subject properties displayed in the [Subjects Selection](#)^[56] view. Click the desired property and overwrite its content.

6.6.3 Settings

6.6.3.1 View Settings Dialog

In the **View Settings** dialog, you can change the bee swarm display to your needs.



1. Right click the [Bee Swarm Main Window](#)^[104] to open a context menu.
2. Select the **Settings** command to open the **View Settings** dialog.
3. Switch to one of the following tabs and change settings:
 - In the [Bee Swarm Tab](#)^[106] you can change the general appearance of the bee swarm display.
 - In the [Cursor Tab](#)^[107] you configure the gaze cursor appearance.
4. Confirm your settings with **OK**.




If you open a second [Bee Swarm](#)^[102] data view, the new data view will inherit the current view settings. If you adapt the view settings of the second data view, you can switch between the two different bee

swarm views very fast.

6.6.3.2 Bee Swarm Tab

In the **Bee Swarm** tab of the [Bee Swarm Settings](#)^[105] dialog, you configure the general appearance of the bee swarm display.




- **Data channel:** Select if you want to view **Left eye** or **Right eye** data. If the currently selected trail data set only has monocular gaze data, the available data channel is selected automatically.
- **Hide 0 Data:** The gaze tracker produces data with position (0,0) if – for some reason – gaze tracking was lost during the recording. Activate the **Hide 0 Data** option to hide these artifacts. This option is enabled by default.
- **Hide toolbar data:** This option applies to web stimuli only. Activate this check box if you want to hide the gaze data which are located on the web toolbar of the stimulus from the bee swarm.
- **Fade out mouse clicks:** Mouse click events are drawn on the screen

at the moment they took place in the recording. This settings enables the drawing to fade out while playing after it first appears.

6.6.3.3 Cursor Tab

In the **Cursor** tab of the [Bee Swarm Settings](#) ^[105] dialog, you configure the gaze cursor appearance.



- **Gaze cursor:** Configures the appearance of the shape that shows the current gaze position. You can switch between a **Crosshair**, a **Circle**, and a **Translucent dot** shape.

It is also possible to use a 64x64 pixel bitmap as customized shape. Switch to **Cursor image** and click the **Choose...** button to select a suitable external bitmap graphics file.

- **Line width** (not used with **Cursor image** setting): Changes the line width of the gaze cursor (in pixels).
- **Size** (not used with **Cursor image** setting): Changes the diameter of


the gaze cursor (in pixels).

- **Color** (not used with **Cursor image** setting): Changes the gaze cursor color. Click the drop-down icon and select the desired color.

6.7 Scan Path

6.7.1 Overview

The **Scan Path** data view shows gaze positions and eye events of the selected trial data set plotted on the stimulus image or video.



Operate the **Scan Path** data view with the following steps:

1. Use the [Stimulus Selection](#) ^[55] to change to the desired stimulus.

The [Subjects Selection](#) ^[56] displays matching subjects together with their trial gaze data sets.

2. In the [Subjects Selection](#) ^[56], activate the desired trial or filter

combination.

The [Scan Path Main Window](#)^[110] is updated and shows the scan path for the activated trial combination.

While selecting trials, the [Events Selection](#)^[61] view and the [Trial Details](#)^[60] view shows information about the currently selected trial or event.

3. If you click on an event in the **Events** selection view, the corresponding event is automatically selected in the main view.
4. Select the scan path time position in the [Thumbnail Control](#)^[70]. Use the [Playback Control](#)^[67] to view an animated scan path.
5. You can export the animated scan path display to an AVI file. From the **Export** menu, select the **Export Scan Path Video** command.

Alternatively, you can export the current view of the scan path to an image file. From the **Export** menu, select the **Save Image to File...** command.




All features of this data view are available with gaze tracking data generated with the iView X™ system. Note that the type of stimuli (still images and/or videos) which can be analyzed depends on the acquired [BeGaze 2.4 program version](#)^[11].



Screen recording experiments and HED videos are only compatible with gaze tracking data which have been generated with the iView X™ version 2.1 or higher.

6.7.2 Main Data View

The **Scan Path** main view visualizes the selected trial data set as a 2D plot over the stimulus image or video. The following image shows an example for a fixation and saccade plot with dynamic fixation radius and AOIs:



Generally, you can select to plot either raw data or to plot fixations and saccades. If you select to plot fixations and saccades, a fixation point is displayed in the center of a circle and the saccades are plotted as connecting lines in-between. It is also possible to configure a fixed circle radius or a circle radius that relates to the fixation duration. A fixation counter can also be displayed in the center of the fixation circle.

You can change the scan path display with the following steps:

1. Right click the scan path display to open a context menu.

2. Select the **Settings** command to display the [Scan Path Settings](#) ^[113] dialog. In the **Scan Path** tab, select between **Fixations** or **Raw data** display. Change other settings as well and confirm with **OK**.

The scan path display is updated.

3. Select the **Show AOIs** command, to toggle the visibility of AOIs in the scan path display.
4. In the **Export** menu, either select the **Save Image to File** ([CTRL] + [S]) or select the **Copy Image to Clipboard** ([CTRL] + [C]) keyboard command to export the current scan path display to a single image. You can also export the scan path to a video file using the **Export Scan Path Video** command from the **Export** menu.

Select Events

If you click on a fixation circle or on a saccade line, the clicked item will be highlighted. At the same time the corresponding subject and event will be highlighted in the [Subjects Selection](#) ^[56] and the [Events Selection](#) ^[61]. The subject and event will be highlighted when clicking on raw data cursors also.

Highlighted event in the **Eye Events** selection:

Eye Events for cv5 - Trial003				
Type	Index	Start (ms)	Duration (ms)	
Fixation	1	4	298	
Saccade	1	302	39	
Fixation	2	342	278	
Saccade	2	620	19	
Fixation	3	640	477	
Saccade	3	1118	19	
Fixation	4	1138	318	

Highlighted fixation in the **Scan Path** display:



i The scan path is drawn in the color of the corresponding subject unless special timers are defined in the [Scan Path Settings](#) ¹¹³.

Modify subject properties

If required, you can edit the subject properties displayed in the [Subjects Selection](#) ⁵⁶ view. Click the desired property and overwrite its content.

6.7.3 Settings

6.7.3.1 View Settings Dialog

In the **View Settings** dialog, you can change the scan path display to your needs.




1. Right click the [Scan Path Main Window](#)^[110] to open a context menu.
2. Select the **Settings** command to open the **View Settings** dialog.
3. Switch to one of the following tabs and change settings:
 - In the [Scan Path Tab](#)^[113] you can change the general appearance of the scan path display.
 - In the [Cursor Tab](#)^[116] you configure the gaze cursor appearance.
 - In the [Fixations Tab](#)^[117] you adapt the fixations display (tab is inactive if "raw data" is selected in the Scan Path Tab).
4. Confirm your settings with **OK**.



If you open a second [Scan Path](#)^[108] data view, the new data view will inherit the current view settings. If you adapt the view settings of the second data view, you can switch between the two different scan path views very fast.

6.7.3.2 Scan Path Tab

In the **Scan Path** tab of the [Scan Path Settings](#)^[113] dialog, you configure the general appearance of the scan path display.



- **Display:** Select if you want to view **Fixations** or **Raw data**. To view saccades as well, enable the **Trailer** option (see below).
- **Data channel:** Select if you want to view **Left eye** or **Right eye** data. If the currently selected trail data set only has monocular gaze data, the available data channel is selected automatically.
- **Hide 0 Data:** The gaze tracker produces data with position (0,0) if – for some reason – gaze tracking was lost during the recording. Activate the **Hide 0 Data** option to hide these artifacts. This option is enabled by default.

- **Hide toolbar data:** This option applies to web stimuli only. Activate this check box if you want to hide the gaze data which are located on the web toolbar of the stimulus from the scan path.
- **Fade out mouse clicks:** Mouse click events are drawn on the screen at the moment they took place in the recording. This settings enables the drawing to fade out while playing after it first appears.
- **Line width:** Select the line widths for the scan path lines (in pixels).
- **Draw connection lines:** Activate this option, if raw data should be connected with lines. This option is enabled by default.
- **Time interval:** You can define two intervals in which the scan path should be plotted in a different color. After these intervals ended, the scan path plot continues with the defined subject color property in the **Subjects** list view. Activate the **Make standard intervals** option if the scan path plot should continue with alternating intervals according to the time interval definition.
- **Trailer:** Determines, how many gaze data is accumulated to display fixations and saccades. Note that the following settings relate to the time window you have set in the [Thumbnail Control](#) ^[70].

From beginning (still image stimulus only): If activated, all gaze data is displayed from the first sample to the current analysis position.


Constant length: If activated, the current analysis position leaves "a trail behind". This means: a certain time window of gaze data – which immediately precedes the current analysis position – is displayed. Use the slider to change the length of time window from 0 seconds up to 10 seconds.



If you display an overlay of the real-time gaze positions of multiple subjects, this is called the "bee swarm" mode. To activate this display mode, enable the Raw Data display and configure the trailer with a Constant length of zero. Select multiple subjects / trials and press play.

6.7.3.3 Cursor Tab

In the **Cursor** tab of the [Scan Path Settings](#) ^[113] dialog, you configure the gaze cursor appearance.




- **Gaze cursor:** Configures the appearance of the shape that shows the current gaze position. You can switch between a **Crosshair**, a **Circle**, and a **Translucent dot** shape.

It is also possible to use a 64x64 pixel bitmap as customized shape. Switch to **Cursor image** and click the **Choose...** button to select a suitable external bitmap graphics file.

- **Line width** (not used with **Cursor image** setting): Changes the line width of the gaze cursor (in pixels).
- **Size** (not used with **Cursor image** setting): Changes the diameter of the gaze cursor (in pixels).
- **Color** (not used with **Cursor image** setting): Changes the gaze cursor color:
 - **Subject**: sets the gaze cursor color to the subject color property in the **Subjects** list view. This is the default selection.
 - **Same for all**: Click the drop-down icon and select the desired color to use for the gaze cursor.

6.7.3.4 Fixations Tab

In the **Fixations** tab of the [Scan Path Settings](#)^[113] dialog, you configure how fixations are plotted on the scan path display. The following settings only apply if you have activated the **Fixations** option in the [Scan Path Settings – Scan Path Tab](#)^[113].



- **Shape:** Selects between a **Crosshair** and a **Circle** shaped fixation display.
- **Size:** Determines the fixation shape size.
Constant size: If checked, the size of the fixation shapes is constant. You can change the shape's size (in pixels).
Raindrop: If checked, the size of the fixation shape is proportional to the fixation duration. On the slider, you can set how many pixels represent a 500 ms fixations.
- **Use raindrops in playback:** If checked, the radius of the fixation

shapes also changes during replay or while moving the current analysis position.

- **Fill Mode:** Selects the fixation shape fill mode: **Hatched**, **Half Transparent** or **Transparent** fills are supported.
- **Highlight:** Selects the highlight color for the fixation shape. Click the drop-down icon and select the desired color.
- **Show fixation counter:** Counts up the fixations and indicates a counter for each fixation.

6.8 Focus Map

6.8.1 Overview

With the **Focus Map** data view, gaze patterns are visualized by altering the transparency of the stimulus display based on the amount of attention received.

The screenshot illustrates the BeGaze 2.4 software interface during a Focus Map analysis. The interface is divided into several key sections:

- Stimulus Selection:** Located at the top left, it shows the current stimulus as 'image13.bmp' and provides a 'Change Stimulus' button.
- Subjects List:** Below the stimulus, a table lists subjects with columns for 'Subjects (2/2)', 'Gender', and 'Age'. The selected subject is 'cv4 - Trial003' with a 'distus' value of 342.
- Eye Events Table:** A table showing gaze events for 'cv4 - Trial003'. The table has columns for 'type', 'fixes', 'start [ms]', and 'duration [ms]'.

type	fixes	start [ms]	duration [ms]
Fixation	1	37	338
Saccade	1	355	39
Fixation	2	395	407
Saccade	2	852	30
Fixation	3	872	159
Saccade	3	1032	19
Fixation	4	1052	170
- User Events Table:** A table showing user events for 'cv4 - Trial003'.

Property	Value
Sample	66-00-04-992
Sampling rate	50 Hz
Fixes	Both
Number of samples	250
Number of fixations	20
Number of saccades	18
Number of blinks	0
- Focus Map Visualization:** The central area displays a heatmap overlay on a stimulus image (a bottle), indicating areas of high attention (red/orange) and low attention (blue).
- Timeline and Playback:** At the bottom, a timeline shows the sequence of frames with a play button and a 'Change Time Position' control.

Operate the **Focus Map** data view with the following steps:

1. Use the [Stimulus Selection](#) ^[55] to change to the desired stimulus.

The [Subjects Selection](#) ^[56] displays matching subjects together with their trial gaze data sets.

2. In the [Subjects Selection](#) ^[56], activate the desired trial or filter combination.

The [Focus Map Main Window](#) ^[121] is updated and shows the focus map for the activated trial combination.

While selecting trials, the [Events Selection](#) ^[61] view and the [Trial Details](#) ^[60] view shows information about the currently selected trial or event.

3. If you click on an event in the **Eye vents** selection view, the corresponding event is automatically selected in the main view.
4. Select the focus map time position in the [Player Control](#) ^[66]. Use the [Playback Control](#) ^[67] to view an animated attention map.
5. You can export the animated focus map display to an AVI file. From the **Export** menu, select the **Export Focus Map Video** command.

Alternatively, you can export the current view of the attention map to an image file. From the **Export** menu, select the **Save Image to File...** command.



All features of this data view are available with gaze tracking data generated with the iView X™ system. Note that the type of stimuli (still images and/or videos) which can be analyzed depends on the acquired [BeGaze 2.4 program version](#) ^[11].




Screen recording experiments and HED videos are only compatible with gaze tracking data which have been generated with the iView X™ version 2.1 or higher.

6.8.2 Main Data View

After selecting the desired trial data, the **Focus Map** main view displays the updated map. Two visualization styles are possible:

- The **Focus map** shows fixation hits related to brightness between darkest (less hits) and normal brightness (most hits).



- The **Custom map** shows fixation hits related to a custom defined color scale.



Note, that the data interpretation differs with the stimulus type. The map displayed for a still image stimulus is based on fixations while the map displayed for a video stimulus is based on raw data.

Change the focus map display

To change the focus map display settings proceed as follows:

1. Right click the map display to open a context menu.
2. Select the **Settings** command to display the [Focus Map Settings](#) dialog. Select the map style and confirm with **OK**.

The focus map display is updated.


3. Select the **Show AOIs** command, to toggle the visibility of AOIs in the map display.
4. In the **Export** menu, either select the **Save Image to File** ([CTRL] + [S]) or select the **Copy Image to Clipboard** ([CTRL] + [C]) keyboard command to export the current focus map display to a single image. You can also export the focus map to a video file using the **Export Focus Map Video** command from the **Export** menu.

Modify subject properties

If required, you can edit the subject properties displayed in the [Subjects Selection](#) view. Click the desired property and overwrite its content.

6.8.3 Settings

In the **View Settings** dialog, you can configure the visualization style and parameters of the **Focus Map**.



General Settings

- **Data channel:** Select if you want to view left or right eye. In case of monocular gaze data files, the available data channel is selected automatically.
- **Saturation:** Select if you want to see a **Focus map** with the maximum saturation value computed dynamically for each trial or a **Custom map**.
- **Fade out mouse clicks:** Mouse click events are drawn on the screen at the moment they took place in the recording. This settings enables the drawing to fade out while playing after it first appears.



Note that the saturation and color parameters settings in the dialog are available for the Custom map only.

Parameters

- **Opaque Color:** The overlay background color used for unfocused areas (default is black)
- **Data Range (min..max):** For every pixel displayed on the map, the fixation duration is counted and integrated over time. For multiple subjects, the sum (over all subjects) of the fixation duration is calculated. The double slider defines the minimum and maximum duration of the scale.

If the maximum value is reached or exceeded the matching image pixels will be drawn with the highest value, which is

- normal brightness for the Focus map,
- a customized color for Custom map style

If the minimum value is not reached, the matching image pixels will be drawn with the lowest value, which is

- no brightness for the Focus Map (or the selected opaque color if changed from black),
- a customized color for the Custom Map.

Changing this parameter is useful if you are interested in fixations that exceed a specific fixation duration.

- Use the **Opacity** double slider to change the opacity level for the corresponding minimum and maximum data range values above.
- **Kernel width:** To calculate the Focus Map, all fixation hits are filtered with a Gaussian filter. This setting defines the width (in pixels) of the Gaussian curve. If you decrease the value, the analysis resolution will increase. At the same time, the hot spots will become smaller and more spread.
- **Trailer:** Determines, how many gaze data is accumulated to display fixations. Note that the following settings relate to the time window you

have set in the [Thumbnail Control](#)⁷⁰.


From Start (still image stimulus only): If selected, all gaze data is displayed from the first sample to the current analysis position.

Constant length: If selected, the current analysis position leaves "a trail behind". This means: a certain time window of gaze data – which immediately precedes the current analysis position – is displayed. Use the slider to change the length of time window from 0 seconds up to 10 seconds.

6.9 Heat Map

6.9.1 Overview

With the **Heat Map** data view, gaze patterns are visualized by altering the color of the stimulus display based on the amount of attention received.



Operate the **Heat Map** data view with the following steps:

1. Use the [Stimulus Selection](#)⁵⁵ to change to the desired stimulus.

The [Subjects Selection](#) ^[56] displays matching subjects together with their trial gaze data sets.

2. In the [Subjects Selection](#) ^[56], activate the desired trial or filter combination.

The [Heat Map Main Window](#) ^[126] is updated and shows the heat map for the activated trial combination.

While selecting trials, the [Events Selection](#) ^[61] view and the [Trial Details](#) ^[60] view shows information about the currently selected trial or event.

3. If you click on an event in the **Events** selection view, the corresponding event is automatically selected in the main view.
4. Select the heat map time position in the [Player Control](#) ^[66]. Use the [Playback Control](#) ^[67] to view an animated heatmap.
5. You can export the animated heat map display to an AVI file. From the **Export** menu, select the **Export Heat Map Video** command.

Alternatively, you can export the current view of the heat map to an image file. From the **Export** menu, select the **Save Image to File...** command.



All features of this data view are available with gaze tracking data generated with the iView X™ system. Note that the type of stimuli (still images and/or videos) which can be analyzed depends on the acquired [BeGaze 2.4 program version](#) ^[11].



Screen recording experiments and HED videos are only compatible with gaze tracking data which have been generated with the iView X™ version 2.1 or higher.

6.9.2 Main Data View

After selecting the desired trial data, the **Heat Map** main view displays the updated map. Two visualization styles are possible:

- The **Heat map** shows fixation hits related to the color scale between

blue (less hits) and red (most hits).



- The **Custom** map shows fixation hits related to a custom defined color scale.



Note, that the data interpretation differs with the stimulus type. The map displayed for a still image stimulus is based on fixations while the map displayed for a video stimulus is based on raw data.

Change the heat map display

To change the heat map display settings proceed as follows:

1. Right click the map display to open a context menu.
2. Select the **Settings** command to display the [Heat Map Settings](#) ¹²⁸ dialog. Select the map style and confirm with **OK**.

The Heat map display is updated.


3. Select the **Show AOIs** command, to toggle the visibility of AOIs in the map display.
4. In the **Export** menu, either select the **Save Image to File** ([CTRL] + [S]) or select the **Copy Image to Clipboard** ([CTRL] + [C]) keyboard command to export the current heat map display to a single image. You can also export the heat map to a video file using the **Export Heat Map Video** command from the **Export** menu.

Modify subject properties

If required, you can edit the subject properties displayed in the [Subjects Selection](#) ^[56] view. Click the desired property and overwrite its content.

6.9.3 Settings

In the **View Settings** dialog, you can configure the visualization style and parameters of the **Heat Map**. The available settings are identical to the ones in the **Focus Map** except for the coloring selection which is described below (and replaces the Opaque color setting in Focus Map). For a detailed description of the common settings see [Focus Map Settings](#) ^[122].



- **Color coding:** select between predefined 3-color and 2-color codings and a user defined 3-color coding for the heat map. The heat map is colored with the selected range of colors starting with the left color for the shortest fixations and ending with the right color for the longest ones.

See also:


[» Heat Map Main Window](#) ^[126]

[» Heat Map Overview](#) ^[125]

6.10 Key Performance Indicators

6.10.1 Overview

With the **Key Performance Indicators** data view, a number of important statistical indicators are visualized in text bubbles associated to each AOI. The statistical data is updated in realtime and reflects the selected subjects in the Subjects list view.



Operate the **Key Performance Indicators** data view with the following steps:

1. Use the [Stimulus Selection](#) ^[55] to change to the desired stimulus.
The [Subjects Selection](#) ^[56] displays matching subjects together with their trial gaze data sets.
2. In the [Subjects Selection](#) ^[56], activate the desired trial or filter combination.
The [Key Performance Indicators Main Window](#) ^[131] is updated and

shows the KPIs for the activated trial combination.

While selecting trials, the [Events Selection](#)^[61] view and the [Trial Details](#)^[60] view shows information about the currently selected trial or event.

3. If you click on an event in the **Events** selection view, the corresponding event is automatically selected in the main view.
4. Select the KPI time position in the [Player Control](#)^[66]. Use the [Playback Control](#)^[67] to view the KPIs in real time.
5. You can export the animated KPI display to an AVI file. From the **Export** menu, select the **Export KPIs Video** command.

Alternatively, you can export the current view of the KPIs to an image file. From the **Export** menu, select the **Save Image to File...** command.



All features of this data view are available with gaze tracking data generated with the iView X™ system. Note that the type of stimuli (still images and/or videos) which can be analyzed depends on the acquired [BeGaze 2.4 program version](#)^[11].



Screen recording experiments and HED videos are only compatible with gaze tracking data which have been generated with the iView X™ version 2.1 or higher.


6.10.2 Main Data View

The **Key Performance Indicators** (KPI) main view gets you immediate responses at a glance:

- Which stimuli elements were the eye catchers?
- How many subjects watched each element?
- In which order?
- How many revisits?
- What is the rank and share of visual attention among all subjects?

- and other indicators

It makes the results quantitative and visible.



KPI functionalities and handling

- Works with still images and video clips, on websites or screen recording videos
- Displayed as overlay on Areas of Interest (AOI) visualization
- Interactive information updated based on selected subjects (individual, groups, all) and time of regard
- Select and deselect KPI windows, move their position freely
- Export visualization as BMP or AVI for your exposé, report, documentation etc.
- A White Space KPI exists for still image stimuli only and shows indicators for the area left outside of the AOIs

Change the KPI display

To change the KPI display settings proceed as follows:

1. Right click the main view to open a context menu.
2. Select the **Settings** command to display the [KPI Settings](#) ¹³³ dialog. Select the indicators to display and confirm with **OK**.

The KPI display is updated.

3. In the **Export** menu, either select the **Save Image to File** ([CTRL] + [S]) or select the **Copy Image to Clipboard** ([CTRL] + [C]) keyboard command to export the current KPI display to a single image. You can also export the KPIs to a video file using the **Export KPIs Video** command from the **Export** menu.

Modify subject properties

If required, you can edit the subject properties displayed in the [Subjects Selection](#) ⁵⁶ view. Click the desired property and overwrite its content.

6.10.3 Settings

In the **View Settings** dialog, you can select which indicators to show in the **Key Performance Indicators** data view.

Data Channel: Left Eye Right Eye

Fade out mouse clicks

Show KPI(s): Sequence

Entry Time

Dwell Time

Hit Ratio

Revisits

Revisitors

Average Fixation

First Fixation

Fixation Count

None

None

None

Font Size: 100%

Fixed Size

OK Cancel

General Settings

- **Data channel:** Select if you want to view left or right eye. In case of monocular gaze data files, the available data channel is selected automatically.
- **Fade out mouse clicks:** Mouse click events are drawn on the screen at the moment they took place in the recording. This settings enables the drawing to fade out while playing after it first appears.

Indicators

The available key performance indicators and their meaning are described in the table below:

KPI Name	unit	Description
Sequence	count	Order of gaze hits into the AOIs based on Entry time, lowest entry time = first in Sequence
Entry time	ms	Average duration for the first fixation into the AOI
Dwell time	ms and %	Dwell time average ms = sum (all fixations and saccades within an AOI for all selected subjects) / by number of selected subject Dwell time average % = dwell time average * 100 / (current time - start time)
Hit ratio	count and %	How many subjects out of the selected subjects looked at least one time into the AOI - "total hit count" / "number of selected subjects"
Revisits	count	Average Revisits = (Number of glances divided by selected subjects with at least one visit) -1 Glances = Increments the counter each time a fixation hits the AOI if not hit before
Revisitors	count	1. Number of subjects with more than one visit in an AOI 2. Total number of subjects with at least one visit into an AOI e.g. 3 revisitors out of 7 visitors

Average Fixation	ms and %	Sum of "average fixation time per subject in an AOI" divided by number of selected subjects
First fixation	ms	Sum of all "first fixations" for selected subjects divided by number of selected subjects
Fixation count	count	Number of all fixations for selected subjects divided by number of selected subjects

Additionally there are three combo-boxes that allow to select two more indicators (one each) to show together with the ones in the table above. For the description of these parameters see the [AOI Summary Statistics](#) ^[175] list. The available extra indicators can be the following:

- AOI Area
- AOI Coverage
- Glance Duration
- Diversion Duration
- Appearance Count
- Visible Time
- Net Dwell Time


Font

- **Font Size:** Selects the size of the KPIs font as a percent of the standard font size used for the main view (the font size used for AOI names in the AOI Editor for example).
- **Fixed Size:** If checked the KPI font size remains the same at all zoom levels, otherwise the font size scales together with the AOIs at different zoom levels. Default is not checked.

6.11 Gridded AOIs

6.11.1 Overview

With the **Gridded AOIs** (aka content independent AOIs) data view, gaze patterns and statistics parameters are visualized by altering the color of a grid of AOIs overlaid over the stimulus based on the amount of attention received. Gridded AOI maps allows complementary interpretation to heat maps – qualitative and quantitative - and allows the comparison of different stimuli independent of their content.



Operate the **Gridded AOIs** data view with the following steps:

1. Use the [Stimulus Selection](#) ^[55] to change to the desired stimulus.

The [Subjects Selection](#) ^[56] displays matching subjects together with their trial gaze data sets.

2. In the [Subjects Selection](#) ^[56], activate the desired trial or filter combination.

The [Gridded AOIs Main Window](#)^[139] is updated and shows the gridded AOIs for the activated trial combination.

While selecting trials, the [Events Selection](#)^[61] view and the [Trial Details](#)^[60] view shows information about the currently selected trial or event.

3. If you click on an event in the **Events** selection view, the corresponding event is automatically selected in the main view.
4. Select the gridded AOIs time position in the [Player Control](#)^[66]. Use the [Playback Control](#)^[67] to view an animated heatmap.
5. You can export the animated gridded AOIs display to an AVI file. From the **Export** menu, select the **Export Gridded AOIs Video** command.

Alternatively, you can export the current view of the gridded AOIs to an image file. From the **Export** menu, select the **Save Image to File...** command.



All features of this data view are available with gaze tracking data generated with the iView X™ system. Note that the type of stimuli (still images and/or videos) which can be analyzed depends on the acquired [BeGaze 2.4 program version](#)^[11].



Screen recording experiments and HED videos are only compatible with gaze tracking data which have been generated with the iView X™ version 2.1 or higher.

6.11.2 Main Data View

The **Gridded AOIs** main view visualizes the selected trial data set as a rectangular AOIs grid over the stimulus image or video. The AOIs in the grid show various statistical values like Entry Time, Dwell Time, Revisits and more. The following image shows an example for an 8x8 grid using the Average Entry Time as parameter in milliseconds:



You can change the gridded AOIs display with the following steps:

1. Right click the gridded AOIs display to open a context menu.
2. Select the **Settings** command to display the [Gridded AOIs Settings](#) ¹⁴³ dialog. Select the number of rows and columns for the AOI grid. Change the displayed statistics parameter as well and confirm with **OK**.

The AOI grid is updated.

- In the **Export** menu, either select the **Save Image to File** ([CTRL] + [S]) or select the **Copy Image to Clipboard** ([CTRL] + [C]) keyboard command to export the current gridded AOIs display to a single image. You can also export the gridded AOIs to a video file using the **Export Gridded AOIs Video** command from the **Export** menu.

Parameters


The **Gridded AOIs** view can display one of the following statistics parameters:

- Entry Time (Average)
- Dwell Time (Total)
- Dwell Time (Average)
- Revisits
- Fixation Count (Total)
- Fixation Count (Average)
- Subject Hit
- Sequence (Average)

The displayed parameter can be changed from the **Parameter** drop-down box in [Gridded AOIs Settings](#)^[143].

Export Statistics

If you right click on the gridded AOIs display the context menu is displayed and the option to **Export Statistics** can be selected. This exports to file or to clipboard all the AOI parameters (name, area) and all the statistics parameters that can be displayed in the gridded AOIs view.



Export Scan Path Strings

Please see [Scanpath String](#) ^[141].

SPSS case format

Checking the **Use SPSS case format** changes the output format so that it contains a single line of data instead of having each trial data on its own line. This is useful to group the data for so called "cases" in SPSS.

Modify subject properties

If required, you can edit the subject properties displayed in the [Subjects Selection](#) ^[56] view. Click the desired property and overwrite its content.


6.11.3 Scan Path Strings

Scanpath strings are used in research to measure scanpath similarities (e.g. Levenshtein distance measure, ClustalG method)

When the scanpath runs over the gridded AOIs, each fixation is replaced by the name of the AOI hit.

Export Scan Path Strings

Selecting the **Export Scanpath Strings...** from the context menu allows to export to file the scanpath string for each trial in the experiment. The scanpath string represents the sequence of AOIs in the grid that the scan path has fixations in. See the [Scan Path](#)^[110] description for more details.



Raw scanpath strings

An AOI in the grid is represented as a letter-number combination representing the row and the column of that particular AOI. The rows are labeled left to right as A, B, C and so on and the columns top to bottom are 1, 2, 3... So a scanpath string can look like this: F5-C5-C4. This shows that the scan path for that trial had fixations in order in AOIs F5, C5 and C4. This string is called the *raw scanpath string*.

Compressed scanpath string


Additionally a *compressed scanpath string* is also exported. The compressed string is obtained by eliminating duplicated consecutive AOIs (A1A1 becomes A1) and duplicated sequences (A1-B1-C1-A1-B1 becomes A1-B1-C1).

The compressed string is obtained by eliminating duplicated consecutive

AOIs (A1-A1 becomes A1) and duplicated sequences (A1-B1-C1-A1-B1 becomes A1-B1-C1). As described in <http://research.chtsai.org/papers/scanpath-compression.html>

6.11.4 Settings

In the **View Settings** dialog, you can select which indicators to show in the **Gridded AOIs** data view.



General Settings

- **Data channel:** Select if you want to view left or right eye. In case of monocular gaze data files, the available data channel is selected automatically.

- **Fade out mouse clicks:** Mouse click events are drawn on the screen at the moment they took place in the recording. This settings enables the drawing to fade out while playing after it first appears.

Grid Configuration

- **Rows:** number of rows for the generated AOI grid
- **Columns:** number of columns for the generated AOI grid
- **Show borders:** display the grid lines between AOIs
- **Show value:** display the values of the selected statistics parameter inside the AOIs
- **Opacity:** selects the opacity level of the AOI grid colors

Parameter

The available parameters to be displayed and their meaning are described in the table below:

KPI Name	unit	Description
Entry time (Average)	ms	Average duration for the first fixation into the AOI
Dwell time (Total)	ms	Dwell time ms = sum (all fixations and saccades within an AOI for all selected subjects)


Dwell time (Average)	ms	Dwell time average ms = sum (all fixations and saccades within an AOI for all selected subjects) / by number of selected subjects
Revisits	count	Average Revisits = (Number of glances divided by selected subjects with at least one visit) -1 Glances = Increments the counter each time a fixation hits the AOI if not hit before
Fixation count (Total)	count	Number of all fixations for selected subjects
Fixation count (Average)	count	Number of all fixations for selected subjects divided by number of selected subjects
Subject Hit	count	Number of subjects that looked into the AOI
Sequence (Average)	count	How many subjects out of the selected subjects looked at least one time into the AOI - "total hit count" / "number of selected subjects"

These parameteres are among those found in the [AOI Summary Statistics](#) ^[175] list.

6.12 AOI Sequence Chart

6.12.1 Overview

The **AOI Sequence Chart** shows the temporal order at which AOIs were hit by a particular subject.



Operate the **AOI Sequence Chart** data tab with the following steps:

1. Use the [Stimulus Selection](#)^[55] to change to the desired stimulus.

The [Subjects Selection](#)^[56] displays matching subjects together with their trial gaze data sets.

2. In the [Subjects Selection](#)^[56], select one or multiple trials.

The [AOI Sequence Chart Main View](#)^[147] is updated and shows the AOI hits for the selected trial.


While selecting trials, the [Trial Details](#)^[60] view shows information about the currently selected trial.

3. Configure the [Chart Display Modes](#) ^[74] to further adapt the display to your needs.

6.12.2 Main Data Tab

Single Subject Selection


After selecting the desired trial data, the **AOI Sequence Chart** main view displays the updated chart. The following image shows the display for a still image stimulus.



The colored bars represent the different AOIs hits. If the AOIs are labeled, their names appear at the y-axis. The x-axis shows the time in milliseconds. If you right click on one of the bars, a tooltip will pop up displaying detailed information on the AOI (name, start / end time of its presentation, and the duration of the AOI presentation).

In the example above the selected subject was looking at the AOI labeled "Cap" (colored in blue), then the gaze switches to the AOI labeled "Corona on bottle" (colored in violet).


For video stimuli, the display also indicates when a specific AOI has the visibility property set. In the example below, the AOI labeled "Bee" is visible ("active") from start until the 24th second while the AOI labeled "SMS" starts invisible ("not active") and gets visible around the 16th second.



You can change the AOIs and also change the AOI colors in the [AOI Editor](#)⁷⁵.

Multiple Subject Selection

After selecting the desired trial data, the **AOI Sequence Chart** main view displays the updated chart. The representation is the same for still images and video stimuli.



The colored bars represent the different AOIs hits. If the AOIs are labeled, their names appear at the Legend. The x-axis shows the time in milliseconds. If you right click on one of the bars, a tooltip will pop up

displaying detailed information on the AOI (name, start / end time of it's presentation, and the duration of the AOI presentation).

In the example above the selected subject was looking at the AOI labeled "Cap" (colored in blue), then the gaze switches to the AOI labeled "Corona on bottle" (colored in violet).

Click the **Reset Scaling** icon in the top left corner to revert display scaling and positioning.

Click the **Legend** button in the top right corner to hide or unhide the legend.


Modify subject properties

If required, you can edit the subject properties displayed in the [Subjects Selection](#) view. Double click the desired property and overwrite its content.

6.13 Binning Chart

6.13.1 Overview

The **Binning Chart** shows a statistical overview of AOI hits for separated time slices (bins). For each time slice, the AOI hit percentages for all selected trials are summarized and displayed as stacked column.



Operate the **Binning Chart** data view with the following steps:

1. Use the [Stimulus Selection](#) ^[55] to change to the desired stimulus.

The [Subjects Selection](#) ^[56] displays matching subjects together with their trial gaze data sets.

2. In the [Subjects Selection](#) ^[56], activate the desired trial or filter combination.

The [Binning Chart Main Window](#) ^[151] is updated and shows the AOI hit percentages for the activated trial combination.

While doing this, the [Trial Details](#) ^[60] view shows information about the

currently selected trial.


- Configure the [Chart Display Modes](#) ^[74] to further adapt the display to your needs.



You can change the time slice granulation in the configuration area available below the main display area. You can change the Bins integration time [ms] setting from sampling frequency (e.g. 20ms for 50Hz data) up to 60 seconds.

6.13.2 Main Data Tab

After selecting the desired trial data, the **Binning Chart** main view displays the updated chart.



The AOI hit percentages are presented using different colors. The legend below the chart shows which colors are used.

In the above example between the 20th and 21st second the "text" AOI was hit at about 14%, whereas all other AOIs were not hit in this time slice. In the next second another AOI ("small logo") was also hit.



You can change the AOIs and also change the AOI colors in the [AOI Editor](#) ¹⁷⁵.

Modify subject properties

If required, you can edit the subject properties displayed in the [Subjects Selection](#) ¹⁵⁶ view. Double click the desired property and overwrite it's content.

6.14 Event and Reading Statistics

6.14.1 Overview

The **Event Statistics** and **Reading Statistics** data tabs presents information and statistics regarding gaze tracking events. The data view's main view consists of different parts identified in the image below.

Trial	Subject	Color	Gender	Stimulus	Start Time [ms]	End Time [ms]	Fixation Start [ms]
Trial002	o5	Coral	male	image11.bmp	0	4980	5
Trial002	o5	Coral	male	image11.bmp	0	4980	293
Trial002	o5	Coral	male	image11.bmp	0	4980	601
Trial002	o5	Coral	male	image11.bmp	0	4980	940
Trial002	o5	Coral	male	image11.bmp	0	4980	1417
Trial002	o5	Coral	male	image11.bmp	0	4980	1955
Trial002	o5	Coral	male	image11.bmp	0	4980	2074
Trial002	o5	Coral	male	image11.bmp	0	4980	2452
Trial002	o5	Coral	male	image11.bmp	0	4980	2751
Trial002	o5	Coral	male	image11.bmp	0	4980	3228
Trial002	o5	Coral	male	image11.bmp	0	4980	3686
Trial002	o5	Coral	male	image11.bmp	0	4980	4144
Trial002	o5	Coral	male	image11.bmp	0	4980	4522
Trial002	o4	CornflowerBlue	female	image11.bmp	0	5014	33
Trial002	o4	CornflowerBlue	female	image11.bmp	0	5014	317
Trial002	o4	CornflowerBlue	female	image11.bmp	0	5014	735
Trial002	o4	CornflowerBlue	female	image11.bmp	0	5014	1352
Trial002	o4	CornflowerBlue	female	image11.bmp	0	5014	1949
Trial002	o4	CornflowerBlue	female	image11.bmp	0	5014	2168

You operate the **Event Statistics** and **Reading Statistics** data views with the following steps. While doing so, the [Results Grid](#) ¹⁵⁸ updates in real-time displaying the outcome of your selections and settings.

1. Use the **Selection Tree** displayed to the lower left to select the stimuli, trials, and areas of interest for statistic analysis. To narrow down or

qualify your selection, enable the **Filter** option to display the **Filter Tree** (upper left). See [Statistics Selection Trees](#)^[153] for an in depth explanation.

2. Choose the desired **Statistics Template** from the Statistics selection box. The list offers both predefined and user defined templates. You may duplicate and change a predefined statistics template. See [Statistics Template](#)^[155] for an in-depth explanation.
3. Press **Settings** button to select or deselect cells from the template, to create own templates and switch between evaluation of **Left eye** or **Right eye** gaze tracking data
4. As an option, you may specify the desired [Time Interval](#)^[157]. Furthermore, it is also possible to re-arrange the columns, sort the data or only show columns of your interest within the [Results Grid](#)^[158].
5. If the display suits your requirements, click **Export** to write the current display to a file. See [Export Statistics](#)^[158] for details.
6. Click on **Copy to Clipboard** button to copy the current shown statistic into the clipboard for further use in other programs, e.g. MS-Excel.



The statistics display is calculated in real-time. Depending on the complexity of the experiment and on the computer performance, the calculation might take some time.



The [Reading Statistics](#)^[199] data view is available when the Reading Package is licensed.

6.14.2 Selection Trees

Selection Tree

The **Selection Tree** is used to select the stimuli, trials and areas of interest for which the **Event Statistics** data view outcome is computed. Using the selection tree is straightforward:

1. The top level (root) nodes selects or de-selects stimuli available in the current experiment. To help in the selection, a thumbnail of the stimulus is displayed as tooltip when you hover the mouse over the respective screen region.
2. If you enable or disable a node, all child nodes follow that selection. For example: to de-select all child entries associated below a specific stimulus, disable the corresponding top level node.
3. On the tree's second level, you select or de-select statistics for all **Areas of interest** or statistic entries for all **Subjects – Trials**. Note, that you can narrow down the selection of subjects and trials with the **Filter Tree** (see below).
4. On the tree's third level, you select a specific combination of AOIs or a specific combination of trials. A "white space" AOI is generated to cover all areas left outside of defined AOIs.

Once a selection is made, the results are computed and displayed in the [Results Grid](#)^[158] immediately.

Filter Tree

With the **Filter Tree**, a specific set of trials / subjects can be selected. This is especially helpful, if you have a large number of trials or if you want to select trials / subjects by additional subject properties collected while running the experiment.

1. Activate the **Filter** option above the **Selection Tree**.

A separate tree view opens. The new tree view lists all **Subjects** as well as customized subject properties as top level experiment. Note, that customized subject properties (for example **Gender** or **Age**) need to be defined when creating the experiment using SMI Experiment Center. When running the experiment, these properties are available for operator input when starting a new trials.

2. Open the available top level nodes and select the desired combination of **Subjects** or customized subject properties. For example: if your experiment includes the subject property **Gender**, you are now able to select trials linked to male or female subjects.

The selected filter combination is applied. The results are computed

and displayed in the [Results Grid](#)^[158] immediately. Note, that the selection in the **Filter Tree** is independent from the selection already made in the **Selection Tree**. For this reason, already de-selected items from the **Selection Tree** may show up in the **Results Grid** now.

3. After doing the selection in the **Filter Tree**, you can de-select items in the **Selection Tree** to temporarily hide specific items from the **Results Grid**.
4. Deactivate the **Filter** option to switch off the settings made in the **Filter Tree**.

Switch between tooltip view of AOI and AOI preview

1. To switch between the tooltip view of an AOI and the AOI preview, press [CTRL] + [T].


6.14.3 Template List

For optimized handling of the large count of statistical data items, BeGaze 2.4 groups them as **Statistic Templates**. Each statistic template covers a specific purpose. For details about the predefined templates see [Statistics Definitions and Examples](#)^[160].

To operate the statistics templates, proceed as follows:

1. Select an item from the **Statistics** list at the top left of the data view.
This will activate a set of statistic items, which are computed and displayed in the [Results Grid](#)^[158] immediately.
2. After activating the desired template, you can modify the **Results Grid** to suit your needs. This can be done by
 - changing the column selection,
 - changing the column sorting, or by
 - changing the column order.

3. Click the **Settings** button to change the columns selection or to copy the modified settings to a new statistic template.



To save the customized **Statistic Templates** press the "Save As..." button in the settings dialog

4. To remove a customized statistic template, open the settings dialog and click the **Delete** button.
5. Optionally, when the settings dialog is closed, you can ...
 - select the **Save Settings for Experiment** menu command or press the [CTRL] + [E] key combination to save the **Statistic Templates** list to the currently opened experiment or

– select the **Save Settings Globally** menu command or press the [CTRL] + [G] key combination, to save the **Statistic Templates** list for use with other experiments. Note that this command will overwrite a previously saved global list.



It is not possible to delete the default statistic templates.

6.14.4 Time Interval

The settings grouped under **Time Window** limit the data to be evaluated while computing the event statistics. The default setting includes all gaze tracking data currently selected for display in the [Statistics Selection Trees](#)^[153]. Both time settings denote a relative time in milliseconds where each trial starts at zero. You can narrow the time window with the following steps:

1. Enter the starting time in the **Start** input. You can enter a number in milliseconds, which is automatically converted to the hh:mm:ss:ms format. You can also enter the time value in the hh:mm:ss:ms format where **hh** denotes a two digit hour value, **mm** denotes minutes, **ss** denotes seconds, and **ms** denote milliseconds.
All gaze tracking data before this time will be filtered out.
2. Enter the ending time in the **End** input. Note, that the **End** time needs to be larger than the **Start** time.
All gaze tracking data after this time will be filtered out.



To revert to the default setting, enter "0" in both the Start and End input fields and select a new trial data set in the selection tree.

6.14.5 Results Grid

The **Result Grid** shows the parameters of the statistics and the computed values. You can customize the results grid view settings and export the current view to a statistics data file (see [Export Statistics](#)^[158]).

To operate the results grid in order to customize the view settings proceed as follows:

1. To resize columns drag a column header's separator.
2. To move columns to another position drag and drop a column header.
3. To sort the results grid click on the desired column header. To reverse the sort order, click the same column header again.
4. To remove columns, click on the Settings button to open the settings dialog
5. To resize all rows hover the mouse over the left border of the results grid. If the mouse cursor changes, drag and drop to indicate the new height.

The results grid view settings are applied temporary for the currently displayed results. The results grid reverts to the former settings, if new results are computed. New results are computed if you change the [Selection Tree](#)^[153] or when you change the [Time Interval](#)^[157] settings. To make the results grid settings permanent, proceed as described under [Statistics Template](#)^[155].

6.14.6 Export Statistics


You can export the current display of the [Results Grid](#)^[158] to an ASCII data file.

Copy to Clipboard

Click on **Copy to Clipboard** button to copy the current shown statistic

into the clipboard for further use in other programs, e.g. MS-Excel.

Export to file



1. Click the **Export...** button available at the top of the [Event Statistics](#) ¹⁵² data view.

The **Export Statistics** dialog opens. The dialog shows a preview of the ASCII data to be exported.

2. Change the exported number precision in the **Decimal places** input.
3. Change the data separator character in the **Decimal Separator** drop-down list. While most applications will import ASCII data separated by the tab character, some applications may require another separator character.
4. If the first two columns of the exported statistics are "Trial" and "Subject" then a checkbox option called **Use SPSS case format** appears in the File Format area. Checking this option changes the output format so that it contains a single line of data instead of having each trial data on its own line. This is useful for certain analysis done

outside the program.

5. Click the **Export** button. Select the storage location and enter a file name in the subsequent **Save as...** dialog.



The first line of the exported data file lists the column header names. If you import the ASCII file to another application, these names are then available for identifying the columns.

6.14.7 Event Statistics - Definitions and Examples

The following tables list details about the default statistic templates that are shipped with the BeGaze 2.4.

Default Statistic Templates

[Fixation Details](#) ¹⁶²

One row per fixation, process all fixations from all selected trials

[Saccade Details](#) ¹⁶³

One row per saccade, process all saccades from all selected trials

[Blink Details](#) ¹⁶⁴

One row per blink, process all blinks from all selected trials

[Trigger Line Details](#) ¹⁶⁵

One row per trigger event, taken from IDF file

[Event Detailed Statistics](#) ¹⁶⁶

One row per trial, process all selected trials

[Event Summary Statistics](#) ¹⁶⁸

One row for all trials, compute values over all selected trials

[AOI Fixations](#) ¹⁷¹

One row for each fixation that hits one AOI, process all selected trials, only on selected AOIs

[AOI Detailed Statistics](#) ¹⁷²

One row for each AOI-trial combination, process all selected trials, only on selected AOIs

[AOI Summary Statistics](#) ¹⁷⁵

One row per AOI, compute values over all selected trials associated with one AOI

[AOI Transition Matrix](#) ¹⁸¹

One row per AOI, number of consecutive fixation transitions inside and between selected AOIs for all selected trials

[User Event Statistics](#) ¹⁸²

One row per recorded user event for all selected trials.

[Noldus Observer Export](#) ¹⁸³

One state change per row

[Questionnaire Statistics](#) ¹⁸³

One questionnaire per line, taken from Experiment Center questionnaires

[Subject Statistics](#) ¹⁸³

One row per subject, shows subject calibration information

[Stimulus Statistics](#) ¹⁸⁴

One stimulus per row, shows stimulus information

Notes and Definitions

All processing is constrained to the selected time interval. All fields without a comment represent information extracted directly from the event properties, with average/max/min as the only statistic measurement done when indicated.

The following table comments terms used in the subsequent table texts.

Name	Definition
Dwell time	Dwell time starts at the moment the AOI is fixated and ends at the moment the last fixation on the AOI ends = sum of durations from all fixations and saccades that hit the AOI

Glance Duration	Saccade duration for entering the object + sum of all fixation durations and saccade durations before the eyes begin to leave the AOI = dwell time + duration of saccade entering AOI
Diversion Duration	Sum of saccade durations for entering and leaving the object + sum of all fixation durations and saccade durations before the eyes begin to leave the AOI = glance duration + duration of saccade leaving AOI
Duration Before	Time until AOI is found = start time of first fixation to enter the AOI
Glances	Number of glances to a target (saccades coming from outside) within a certain period (increment the counter each time a fixation hits the AOI, if not hit before)
Saccade latency	Duration between consecutive saccades = average of the time difference between the end of a saccade and the start of the consecutive one

The following color codes denote the parameter origin:

- parameters
- event properties
- computed values

Fixation Details

This template shows one row per fixation, process all fixations from all selected trials.

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code

Stimulus		Stimulus Name
Start Time	[ms]	Trial Start Time, normally zero
End Time	[ms]	Trial End Time
Fixation Start	[ms]	Beginning of a fixation.
Fixation Duration	[ms]	Duration of a fixation.
Fixation End	[ms]	End of a fixation.
Position XY		Geographical position of a fixation .
Average pupil size	[px]	Average size of a pupil.
Dispersion	[px]	Dispersion of a fixation.
Eye L/R		Which eye fixated
Number		Number of the fixation.

Saccade Details

This template shows one row per saccade, process all saccades from all selected trials.

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Start Time	[ms]	Trial Start Time, normally zero
End Time	[ms]	Trial End Time
Saccade Start	[ms]	Beginning of a saccade.
Saccade Duration	[ms]	Duration of a saccade.
Fixation End	[ms]	End of a saccade.

Parameter	Dimension unit	Description
Start Position XY		Geographical position where the saccade begins.
End Position XY		Geographical position where the saccade ends.
Amplitude	[°]	Max. oscillation from the rest position of a saccade.
Acceleration average	[°/s ²]	Average acceleration of a saccade in x.
Acceleration peak	[°/s ²]	Peak value of acceleration of gaze during a saccade.
Deceleration peak	[°/s ²]	Peak value of deceleration of gaze during a saccade.
Velocity average	[°/s]	Average velocity of gaze during a saccade.
Velocity peak	[°/s]	Peak value of velocity of gaze during a saccade.
Peak velocity at	[%]	Position of the peak velocity within the saccade.
Eye L/R		Which eye does a saccade
Number		Number of the saccade.

Blink Details

This template shows one row per blink, process all blinks from all selected trials.

Parameter	Dimension unit	Description
Trial		Trial Number

Subject		Subject Code
Stimulus		Stimulus Name
Start Time	[ms]	Trial Start Time, normally zero
End Time	[ms]	Trial End Time
Blink Start	[ms]	Beginning of a blink.
Blink Duration	[ms]	Duration of a blink.
Blink End	[ms]	End of a blink.
Eye L/R		Which eye blinked
Number		Number of the blinks.

Trigger Line Details

Parameter	Dimension unit	Description
Trigger Line start	[ms]	Start time of the trigger event.
Trigger Line duration	[ms]	Duration of the trigger event.
Trigger Line end	[ms]	End time of the trigger event.
Number		Trigger event count.
Port Status		Hardware port ID from where the event was triggered.

Event Detailed Statistics

This template shows one row per trial, process all selected trials.

Parameter	Dimension unit	Description
Trial		Trial Number

Parameter	Dimension unit	Description
Subject		Subject Code
Stimulus		Stimulus Name
Start Time	[ms]	Trial Start Time, normally zero
End Time	[ms]	Trial End Time
Blink count		Number of blinks in the trial.
Blink frequency	[count/s]	Number of blinks per second.
Blink duration total	[ms]	Sum of duration of all blinks.
Blink duration average	[ms]	Sum of duration of all blinks divided by number of blinks in the trial.
Blink duration maximum	[ms]	Longest blink duration.
Blink duration minimum	[ms]	Shortest blink duration.
Fixation count		Number of fixations in the trial.
Fixation frequency	[count/s]	Number of fixations per second.
Fixation duration total	[ms]	Sum of duration of all fixations.
Fixation duration average	[ms]	Sum of duration of all fixations divided by number of fixations in the trial.
Fixation duration maximum	[ms]	Longest fixation duration.
Fixation duration minimum	[ms]	Shortest fixation duration.
Fixation dispersion total	[px]	Sum of all fixation dispersions on X and Y

Parameter	Dimension unit	Description
Fixation dispersion average	[px]	Sum of all fixation dispersions on X and Y divided by number of fixations in the trial.
Fixation dispersion maximum	[px]	Largest value for the sum of X and Y dispersions of one fixation.
Fixation dispersion minimum	[px]	Smallest value for the sum of X and Y dispersions of one fixation.
Saccade count		Number of saccades in the trial.
Saccade frequency	[count/s]	Number of saccade per second.
Saccade duration total	[ms]	Sum of duration of all saccades..
Saccade duration average	[ms]	Sum of duration of all saccades divided by number of saccades in the trial.
Saccade duration maximum	[ms]	Longest saccade duration.
Saccade duration minimum	[ms]	Shortest saccade duration.
Saccade amplitude total	[°]	Sum of all saccades amplitude.
Saccade amplitude average	[°]	Sum of all saccades amplitude divided by number of saccades in the trial.
Saccade amplitude maximum	[°]	Max. saccade amplitude

Parameter	Dimension unit	Description
Saccade amplitude minimum	[°]	Min. saccade amplitude
Saccade velocity total	[°/s]	Sum of all saccades velocities.
Saccade velocity average	[°/s]	Sum of all saccades velocities divided by number of saccades in the trial.
Saccade velocity maximum	[°/s]	Max. value of the saccade velocity.
Saccade velocity minimum	[°/s]	Min. value of the saccade velocity.
Saccade latency average	[°/s]	saccade latency = time between the end of a saccade and the start of the next saccade. Saccade latency average = total saccade latency for all saccades / saccade count

Event Summary Statistics

This template shows one row for all trials, compute values over all selected trials.

Parameter	Dimension unit	Description
Start Time	[ms]	Trial Start Time, normally zero
End Time	[ms]	Trial End Time
Blink count		Number of blinks of all selected trials.

Parameter	Dimension unit	Description
Blink frequency	[count/s]	Number of blinks of all selected trials per second divided by the number of selected trials.
Blink duration total	[ms]	Sum of duration of all blinks of all selected trials.
Blink duration average	[ms]	Sum of duration of all blinks of all selected trials divided by the number of selected trials.
Blink duration maximum	[ms]	Longest blink duration of all selected trials.
Blink duration minimum	[ms]	Shortest blink duration of all selected trials.
Fixation count		Number of fixations of all selected trials.
Fixation frequency	[count/s]	Number of fixations of all selected trials per second divided by the number of selected trials.
Fixation duration total	[ms]	Sum of duration of all fixations of all selected trials.
Fixation duration average	[ms]	Sum of duration of all fixations of all selected trials divided by the number of selected trials.
Fixation duration maximum	[ms]	Longest fixation duration of all selected trials.
Fixation duration minimum	[ms]	Shortest fixation duration of all selected trials.
Fixation dispersion total	[px]	Sum of all fixation dispersions on X and Y of all selected trials.
Fixation dispersion average	[px]	Sum of dispersion of all fixations of all selected trials divided by the number of selected trials.

Parameter	Dimension unit	Description
Fixation dispersion maximum	[px]	Largest value for the sum of X and Y dispersions of fixation of all selected trials.
Fixation dispersion minimum	[px]	Smallest value for the sum of X and Y dispersions of fixation of all selected trials.
Saccade count		Number of saccades of all selected trials.
Saccade frequency	[count/s]	Number of saccades per second of all selected trials divided by the number of selected trials.
Saccade duration total	[ms]	Sum of all saccade duration of all selected trials.
Saccade duration average	[ms]	Sum of all saccade duration of all selected trials divided by the number of selected trials.
Saccade duration maximum	[ms]	Longest saccade duration of all selected trials.
Saccade duration minimum	[ms]	Shortest saccade duration of all selected trials.
Saccade amplitude total	[°]	Sum of all saccades amplitude of all selected trials.
Saccade amplitude average	[°]	Sum of all saccades amplitude of all selected trials divided by the number of saccades in the trial.
Saccade amplitude maximum	[°]	Max. saccade amplitude of all selected trials.

Parameter	Dimension unit	Description
Saccade amplitude minimum	[°]	Min. saccade amplitude of all selected trials.
Saccade velocity total	[°/s]	Sum of all saccades velocities of all selected trials.
Saccade velocity average	[°/s]	Sum of all saccades velocities of all selected trials divided by the number of saccades in the trial.
Saccade velocity maximum	[°/s]	Max. value of the saccade velocity of all selected trials.
Saccade velocity minimum	[°/s]	Min. value of the saccade velocity of all selected trials.
Saccade latency average	[°/s]	saccade latency = time between the end of a saccade and the start of the next saccade. Saccade latency average = total saccade latency for all saccades / saccade count

AOI Fixations

This template shows one row for each fixation that hits one AOI, process all selected trials, only on selected AOIs.

Parameter	Dimension unit	Description
Subject		Subject Code
Area of Interest		AOI Name
AOI Order		AOI Order Number
Fixation Start	[ms]	Beginning of a fixation in an AOI.

Parameter	Dimension unit	Description
Fixation Duration	[ms]	Duration of a fixation in an AOI.
Fixation End	[ms]	End of a fixation in an AOI.
Position XY		Geographical position of a fixation inside an AOI.
Average pupil size	[px]	Average size of a pupil inside an AOI.
Dispersion	[px]	Dispersion of a fixation inside an AOI.
Eye L/R		Which eye fixated inside an AOI.
Number		Number of the fixation.

AOI Detailed Statistics

This template shows one row for each AOI-trial combination, process all selected trials, only on selected AOIs.

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Area of Interest		AOI Name
AOI Order		AOI Order Number
AOI Size	[px]	Size of AOI in pixel - the part overlaying the stimulus is taking into consideration, parts outside of the stimuli area are ignored
AOI Coverage	[%]	AOI size in comparison to Stimulus size
Start Time	[ms]	Trial Start Time, normally zero

Parameter	Dimension unit	Description
End Time	[ms]	Trial End Time
Duration before	[ms]	Duration from start of the trial to the first hit of the AOI.
Sequence		Order of gaze hits into the AOIs based on Entry time (Duration before), lowest entry time = first in sequence.
Net dwell time	[ms]	Sum of sample durations for all gaze data samples that hit the AOI
Dwell time	[ms]	Starts at the moment the AOI is fixated and ends at the moment the last fixation on the AOI ends = sum of durations from all fixations and saccades that hit the AOI.
Glance duration	[ms]	Saccade duration for entering the object + sum of all fixation durations and saccade durations before the eyes begin to leave the AOI = dwell time + duration of saccade entering AOI.
Diversion duration	[ms]	Sum of saccade durations for entering and leaving the object + sum of all fixation durations and saccade durations before the eyes begin to leave the AOI = glance duration + duration of saccade leaving AOI.
First fixation duration	[ms]	Duration of the first fixation to hit the AOI.
Glances count		Number of glances to a target (saccades coming from outside) within a certain period (increment the counter each time a fixation hits the AOI, if not hit before). [both eyes]

Parameter	Dimension unit	Description
Fixation count		Number of fixations inside the AOI.
Fixation count		Number of fixations inside the AOI.
Appearance count		Sum of all appearances of one AOI within one trial: <ul style="list-style-type: none"> – For static AOIs on still images it is always 1 – For dynamic AOIs it is the number of slices where the AOI was visible
Visible time	[ms]	Sum of AOI duration within one trial <ul style="list-style-type: none"> – For static AOI it is end time – start time – For dynamic AOI it is the sum of all durations where the AOI was visible within start and end time
Net dwell time	[%]	Value is calculated with: $\text{Net dwell time (ms)} / (\text{end time} - \text{start time})$
Dwell time	[%]	Value is calculated with: $\text{Dwell time (ms)} / (\text{end time} - \text{start time})$
Fixation time (ms)	[ms]	Adds the fixations times
Fixation time (%)	[%]	Value is calculated with: $\text{Fixation time (ms)} / (\text{end time} - \text{start time})$
Time to first mouse click	[ms]	Time of first mouse click into the AOI, similar to “Duration before” for gaze data.

The Duration before cell contains "-" if the corresponding AOI is not hit by any fixation during the selected period of time.

AOI Summary Statistics

This template shows one row per AOI, compute values over all selected trials associated with one AOI.

Parameter	Dimension unit	Description
Stimulus		Stimulus Name
Area of Interest		AOI Name
AOI Group		AOI Group Name
AOI Scope		Scope of AOI - local or global
AOI Order		AOI Order Number
AOI Size	[px]	Size of AOI in pixel - the part overlaying the stimulus is taking into consideration, parts outside of the stimuli area are ignored
AOI Coverage	[%]	AOI size in comparison to Stimulus size
Start Time	[ms]	Trial Start Time, normally zero
End Time	[ms]	Trial End Time
Duration before total	[ms]	Sum of duration before of all subjects.
Duration before average	[ms]	Sum of duration before of all subjects divided by number of the subjects.
Duration before maximum	[ms]	Max. duration before of all subjects.
Duration before minimum	[ms]	Min. duration before of all subjects.
Sequence		The order in which the AOIs were fixated.

Parameter	Dimension unit	Description
Stimulus		Stimulus Name
Area of Interest		AOI Name
AOI Group		AOI Group Name
AOI Scope		Scope of AOI - local or global
AOI Order		AOI Order Number
AOI Size	[px]	Size of AOI in pixel - the part overlaying the stimulus is taking into consideration, parts outside of the stimuli area are ignored
Net dwell time total	[ms] / [%]	Sum of net dwell time of all subjects.
Net dwell time average	[ms] / [%]	Sum of net dwell time of all subjects divided by number of the subjects.
Net dwell time maximum	[ms] / [%]	Max. net dwell time of all subjects.
Net dwell time minimum	[ms] / [%]	Min. net dwell time of all subjects.
Dwell time total	[ms] / [%]	Sum of dwell time of all subjects.
Dwell time average	[ms] / [%]	Sum of dwell time of all subjects divided by number of the subjects.
Dwell time maximum	[ms] / [%]	Max. dwell time of all subjects.
Dwell time minimum	[ms] / [%]	Min. dwell time of all subjects.
Glance duration total	[ms]	Sum of glance duration of all subjects.
Glance duration average	[ms]	Sum of glance duration of all subjects divided by number of the subjects.

Parameter	Dimension unit	Description
Stimulus		Stimulus Name
Area of Interest		AOI Name
AOI Group		AOI Group Name
AOI Scope		Scope of AOI - local or global
AOI Order		AOI Order Number
AOI Size	[px]	Size of AOI in pixel - the part overlaying the stimulus is taking into consideration, parts outside of the stimuli area are ignored
Glance duration maximum	[ms]	Max. glance duration of all subjects.
Glance duration minimum	[ms]	Min. glance duration of all subjects.
Diversion duration total	[ms]	Sum of diversion duration of all subjects.
Diversion duration average	[ms]	Sum of diversion duration of all subjects divided by number of the subjects.
Diversion duration maximum	[ms]	Max. diversion duration of all subjects.
Diversion duration minimum	[ms]	Min. diversion duration of all subjects.
First fixation duration total	[ms]	Sum of first fixation duration of all subjects.
First fixation duration average	[ms]	Sum of first fixation duration of all subjects by number of the subjects.
First fixation duration maximum	[ms]	Max. first fixation duration of all subjects.

Parameter	Dimension unit	Description
Stimulus		Stimulus Name
Area of Interest		AOI Name
AOI Group		AOI Group Name
AOI Scope		Scope of AOI - local or global
AOI Order		AOI Order Number
AOI Size	[px]	Size of AOI in pixel - the part overlaying the stimulus is taking into consideration, parts outside of the stimuli area are ignored
First fixation duration minimum	[ms]	Min. first fixation duration of all subjects.
Glances count total		Sum of first glances count of all subjects.
Glances count average		Sum of first glances count of all subjects by number of the subjects.
Glances count maximum		Max. first glances count of all subjects.
Glances count minimum		Min. first glances count of all subjects.
Fixation count total		Sum of first fixation count of all subjects.
Fixation count average		Sum of first fixation count of all subjects by number of the subjects.
Fixation count maximum		Max. first fixation count of all subjects.
Fixation count minimum		Min. first fixation count of all subjects.

Parameter	Dimension unit	Description
Stimulus		Stimulus Name
Area of Interest		AOI Name
AOI Group		AOI Group Name
AOI Scope		Scope of AOI - local or global
AOI Order		AOI Order Number
AOI Size	[px]	Size of AOI in pixel - the part overlaying the stimulus is taking into consideration, parts outside of the stimuli area are ignored
Appearance count total		Sum of all appearances of one AOI within one trial of all subjects.
Appearance count average		Sum of all appearances of one AOI within one trial of all subjects by number of the subjects.
Appearance count maximum		Max. sum of all appearances of one AOI within one trial of all subjects.
Appearance count minimum		Min. sum of all appearances of one AOI within one trial of all subjects.
Visible time total	[ms] / [%]	Sum of AOI duration within one trial of all subjects.
Visible time average	[ms] / [%]	Sum of AOI duration within one trial of all subjects by number of the subjects.
Visible time maximum	[ms] / [%]	Max. sum of AOI duration within one trial of all subjects.
Visible time minimum	[ms] / [%]	Min. sum of AOI duration within one trial of all subjects.
Fixation time total	[ms] / [%]	Added fixations times of all subjects.

Parameter	Dimension unit	Description
Stimulus		Stimulus Name
Area of Interest		AOI Name
AOI Group		AOI Group Name
AOI Scope		Scope of AOI - local or global
AOI Order		AOI Order Number
AOI Size	[px]	Size of AOI in pixel - the part overlaying the stimulus is taking into consideration, parts outside of the stimuli area are ignored
Fixation time average	[ms] / [%]	Added fixations times of all subjects by number of the subjects.
Fixation time maximum	[ms] / [%]	Max. added fixations times of all subjects.
Fixation time minimum	[ms] / [%]	Min. added fixations times of all subjects.
Subject Hit Count		Number of subjects that looked into the AOI
Subject Hit Count	[%]	Number of subjects that looking into the AOI in comparison to all selected subjects
Revisitors count		Number of subjects that looked into the AOI at least 2 times.
Time to first mouse click total	[ms]	Sum of the times of first mouse click into the AOI of all subjects.
Time to first mouse click average	[ms]	Time of first mouse click into the AOI by number of subjects.

Parameter	Dimension unit	Description
Stimulus		Stimulus Name
Area of Interest		AOI Name
AOI Group		AOI Group Name
AOI Scope		Scope of AOI - local or global
AOI Order		AOI Order Number
AOI Size	[px]	Size of AOI in pixel - the part overlaying the stimulus is taking into consideration, parts outside of the stimuli area are ignored
Time to first mouse click maximum	[ms]	Max. time to first mouse click of all subjects.
Time to first mouse click minimum	[ms]	Min. time to first mouse click of all subjects.

The duration before values are computed only on valid trials which associated with a stimulus (the ones that contain at least one fixation inside the corresponding AOI during the selected period of time). The other values are computed on all selected trials associated with the stimulus.

Transition Matrix (Stacking Order, All)

This template shows one row per AOI, number of consecutive fixation transitions inside and between selected AOIs for all selected trials.

Stacking Order: In case of overlapping AOI the most top AOI is taken into consideration

All: All AOI are taken into consideration, even though when they are overlapping

Parameter	Dimension unit	Description
Stimulus		Stimulus Name
from \ to (count)		Column lists all AOI names
Area of Interest		One column for each AOI, all columns for a matrix
[Matrix cells]		Number of transitions from AOI to AOI

User Event Statistics

This template shows one row per recorded user event for all selected trials.

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Time Trial	[ms]	Time, relative to the start of the trial
Time Run	[ms]	Time, relative to the start of the run
Type		User Action/Experiment Event
Event		Scroll/URL/mouse click/key pressed
Content		Content of the message
Content 2		Extra content, e.g. mouse click position

Noldus Observer Export

Parameter	Dimension unit	Description
Time	[ms]	Time of the event
Type		State start/State stop/Point
AOI Name		Name of the AOI

Questionnaire Statistics

Parameter	Dimension unit	Description
Subject		Subject Code
Question		Question text
Answer		User selected answer

Subject Statistics

The subject statistics is independent of the subject, trial and stimuli filtering/selection and shows the general statistics for the subjects.

Parameter	Dimension unit	Description
Subject		Subject Code
<i>Property 1..n</i>		Subject properties
Calibration Deviation X	[°]	Calibration deviation on X
Calibration Deviation Y	[°]	Calibration deviation on Y

Parameter	Dimension unit	Description
Tracking Quality	[%]	Number of non-zero gaze positions divided by sampling frequency multiplied by run duration, expressed in percent.

Stimulus Statistics

Parameter	Dimension unit	Description
Stimulus		Stimulus Name
Subject		Subject Code
Order		Position of the associated trial inside the run
Duration	[ms]	Duration of the associated trial
Width	[px]	Stimulus width
Height	[px]	Stimulus height

6.14.8 Reading Statistics - Definitions and Examples

The following tables list details about the reading statistic templates that are shipped with the BeGaze 2.4 when the reading package is licensed.

Default Statistic Templates

[Fixation Duration](#) ¹⁸⁶

[Saccadic Amplitude](#) ¹⁸⁸

[AOI Statistics](#) ¹⁸⁹

[Landing Position AOI](#) ^[192]

[Pause Duration](#) ^[192]

[First Pass Regression](#)

[Scanpath](#) ^[193]

[Return Sweep](#) ^[194]

[Inner-AOI Regressions](#) ^[196]

[Between AOI Regressions](#)

^[197]

[AOI Hits per Minute](#) ^[199]

Notes and Definitions

All processing is constrained to the selected time interval. All fields without a comment represent information extracted directly from the event properties, with average/max/min as the only statistic measurement done when indicated.

Reading AOI's are generated for

- Paragraphs
- Words
- Sentences
- Characters



Reading AOIs are automatically generated and cannot be self defined but modified in size and position in the AOI editor.



Please note, that character AOIs are disabled by default. When character AOIs are enabled, please be aware that this creates a huge amount of additional data (several thousands of additional AOIs) and will slow down the calculation process for statistics and other computations. It is strongly recommended to leave the character AOIs disabled until they are really needed.

The following color codes denote the parameter origin:

- parameters
- event properties
- computed values

Fixation Duration

This template shows one row per fixation, process all fixations from all selected trials.

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Fixation Start	[ms]	Beginning of a fixation

Fixation Duration	[ms]	Duration of a fixation <i>Note:</i> A longer fixation duration is often associated with a deeper and more effortful cognitive processing. Just and Carpenter (1980) formulated this relation in the influential Strong Eye-Mind Hypothesis, which claims that there is no appreciable temporal lag between what is fixated and what is processed. In reading research, words that are less frequent, and would therefore require a longer lexical activation process, generally get longer fixation durations (Rayner 1998). More complicated texts give rise to longer average fixation durations, ranging from around 200 ms in light fiction to around 260 ms for physics and biology texts (Rayner and Pollatsek, 1989). More complicated grammatical structures give rise to longer fixation durations (Rayner 1978, 1982). Note that fixation duration is an idiosyncratic measure.
Fixation End	[ms]	End of a fixation
Fixation Position XY		Geographical position of a fixation
Word		Fixated word
Reading AOI number		Fixated AOI number
Reading direction		Reading direction (Left to Right or Right to Left)
Eye		Which eye fixated

Saccadic Amplitude

This template shows one row per saccade, process all saccades from all selected trials.

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Saccade start	[ms]	Beginning of a saccade
Saccade duration	[ms]	Duration of a saccade
Saccade end	[ms]	End of a saccade
Saccade startPosition XY		Geographical position where the saccade begins
Saccade endPosition XY		Geographical position where the saccade ends
Saccade amplitude	[px]	<p>Max. oscillation from the rest position of a saccade</p> <p><i>Note:</i> The same effect on saccadic amplitude (and fixation duration) can be found when subject read texts of varying difficulty (Rayner and Pollatsek 1989). Beginning, poor and dyslectic readers have shorter saccadic amplitudes. In oral reading, average saccadic amplitude falls to around 6 letters (1:5), while during music reading and typing, saccades are a mere 1 on average. For subjects reading musical scores, Kinsler and Carpenter (1995) found that the mean saccadic amplitude increased as the tempo of the performance increased.</p>

Start word		Fixated word before saccade started
Start reading AOI number		Fixated AOI before saccade started
End word		Fixated word after saccade ended
End reading AOI number		Fixated AOI after saccade ended
Reading direction		Reading direction (Left to Right or Right to Left)
Eye		Which eye does a saccade

AOI Statistics

This template shows one row for each AOI-trial combination, process all selected trials, only on selected AOIs.

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Area of Interest		AOI name
Reading AOI Type		AOI type
Reading AOI number		AOI number
Fixation count		Number of fixations inside an AOI
Progressive fixations		Number of progressive fixations (preceded by progressive saccades)
Regressions into AOI		Number of regressions into an AOI

Regressions out of AOI		<p>Number of regressions out of an AOI</p> <p><i>Note:</i> While regressions inside words are thought to reflect lexical activation processes (understanding the word), regressions between word reflect sentence integration processes (understanding how several words relate), see chapters 4 and 5 in Underwood (1998).</p>
Regressive fixations		<p>Number of regressive fixations (preceded by regressive saccades)</p>
Single fixation duration	[ms]	<p>The fixation duration of the fixation on a word, for AOIs in which only one fixation has been made</p> <p><i>Note:</i> Single fixation duration is one of the measures for studying lexical activation; known as <i>early processes</i>.</p>
First fixation duration	[ms]	<p>The duration of the first fixation in an AOI (if any)</p> <p><i>Note:</i> Generally, Rayner and Pollatsek (1989) argue that very fast cognitive operation (like lexical activation and recognition) can be measured with first fixation duration, while slower cognitive processes affect gaze duration (=dwell time). The word properties that affect first fixation duration include word frequency, morphological complexity, metaphorical status, orthographic properties, the degree of polysemy and other linguistic computations.</p>

First pass duration	[ms]	<p>Sum of fixation durations from the first entry into an AOI until the eye leaves it in any direction</p> <p><i>Note:</i> First pass gaze duration is considered a measure of linguistic processes slower than lexical activation. Rayner (1998), reviewing reading research using the fixation based gaze duration measure, concludes that gaze duration is indicative both of word frequency and of comprehension processes integrating several words. Gaze duration on a word thus contrasts to first fixation duration, the other major reading measure, which is used as an index on word frequency. "Gaze duration" is a reading research term. It is defined exactly as dwell time.</p>
First return to AOI	[ms]	Time of occurrence for the first re-entry into an AOI
Second pass duration	[ms]	<p>Sum of fixation durations from the second entry into an AOI until the eye leaves it in any direction</p> <p><i>Note:</i> Second pass gaze duration on a word is assumed to reflect late effects (word integration processes).</p>
Ratio saccade / next fixation	[%]	Saccade time divided by next fixation time
Ratio saccade / prev fixation	[%]	Saccade time divided by previous fixation time
Is first skip		<p>AOIs (words) that are not fixated during first pass reading (although they may be fixated during later regressions)</p> <p><i>Note:</i> Readers skip over high predictable words more frequently than low predictable words (Rayner & Well 1996).</p>

Is total skip		AOIs (words) that are never fixated
Eye		Which eye fixated inside an AOI

Landing Position AOI

This template shows one row for each AOI-trial combination, process all selected trials, only on selected AOIs.

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Area of Interest		AOI name
Reading AOI Type		AOI type
Reading AOI number		AOI number
Reading AOI landing position	[%]	Quotient between AOI length and fixation position inside the AOI
Eye		Which eye fixated inside an AOI

Pause Duration

This template shows one row for each AOI-trial combination, process all selected trials, only on selected AOIs.

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code

Stimulus		Stimulus Name
Fixation Start	[ms]	Beginning of a fixation
Fixation Duration	[ms]	Duration of a fixation
Fixation End	[ms]	End of a fixation
Fixation Position XY		Geographical position of a fixation
Word		Fixated word
Reading AOI number		AOI number
Fixation pause	[ms]	Fixation duration + the duration of the subsequent saccade
Eye		Which eye fixated

First Pass Regression Scanpath

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Event type		Type of user event
Start	[ms]	First Pass Regression start time
Duration	[ms]	First Pass Regression duration <i>Note: The duration of the regression scanpath is a measure of sentence integration processes.</i>
End	[ms]	First Pass Regression end time
StartPosition XY		Position when first pass regression started

EndPosition XY		Position when first pass regression ended
Start word		Fixated word when first pass regression started
Start reading AOI number		AOI number when first pass regression started
End word		Fixated word when first pass regression ended
End reading AOI number		AOI number when first pass regression ended
Number		Number of events durring first pass regression
Eye		Which eye fixated

Return Sweep

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Saccade return sweep start	[ms]	Return sweep start time
Saccade return sweep duration	[ms]	Return sweep duration
Saccade return sweep end	[ms]	Return sweep end time
Saccade return sweep startPosition XY		Start position for return sweep
Saccade return sweep endPosition XY		End position for return sweep

Saccade correction start	[ms]	Correction saccade start time
Saccade correction duration	[ms]	Correction saccade duration
Saccade correction end	[ms]	Correction saccade end time
Saccade correction startPosition XY		Start position for correction saccade
Saccade correction endPosition XY		End position for correction saccade
Saccade return sweep start word		Fixated word before return sweep
Saccade return sweep start reading AOI number		Fixated AOI number before return sweep
Saccade return sweep end word		Fixated word after return sweep
Saccade return sweep end reading AOI number		Fixated AOI number after return sweep
Saccade correction end word		Fixated word after correction saccade
Saccade correction end reading AOI number		Fixated AOI after correction saccade
Fixation intermediate start	[ms]	Intermediate fixation start time
Fixation intermediate duration	[ms]	Intermediate fixation duration

Fixation intermediate end	[ms]	Intermediate fixation end time
Fixation intermediate Position XY		Position for intermediate fixation
Fixation intermediate word		Fixated word in intermediate fixation
Fixation intermediate reading AOI number		AOI number in intermediate fixation

Inner-AOI Regressions

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Prev. Fixation start	[ms]	Previous fixation start time
Prev. Fixation duration	[ms]	Previous fixation duration
Prev. Fixation end	[ms]	Previous fixation end time
Prev. Fixation Position XY		Previous fixation position
Next Fixation start	[ms]	Next fixation start time
Next Fixation duration	[ms]	Next fixation duration
Next Fixation end	[ms]	Next fixation end time

Next FixationPosition XY		Next fixation position
Regressive Saccade start	[ms]	Intermediate regressive saccade start time
Regressive Saccade duration	[ms]	Intermediate regressive saccade duration
Regressive Saccade end	[ms]	Intermediate regressive saccade end time
Regressive Saccade startPosition XY		Intermediate regressive saccade start position
Regressive Saccade endPosition XY		Intermediate regressive saccade end position
Area of Interest		AOI name
Reading AOI number		AOI number
Eye		Which eye fixated inside an AOI

Between AOI Regressions

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Prev. Fixation start	[ms]	Previous fixation start time
Prev. Fixation duration	[ms]	Previous fixation duration
Prev. Fixation end	[ms]	Previous fixation end time

Prev. FixationPosition XY		Previous fixation position
Next Fixation start	[ms]	Next fixation start time
Next Fixation duration	[ms]	Next fixation duration
Next Fixation end	[ms]	Next fixation end time
Next FixationPosition XY		Next fixation position
Regressive Saccade start	[ms]	Intermediate regressive saccade start time
Regressive Saccade duration	[ms]	Intermediate regressive saccade duration
Regressive Saccade end	[ms]	Intermediate regressive saccade end time
Regressive Saccade startPosition XY		Intermediate regressive saccade start position
Regressive Saccade endPosition XY		Intermediate regressive saccade end position
Area of Interest start		Previous AOI name
Reading AOI number start		Previous AOI number
Area of Interest end		Next AOI name
Reading AOI number end		Next AOI number
Eye		Which eye fixated inside an AOI

AOI Hits per Minute

This template shows one row per selected trials.

Parameter	Dimension unit	Description
Trial		Trial Number
Subject		Subject Code
Stimulus		Stimulus Name
Reading AOI Hits character		Character AOI hits per minute
Reading AOI Hits word		Word AOI hits per minute <i>Note: This is the word-per-minute (WPM) measure, a classical measure for reading speed. In the eye-tracking version, WPM can be made a continuous measure that varies along the text.</i>
Reading AOI Hits sentence		Sentence AOI hits per minute
Reading AOI Hits paragraph		Paragraph AOI hits per minute
Eye		Which eye fixated inside an AOI

6.14.9 Reading Statistics - References

Berséus, P. (2002). Eye movement in prima vista singing and vocal text reading.

Unpublished master's thesis, Lund University.

Clifton, C., Staub, A., & Rayner, K. (2006). Eye movements: A window on mind and brain. In R. P. G. van Gompel, M. H. Fischer, W. S. Murray, & R. L. Hill (Eds.), (chap. Eye movements in reading words and sentences.). Elsevier, Oxford.

Diepen P.M.J. van, M., W., & G. d'Ydewalle. (1998). Eye guidance in reading, driving, and scene perception. In G. Underwood (Ed.), (chap. Functional division of the visual field: moving masks and moving windows.). Oxford: Elsevier.

Holmqvist, K., Holsanova, J., Barthelson, M., & Lundqvist, D. (2003). Reading or scanning? a study of newspaper and net paper reading. In *The mind's eye: Cognitive and applied aspects of eye movement research* (p. 657-670). Elsevier Science Ltd.

Holmqvist, K., & Wartenberg, C. (2005). The role of local design factors for newspaper reading behaviour – an eye-tracking perspective (Tech. Rep. No. 127). Lund University Cognitive Studies.

Holsanova, J., Holmberg, N., & Holmqvist, K. (2008). Reading information graphics: The role of spatial contiguity and dual attentional guidance. *Appl. Cognit. Psychol.*

Holsanova, J., Rahm, H., & Holmqvist, K. (2006). Entry points and reading paths on newspaper spreads: comparing a semiotic analysis with eyetracking measurements. *Visual Communication*, 5 (1), 65-93. Available from <http://vcj.sagepub.com/cgi/content/abstract/5/1/65>

Huey, E. B. (1898). Preliminary experiments in the physiology and psychology of reading. *American Journal of Psychology*, 9 , 575-586.

Hyrskykari, A. (2006). Utilizing eye movements: Overcoming inaccuracy while tracking the focus of attention during reading. *Computers in Human Behavior*, 22 (4), 657 - 671. Available from <http://www.sciencedirect.com/science/article/B6VDC-4J72YW1-2/2/8ba5466cc55abcaaf639bc20d15314c3>

(Attention aware systems - Special issue: Attention aware systems)

Inhoff, A. W., & Radach, R. (1998). Eye guidance in reading, driving and scene perception. In G. Underwood (Ed.), (p. 29-53). New York: Elsevier.

Johansson, R., Holmqvist, K., Mossberg, F., & Lindgren, M. (2009). Music you prefer to study to impairs reading comprehension.

- Johansson, R., Johansson, V., Wengelin, Å., & Holmqvist, K. (2008).** Reading during writing: Four groups of writers. (Working Papers No. 53). Department of Linguistics, Lund University.
- Johansson, R., Johansson, V., Wengelin, Å., & Holmqvist, K. (in press).** Looking at the keyboard or the monitor: Relationship with text production processes. Reading and Writing.
- Kollmorgen, S., & Holmqvist, K. (2009).** Automatically detecting reading in eye tracking data (LUCS No. 144). Department of Cognitive Science.
- Krischer, C., & Zangemeister, W. (2007).** Scanpaths in reading and picture viewing: Computer-assisted optimization of display conditions. *Computers in Biology and Medicine*, 37 (7), 947 - 956. Available from <http://www.sciencedirect.com/science/article/B6T5N-4NB2SN9-2/2/73537cff18bfb34e14d8e45538be90a2>
- (Vision and Movement in Man and Machines)
- Liversedge, S., Paterson, K., & M.J., P. (1998).** Eye guidance in reading, driving, and scene perception. In G. Underwood (Ed.), (chap. Eye movements and measures of reading time.). Oxford: Elsevier.
- Liversedge, S., White, S., Findlay, J., & Rayner, K. (2006).** Binocular coordination of eye movements during reading. *Vision research*, 46 (15), 2363–2374.
- McConkie, G. W., Kerr, P. W., Reddix, M. D., & Zola, D. (1988).** Eye movement control during reading: I the location of initial eye fixations on words. *Vision Research*, 28 , 1107-1118.
- McConkie, G. W., Kerr, P. W., Reddix, M. D., Zola, D., & Jacobs, A. M. (1989).** Eye movement control during reading: II frequency of refixating a word. *Perception and Psychophysics*, 46 , 245-253.
- Rayner, K. (1998).** Eye movements in reading and information processing: 20 years of research. *Psychological Bulletin*, 124 (3), 372–422.
- Rayner, K., Li, X., Williams, C., Cave, K., & Well, A. (2007).** Eye movements

during information processing tasks: Individual differences and cultural effects. *Vision Research*, 47 (21), 2714–2726.

Rayner, K., & Pollatsek, A. (1989a). *The psychology of reading.*

Englewood

Cliffs, NJ: Prentice Hall.

Rayner, K., & Pollatsek, A. (1989b). *The psychology of reading.*

Renshaw, J. A., Finlay, J. E., Tyfa, D., & Ward, R. D. (2004).

Regressions

re-visited: a new definition for the visual display paradigm. In Chi '04: Chi '04 extended abstracts on human factors in computing systems (pp. 1437–1440). New York, NY, USA: ACM.

Sullivan, B., Jovancevic, J., Hayhoe, M., & Sterns, G. (2005). Use of gaze in natural tasks in stargardt's disease: A preferred retinal region.

International Congress Series, 1282, 608 - 612. Available from

<http://www.sciencedirect.com/science/article/B7581-4H8DPNP-4W/2/e4819ada0573583ad3ed7a0277723bbf>

(Vision 2005 - Proceedings of the International Congress held between 4 and 7 April 2005 in London, UK)

Takeda, Y., Sugai, M., & Yagi, A. (2001). Eye fixation related

potentials in a proof reading task. *International Journal of Psychophysiology*, 40 (3), 181 - 186. Available from

<http://www.sciencedirect.com/science/article/B6T3M-42DP0HD-2/2/586b0c37ba4efcd04248bea842b36f2c>

Underwood, G. (1998). *Eye guidance in reading, driving, and scene perception.*

Elsevier.

6.15 Line Graph

6.15.1 Overview

The **Line Graph** data view displays un-interpreted eye tracking data and gaze events for scientific or informal purposes. Data and events are plotted as lines on a timeline diagram.



Operate the **Line Graph** data view with the following steps:

1. In the **Subjects Selection** ^[56], select a single trial.

The **Line Graph Main Window** ^[206] and **Line Graph Data Table** ^[208] are updated for the selected trial.

While selecting trials, the **Events Selection** ^[61] view and the **Trial Details** ^[60] view shows information about the currently selected trial or event.

2. In the **Line Graph Miniview** ^[209], change to the desired view tab.

The **Miniview** displays the selected stimulus correlated with the gaze

position of the current [Diagram Cursors](#) ^[206].

6.15.2 Events List

The general functionality of this view is described in [Events List](#) ^[61]. The blue data cursor and the red auxiliary cursor will frame the marked event in the [Line Graph Main view](#) ^[206]. The gaze cursor in the [Line Graph Miniview](#) ^[209] will jump to the position at the start time of the event.




Type	Eye	Index	Start (ms)	Duration (ms)
User Message		1		0 0
Fixation	Right	1		2 134
Fixation	Left	1		2 134
Saccade	Left	1	137	19
Saccade	Right	1	137	19
Fixation	Right	2	156	659
Fixation	Left	2	156	483
Fixation	Left	3	676	139
Blink	Left	1	1353	99
Blink	Right	1	1353	119
Blink	Left	2	14051	119
Blink	Right	2	15564	119
Blink	Left	3	15564	119

A highlighted event in the **Events** list. The marked event is framed by two cursors in the Graph Area:



The gaze cursor (blue dot in the full view, a cross in the zoomed view) is at the start time of the event in the Miniview:




[Line Graph Overview](#) ²⁰³

6.15.3 Graph Area

In the **Line Graph** main view, the following gaze data will be visualized over the timeline:

- **Gaze parameters:** The Y-axis at the left displays the gaze position in the stimulus (x- and y-direction) as well as angular velocity and acceleration of the eye.
- **Pupil diameter:** The Y-axis at the right displays the pupil diameter.
- **Time [ms]:** The X-axis at the bottom displays fixations, saccades, blink, and user events.




You can customize the line graph display with the following steps:

1. Right click the line graph display to open a context menu. Select the **Settings** command and change line colors and display in the [Line Graph Settings Dialog](#) ^[210].

2. Click the **Reset Scaling** icon in the top left corner to revert display scaling and positioning.
3. Click the **Legend** button in the top right corner to hide or unhide the legend.

If the legend is displayed, activate or deactivate the options next to the labels. This selects the desired combination of lines to draw.
4. To shift the line graph display scales, drag the left or right Y-axis or drag the bottom X-axis using the finger mouse cursor. To shift the line graph position, hold down the [SPACE] key and drag the display using the hand mouse cursor.
5. To zoom in, simply click into the display. To zoom an arbitrary display portion, click and drag to span a dotted zoom box. If you release the mouse button, the display is zoomed accordingly.
6. To zoom out, hold down the [CTRL] key and click into the display.
7. Click the colored event bar displayed at the bottom of the line graph display. This selects a single event and moves the [Line Graph Diagram Cursors](#)^[208] as well. The respective event is highlighted in the [Events Selection](#)^[61] view, which in turn also updates the [Trial Details](#)^[60] view and the [Line Graph Miniview](#)^[209]. Note, that you can manually drag the diagram cursors using the drag mouse cursor.



Click on the bar to find out more about the event.

8. In the **Export** menu, either select the **Save image to file** ([CTRL] + [S]) or select the **Copy image to clipboard** ([CTRL] + [C]) keyboard command to export the current line graph display to a single image.

6.15.4 Diagram Cursors

If you create a new Line Graph, you will find a blue vertical line in the middle of the Graph Area, the main data cursor. The data cursor is movable, you can drag it to every time in the Graph Area. Simply hover with the mouse over the data cursor until a double slider becomes visible, then click the left mouse button and drag the data cursor to the desired position. Alternatively, you can use the arrow left / arrow right keys on the keyboard.



The data cursor can be used to find out the exact values for the gaze data at a particular time. The gaze data values are displayed in the [Data Table](#) [208] and are immediately updated for the current data cursor position. Furthermore, the gaze point at this time on the stimulus image is displayed in the [Miniview](#) [209] below the Graph Area.

Apart from the data cursor a red auxiliary cursor is visible.

6.15.5 Data Table

In the data table, the data values are displayed numerically for the current [Line Graph Diagram Cursor](#) [208] positions. You will find information about:

- the exact time for the time cursor positions.
- the pupil diameter at this time


- the gaze point in x- and y-direction in [pixels]. (0,0) is the upper left corner of the stimulus image.
- the angular speed of the eye
- the angular acceleration of the eye

If you work with binocular data, the values for both eyes will be displayed.

6.15.6 Miniview

In the **Miniview**, the gaze position at the current [data cursor](#)^[208] is displayed in the stimulus. The Miniview offers two display tabs:

- **Full** tab: shows the complete and resized stimulus.
- **Zoom** tab: shows a magnified area around the gaze position.



6.15.7 Settings

In the **Linegraph Settings** dialog, you select line colors, event colors and customize the line graph chart settings.


1. Right click into the [Line Graph Main Window](#)^[206] to open a context menu. Select the **Settings** command.

The **Linegraph Settings** dialog opens.

2. Switch to one of the available dialog tabs and change settings.
3. Confirm your settings with **OK**.

The following dialog tabs are available.

Left Eye




In this tab you can configure, for left data channel the color and the visibility of:

- gaze data on X
- gaze data on Y
- pupil diameter

- angular velocity
- angular acceleration

Right Eye



In this tab you can configure, for right data channel the color and the visibility of:

- gaze data on X
- gaze data on Y
- pupil diameter
- angular velocity
- angular acceleration


Events



In this tab you can configure the color for the following types of events:

- fixation
- saccade
- blink
- user event

Chart



In this tab you can configure:

- the time range in [s]
- the color of the first cursor line
- the color of the second cursor line
- whether to show values in steps
- whether to show the legend

Event Detection

Chapter



VII

7 Event Detection

7.1 Built-In Event Detector

BeGaze 2.4 has a built-in saccade, fixation and blink detector. A saccade is defined as a rapid change in gaze location, and a fixation is regarded as being bordered by two saccades. A blink can be considered a special case of a fixation, where eye data is not present.

In general, there are two approaches for the built-in detector: Either it can first look for fixations and the other events are derived from them, or it can first look for saccades, followed by the computation of the other events.

Which event the detector searches first, we call *primary event*. If the primary event is *fixation*, the detector uses a *dispersion* based algorithm. If the primary event is saccade, a *velocity* based algorithm is used.

For low speed eye tracking data (< 200 Hz), choosing fixations as primary event achieves the best results, whereas primary looking for saccades is sensible for high speed eye tracking data.

Depending on the sample rate the built-in detector selects the detection method:

sample rate	detection method	primary event	algorithm based on
all data rates	low speed event detection ²¹⁹	fixation	dispersion
200 Hz and above	high speed event detection ²²¹	saccade	velocity

Please note, that none of the algorithms are currently well suited to detect fixations on moving targets in videos where the eyes are following with a smooth pursuit. This issue is currently addressed in ongoing research work.

7.2 Adjust Event Detection

In the **Adjust Event Detection** dialog, you can change the event detection parameters as well as the stimulus geometry for one or more trials.

1. In the [File menu](#)^[246] select the **Adjust Event Detection** command.

The **Adjust Event Detection** dialog opens.

2. In the **Fixation detection parameters** section of the dialog, you can change settings for low speed event detection or for high speed event detection. Which type of settings are available, depends on the gaze tracking device used.
3. In the **Geometry** section of the dialog, you can adapt resolution and dimension of the presented stimuli.
4. Confirm you settings with **OK**.

When creating an experiment, you can adjust these parameters in the [Event Detection](#)^[37] tab of the **Create Experiment** wizard.


Exclude first fixation

The first fixation can be deleted from all datasets in the experiment if required.

Exclude first fixation

Low Speed Event Detection Settings

For [Low Speed Event Detection](#)^[219] the following parameters are displayed and can be changed:



Fixation detection parameters

Min. duration:	<input type="text" value="80"/>	ms
Max. dispersion:	<input type="text" value="100"/>	px


Min. duration: minimum fixation duration in [ms]

Max. dispersion: maximum dispersion value. The unit depends on the [experiment type](#)^[256]:

	Unit
standard data	pixels
data with head tracking	degrees

High Speed Event Detection Settings

For [High Speed Event Detection](#)^[221] the following parameters are displayed and can be changed:



Min. duration: minimum saccade duration in [ms]. If the Auto option is checked, the minimum duration varies and is automatically set dependent on the peak threshold.

Peak velocity threshold: peak velocity threshold in [°/s]

Min. fixation duration: minimum fixation duration in [ms]. All fixations below the threshold are rejected.

Peak velocity window: The single peak value has to lie in this window. Start and end is given in % of the saccade length.

For more information see [Built-In Event Detector](#)^[215].

If you click on **Adjust**, the saccades, fixations and blinks will be

recalculated for all the trials in the experiment, using the displayed detection parameters. The changes are persistent for each trial.

Geometry Settings

Geometry

Use geometry from file

Stimulus screen resolution: 1280x1024

Horizontal: 1280 px Vertical: 1024 px

Phys. stimulus dimensions: Horizontal: 340 mm Vertical: 270 mm

Distance monitor-head: 680 mm

Activate the **Use geometry from file** option to read in the screen resolution and physical stimulus dimension settings from the gaze tracking data file. Activating this option immediately reads in the settings from the gaze tracking data file and disables the respective controls.

Deactivate the **Use geometry from file** option if you want to overwrite these settings manually. Deactivating this option enables the following settings:

Stimulus screen resolution: Enter the horizontal and vertical resolution (in pixels) of the monitor which originally displayed all visual stimuli. A list of typical screen resolutions is offered in a drop-down list for selection. To enter a screen resolution not available in the list, select the **user defined** entry or simply enter the desired resolution in the respective text input controls. Note, that all visual stimuli attributed with the **Fit to Screen** option will be recalculated and scaled to this resolution.

Phys. stimulus dimensions: Enter the horizontal and vertical display dimensions in millimeters. Note, that a typical CRT or LCD computer monitor has a display resolution between 72 dpi and 120 dpi with the same horizontal and vertical dpi resolution. Example: a 96 dpi LCD

monitor displaying 1280 horizontal pixels should have a width of 338 mm (1280 px / 96 dpi * 25,4 mm per inch). Note also that other displays such as a video beamer emitting camcorder material typically use a different dpi resolution for horizontal and vertical display.

Distance monitor-head: If you change the **Phys. stimulus dimensions** settings, you need to adapt the approximate distance between the displaying monitor and the subjects head accordingly. Note that during calibration the individual relation between the gaze tracking system and the subject is established. The calibration outcome is not changed nor invalidated with this setting.



Overwriting and changing the geometry settings requires BeGaze 2.4 to re-calculate the gaze tracking data in order to adapt to the new settings. For longer experiments, the recalculation may require some time with the progress indicated by a status dialog.

7.3 Low Speed Event Detection

In the Low Speed Event Detection method, Fixation is selected as primary event. The [Built-In Detector](#)^[215] will first search for fixation events, using a dispersion based algorithm, after which saccade events are computed and derived from the primary fixation events.

Blink Detection

A blink can be regarded as a special case of a fixation, where the horizontal and vertical gaze position equals 0. If this is the case, we create a blink event. However, the duration of the blink event is expanded in order to include the transition period between valid gaze data and the blink.

Fixation Detection

The Minimum Fixation Duration defines the minimum time window in which the gaze data is analyzed. Fixations smaller than the time window will not be caught.

The algorithm identifies fixations as groups of consecutive points within a particular dispersion, or maximum separation. It uses a moving window that spans consecutive data points checking for potential fixations. The moving window begins at the start of the protocol and initially spans a minimum number of points, determined by the given Minimum Fixation Duration and sampling frequency.

The algorithm then checks the dispersion of the points in the window by summing the differences between the points' maximum and minimum x and y values; in other words, dispersion $D = [\max(x) - \min(x)] + [\max(y) - \min(y)]$. If the dispersion is above the Maximum Dispersion Value, the window does not represent a fixation, and the window moves one point to the right. If the dispersion is below the Maximum Dispersion Value, the window represents a fixation. In this case, the window is expanded to the right until the window's dispersion is above threshold. The final window is registered as a fixation at the centroid of the window points with the given onset time and duration.

Saccade Detection

At the end a saccade event is created between the newly and the previously created blink or fixation.

Parameters

The parameters can be changed in the [Adjust Event Detection](#)^[216] dialog.



Min. duration: minimum fixation duration in [ms]

Max. dispersion: maximum dispersion value. The unit depends on the [experiment type](#)^[258]:

	Unit
standard data	pixels

	Unit
data with head tracking	degrees

Further Reading:

Dario D. Salvucci & Joseph H. Goldberg:

[Identifying Fixations and Saccades in Eye-Tracking Protocols](#)

In: Proceedings of the Eye Tracking Research and Applications Symposium (pp. 71-78). New York, 2000

7.4 High Speed Event Detection

In the High Speed Event Detection method, Saccade is selected as primary event. The [Built-In Detector](#)^[215] will first search for saccade events, using a velocity based algorithm. Blinks and fixations are computed and derived from the primary saccade events.

Saccade Detection

From the gaze stream all velocities are calculated. From all velocities the peaks are detected. A peak is defined as the peak value of velocities above the Peak Threshold [°/s]. The peak could indicate a saccade, but as we are not sure, yet, we call it saccade-like event. To detect the start of the saccade-like event, we search for the first velocity to the left which is lower than the fixation velocity threshold. To detect the end of the saccade-like event, we search for the first velocity to the right which is lower than the fixation velocity threshold. The fixation velocity threshold is an internal value calculated from the first peak less velocities of the velocity stream. We assume the saccade-like event a real saccade, if

- the distance between start and end exceeds the Minimum Saccade Duration [ms] and
- the single peak value lies in the range of 20% to 80% of the distance between start and end

Blink Detection


However, the saccade we have found could still be an artifact as a result of a start or end of a blink. If so, we discard the saccade event and assign the artificial saccade to a blink. To determine, if this is the case we evaluate the pupil diameter during the saccade period. If the speed of the pupil diameter change exceeds an internal threshold value, the saccade is assumed artificial and part of the blink.

Fixation Detection

Finally, we create a fixation event between the newly and the previously created blink or saccade.

Parameters

The parameters can be changed in the [Adjust Event Detection](#) ²¹⁶ dialog.



Min. duration: minimum saccade duration in [ms]. If the Auto option is clicked, the minimum duration varies and is automatically set dependent on the peak threshold.

Peak threshold: peak velocity threshold in [°/s]

Min. fixation duration: minimum fixation duration in [ms]. All fixations below the threshold are rejected. The default value is 50 ms.


Peak Velocity Window

The velocity curve must resemble a certain pattern to be regarded as the velocity of a saccade. In a typical saccade the velocity of the eye movement increases, reaches a peak and decreases. At first, the detector assumes this kind of movement to be a saccade. The time between start and end of movement is called saccade length. Then the detector searches, if the velocity peak lies within a certain time window inside of the saccade. If the peak lies outside, the assumed saccade is discarded. The start and end of the time window is given in % of the saccade length.

Default values:

Start: 20% of saccade length

End: 80% of saccade length



Export and Conversions

Chapter



VIII

8 Export and Conversions

8.1 Overview

BeGaze 2.4 allows [events export](#)^[225] and [raw data export](#)^[232]. Furthermore, you can record the replay of the scan path, attention map or key performance indicators to an AVI file (see [Video Export](#)^[239]).

8.2 Export Events

8.2.1 Export Events

In case you want to perform further evaluation with third party software, it is possible to export the events to a custom delimited table in ASCII text format.

If you click on the toolbar item or go to the Export menu and select Export event data to file, a window will be displayed, containing the following tabs:

- General
- Preview

Trial selection

Select the Trials from the Experiment, whose Events should be exported. For each Trial a separated file will be created.

Events to Export

Select from the available events the ones that should be contained in [export file](#)^[228].

Events to Export

Events

Fixations User Events
 Saccades Trigger Line
 Blinks

Channel

Left Eye
 Right Eye
 Binocular

Misc

Include object hit information

File Format

Configure the format of the [export file](#) ²²⁸.

File Format

Write header


Decimal places:

Decimal separator:

Separator:

Export location:

Write Header

Select whether the [Header](#)  will be written in the file.

Decimal Places


Configure the format of the numerical values.

Separator

The separator between values can be one of the following:


- Tab
- Space
- Comma
- Semicolon

Export Location

Click on  to browse for the folder or to create a new folder. BeGaze 2.4 will create the file names automatically.

Preview

You can preview the exact format of the export file. Note: in trial section, only a few data lines are shown.



The Export file may include information about:

- the start and the end time of the fixation, the fixation duration.
- the gaze coordinates at the beginning of the fixation.
- the dispersion during the fixation in [pixels]
- the object hit during the fixation
- the amplitude of a saccade
- the maximum speed and acceleration of the saccade and the time when these maxima occurred

In case the [experiment](#) ^[258] contains head tracking data, additionally will be exported:

- the image name connected to a plane during a fixation on this plane
- the plane number during a fixation on it

8.2.2 Export File Format

8.2.2.1 Export File Format

The BeGaze 2.4 export file starts with a [short header](#) ^[229] section, followed by the [trial section](#) ^[229].

The file can be opened and read with any text editor, but as the entries are tab limited, it will be best read with a spreadsheet program like Microsoft Excel or similar.

8.2.2.2 Header

The header consists of the following few lines:

Converted from:	Complete path of the IDF file.
Date:	Date and time of the export.
Version:	Version, with which the export file is created.
Sample Rate:	Sample rate of the recording.
Subject:	Subject as written to IDF file or modified in experiment creation.
Description:	Description of Run as written to IDF file or modified in experiment creation.

8.2.2.3 Trial Section

The table header description is followed by the list of events.

Every event type has a different table header.

Event Export Fixations

The table header for fixations applies for all lines starting with the word Fixation.

The table headers mean the following:

Event Type:	fixation, L for left or R for right
Trial:	number of current trial
Number:	index of current fixation

Start:	start time in microseconds, relative to start time of beginning of the current trial
End:	end time in microseconds, relative to start time of beginning of the current trial
Duration:	duration of fixation in microseconds
Location X:	horizontal location of fixation in pixel on calibration area
Location Y:	vertical location of fixation in pixel on calibration area
Dispersion X:	horizontal dispersion of fixation in pixel
Dispersion Y:	vertical dispersion of fixation in pixel
Object hit:	name of area of interest (AOI) that is hit by current fixation. The field could be '-', if no AOI is hit.

Event Export Saccades

The table header for saccades applies for all lines starting with the word Saccade.

The table headers mean the following:

Event Type:	saccade, L for left or R for right
Trial:	number of current trial
Number:	index of current saccade
Start:	start time in microseconds, relative to start time of beginning of the current trial
End:	end time in microseconds, relative to start time of beginning of the current trial
Duration:	duration of saccade in microseconds
Start Pos X:	horizontal start position of saccade in pixel on calibration area

Start Pos Y:	vertical start position of saccade in pixel on calibration area
End Pos X:	horizontal end position of saccade in pixel on calibration area
End Pos Y:	vertical end position of saccade in pixel on calibration area
Amplitude:	length of saccade in degrees
Peak Speed:	maximum speed of eye movement during current saccade
Peak Speed At:	location of speed maximum in parts of complete amplitude (a value of 0.416 means peak speed reached at 41.6% of amplitude)
Average Speed:	average velocity of current saccade in degrees per second
Peak Accel.	maximum acceleration of current saccade in deg/s ²
Peak Decel.:	maximum deceleration of current saccade in deg/s ²
Average Accel.	average acceleration of current saccade in deg/s ²

Event Export Blinks

The table header for blinks applies for all lines starting with the word Blink.

The table headers mean the following:

Event Type:	blink, L for left or R for right
Trial:	number of current trial
Number:	index of current blink
Start:	start time in microseconds, relative to start time of beginning of the current trial

End:	end time in microseconds, relative to start time of beginning of the current trial
Duration:	duration of blink in microseconds

Event Export User Messages

The table header for user messages applies for all lines starting with the word Blink.

The table headers mean the following:

Event Type:	user message
Trial:	number of current trial
Number:	index of current user message
Start:	start time in microseconds, relative to start time of beginning of the current trial
Description:	content of the message



Note, that the origin of the calibration area is always in the upper left corner.

8.3 Export Raw Data

8.3.1 Export Raw Data

In case you want to perform further evaluation with third party software, it is possible to export the raw data to a custom delimited table in ASCII text format.

If you click on the toolbar item or go to the Export menu and select Export raw data to file, a window will be displayed, containing the following tabs:

- General
- Preview

Trial selection

Select the Trials from the Experiment, whose Raw Data should be exported. For each Trial a separated file will be created.

Fields to Export

Select from the available events the ones that should be contained in [export file](#) ^[236].

Fields to Export

Raw Data

- Pupil position
- Corneal reflex (CR)
- Pupil diameter
- Head Position
- Head Rotation

Channel

- Left Eye
- Right Eye
- Binocular

Points of Regard (POR)


- Gaze position
- Plane number hit
- Object hit
- Eye position
- Gaze vector
- Quality values
- Quality values

Misc. Data

- Messages
- Frame counter
- Trigger
- Hexadecimal
- Decimal
- Event info
- Stimulus

File Format

Configure the format of the [export file](#) ^[236].



Write Header

Select whether the [Header](#)^[236] will be written in the file.

Decimal Places


Configure the format of the numerical values.

Separator

The separator between values can be one of the following:

- Tab
- Space
- Comma
- Semicolon

Export Location

Click on  to browse for the folder or to create a new folder. BeGaze 2.4 will create the file names automatically.

Preview

You can preview the exact format of the export file. Note: in trial section, only a few data lines are shown.

Preview

```
## [BeGaze]
## Converted from: W:\Research Systems\BeGaze2\test use cases\Demo Cases\Ads light\es-cv5-1.idf
## Date: 04.10.2008 13:28:39
## Version: BeGaze 2.1.30
## Sample Rate: 50
## [Run]
## Subject: cv5
## Description: Run1
## [Calibration]
## Calibration Type: 9-point
## Calibration Area: 1280 1024
## [Geometry]
## Stimulus Dimension [mm]: 376 301
## Head Distance [mm]: 700
## [Hardware Setup]
## [Presentation]
## Number of Samples: 250
## Reversed: none
## Format: MSG
##
Time Type Trial
6961867180 MSG 1 # Message: image11.bmp
6961872225 SMP 1
6961898994 SMP 1
6961919298 SMP 1
6961942882 SMP 1
6961967586 SMP 1
```

8.3.2 Export Raw File Format

8.3.2.1 Export Raw File Format

The BeGaze 2.4 export file starts with a [short header](#)^[236] section, followed by the [trial section](#)^[237].

The file can be opened and read with any text editor, but as the entries are tab limited, it will be best read with a spreadsheet program like Microsoft Excel or similar.

8.3.2.2 Header

The header consists of the following few lines:

Converted from:	Complete path of the IDF file.
Date:	Date and time of the export.
Version:	Version, with which the export file is created.
Sample Rate:	Sample rate of the recording.
Subject:	Subject as written to IDF file or modified in experiment creation.
Description:	Description of Run as written to IDF file or modified in experiment creation.
Calibration Type:	Type of calibration used during recording.
Calibration Area:	Width and height of the calibration area.
Stimulus Dimension:	Width and height of the stimulus.
Head Distance:	Distance between subject and stimulus during recording.
Number of Samples:	Number of samples in the exported trial.
Reversed:	Specifies whether the recorded values were reversed on horizontal and/or vertical axis.
Format:	Format of the exported fields.

8.3.2.3 Trial Section

The table header description is followed by the list of samples and messages.

Raw Data Export Samples

The following fields can be exported for one sample (if available):

Time:	Timestamp of the sample.
Type:	The type is SMP.
Trial:	Number of current trial.
L Raw X [px]:	Horizontal pupil position.
L Raw Y [px]:	Vertical pupil position.
L Dia X [px]:	Horizontal pupil diameter.
L Dia Y [px]:	Vertical pupil diameter.
L CR1 X [px]:	Horizontal corneal reflex position. One or two CRs can be present.
L CR1 Y [px]:	Vertical corneal reflex position.
L POR X [px]:	Horizontal gaze position
L POR Y [px]:	Vertical gaze position
Timing, Latency:	Quality values
L Plane:	Plane number
L Object Hit:	Name of area of interest (AOI) that is hit by current sample.
H POS X [mm]:	Head position on X
H POS Y [mm]:	Head position on Y
H POS Z [mm]:	Head position on Z

H ROT X [°]:	Head rotation on X
H ROT Y [°]:	Head rotation on Y
H ROT Z [°]:	Head rotation on Z
L EPOS X:	Eye position on X
L EPOS Y:	Eye position on Y
L EPOS Z:	Eye position on Z
L GVEC X:	Gaze vector on X
L GVEC Y:	Gaze vector on Y
L GVEC Z:	Gaze vector on Z
Frame:	Frame counter
L Event Info:	Type of event detected for the interval containing this sample (fixation, saccade, blink)
Stimulus:	Stimulus associated with this sample

In case of binocular recordings, data from both channel (named L and R) can be exported.

Raw Data Export Messages

The following fields are exported for one message, along with the actual message:

Time:	Timestamp of the sample.
Type:	The type is MSG
Trial:	Number of current trial



Note, that the origin of the calibration area is always in the upper left corner.

8.4 Video Files


8.4.1 Video Export

You can record the animated **Attention Map** as well as the **Scan Path** or the **Key Performance Indicators** replays to an AVI file.

1. From the **Export** menu, select the **Export Scan Path Video**, **Export Attention Map Video** or **Export KPIs Video** command.

The **Export to File** dialog opens, where you can set the recording options and start the export.

2. Press **Export....**



3. A popup dialog appears allowing you to select the desired video file name and location. Click "Save" to finish.

Dialog Settings

- **Video Compression:** Shows that the “XMP4” video encoder is used. Note, that you need to install this codec from the product installation CD.
- **Frames per second:** This setting applies to a still image stimulus. In case of a video stimulus, the stimulus’ frame rate will be adopted. Select the number of frames per second for the exported AVI video. You can select 10, 25 or 50 frames per second or the eye tracking sampling rate. Higher framerates results in longer export times.
- **Apply Watermark:** Overlay a watermark image over the exported video. The overlay can be Solid or Half Transparent. You can also select a custom image by pressing the button "...".
- **Show Legend:** For plugins that can show a color legend (Heat Map, gridded AOIs) this setting toggles the visibility of such legend in the exported video.

For experiments that contain user videos (user data recorded with a webcam) several other options are available. If no used data exists the **User Video** options are grayed out.

- **User Video Overlay:** If checked the user video is overlaid as a smaller image (picture-in-picture style) inside the animated data visualization.
 - **Opacity:** Selects the opacity level of the user video. Moving the slider to the left fades out the user video more.
 - **Dimensions:** Size of the user video to embed in the main video.
- **User Audio:** If checked the sound from the user video is used as the sound for the exported AVI (if the stimulus is a video with sound then this setting replaces the stimulus sound with the user sound)
- **User Video Location:** The yellow rectangle can be dragged on the gray surface to set the position of the user video relative to the main video in the exported AVI.

8.4.2 Optimizing AVI Videos

The real-time video display and edit functions require appropriate computing resources. While it is necessary to use a modern and powerful PC, it is possible to optimize video data for use with BeGaze 2.4. The video file conversion described below will give a faster response while editing AOIs and working with the video data during analysis.

All video streams are stored as a sequence of single images. To save disk space or transport bandwidth, the following techniques are used:

- The stored image frames are compressed, which normally means that an algorithm is used to encode and decode the single image frames. Most of the image codecs (“Coder/Decoder”) will discard visible information for better compression. There is a tradeoff between file size and visible details.
- If you store images frame after frame, the resulting file size is huge even if the frames are compressed. For this reason, only some frames are stored completely – as “key frames”. All frames following a key frame are generated based on the key frame with additional transformations applied. A high compression video codec will insert key frames only, if it detects major scene changes in the base material. While this is fine for sequential watching, stepping some frames backward requires a lot of calculation. There is also a tradeoff between file size and necessary CPU resources.
- To optimize the user experience for the standard use case “watching the video”, post-processing is applied while reading the video file and displaying it’s contents on the screen. This includes for example to sharpen the video, video scaling or de-interlacing TV material for a non-interlaced computer monitor. There is a tradeoff between screen rendering quality and CPU resources.

BeGaze 2.4 works best with the customized Xvid Solutions MPEG-4 codec (XMP-4) installed during BeGaze 2.4 setup. The post-processing configuration for this codec, which is also applied during setup, is optimized for editing and analyzing purposes. You should convert your video material to this codec and insert more key frames while doing so.



The XMP-4 codec is compatible to standard Xvid and DivX codecs

for playback.

8.4.3 Converting Videos with SMI Video Optimizer

Videos are automatically converted while creating experiment in Experiment Center V2.4.

As an alternative, the SMI Video optimizer can be used to convert videos as well (**Start->SMI->Experiment Suite 360°->Tools**)

The SMI Video Optimizer converts (re-encodes) nearly every kind of video into our recommended XMP4 avi format with optimal codec settings.

The XMP4 codec is automatically installed and configured on the PC during the installation of Experiment Center and BeGaze.

Supported Video Formats


The video optimizer has been successfully tested with a huge variety of video formats and codecs, including DVD (vob), MPEG and Flash (flv) videos.

Nevertheless it depends on the installed and licensed codecs if the selected videos can be converted, which is the responsibility of the user.

The Video Optimizer is using Microsofts Direct X interface to read and convert the selected videos. Please ensure that you have all codecs licensed and installed that your original video needs in order to be read.

Open Video Optimizer

Click on the Video Optimizer entry in the start menu under All Programs -> SMI -> Experiment Suite 360° -> Tools -> Video Optimizer to execute the program.



Adding Videos

1. Click on the buttons **Add file(s)** or **Add folder** to add the videos you want to re-encode.
2. You can also add files by drag&drop of video files from programs like the Windows Explorer into the **Input file** area of the Video Optimizer.

Convert Videos

Press the **Convert Movies** button to start the re-encoding of your videos.

Please note that the new video files are renamed. The re-encoded videos are saved as Originalname + "(optimized).avi"



Original AVI file are not overwritten, please rename the converted video if necessary.

8.4.4 Background Information

The AVI (“Audio Video Interleaved”) container file format is highly suitable for editing purposes. The file format was invented in the 1990’s, with the developing focus on CPU resources with no copy/edit protection nor internet distribution in mind. One of the major drawbacks of this format is the CBR (“Constant Bit Rate”) audio support. It is possible to add VBR (“Variable Bit Rate”) audio material – but this violates the original format specification which may trigger viewer incompatibilities. VBR audio is used most likely for internet video or converted DVD material while self-recorded material usually has CBR audio. If you experience audio dropouts or audio-lag, you can extract the audio file from the AVI file, convert the audio using a CBR codec and re-include the CBR audio to a new AVI file. Another option is to use a special version of VirtualDub called “Nandub” for writing an AVI with VBR audio.

Workspace Reference

Chapter



IX

9 Workspace Reference

9.1 Menu Commands

The following gives an overview of the menu commands:

File

New Experiment...	Starts the Create Experiment wizard ^[27] to create a new experiment
New Experiment from Folder...	Creates an experiment on the basis of a results folder which has been created with the SMI Experiment Center
Open Experiment...	Opens a dialog box to select a saved experiment from the database ^[259]
Close Experiment	Closes the current experiment
Save Experiment	Saves the current experiment to the database ^[259]
Save Experiment As...	Saves the current experiment as a new experiment in the database ^[259]
Define Annotations...	Opens the Define Annotations ^[46] dialog where new annotation types can be defined
Modify Experiment...	Opens the Modify Experiment wizard ^[41] , where all parameters used to create an experiment can be changed
Adjust Event Detection...	Opens the dialog to change and edit the event detection parameters
Delete Experiment from Database...	Opens a dialog to delete a saved experiment from the database ^[259]

Backup Experiment to File...	Opens a dialog to select a saved experiment from the database ^[259] . A backup of the selected experiment will be created in a file.
Restore Experiment from File...	Opens a file selection dialog to select and restore an experiment from file
Print Preview	Opens the print preview
Print...	Opens the printing dialog
Change Data Storage Location...	Opens a folder selection dialog to select another location for the database ^[259]
Reset Plugin Detection	On the next run of BeGaze 2.4, the available data views will be dynamically detected
Recent Experiments	Opens a sub menu with the last opened experiments
Quit	Closes BeGaze 2.4

View

Close Selected View	Closes the selected view
Close All	Closes all opened views
Close All but Selected View	Closes all the views except selected one
Toolbar	Toggles activation/deactivation of the toolbar ^[249]

Analysis

AOI Editor	Opens the AOI Editor ^[75] data view
Gaze Replay	Opens the Gaze Replay ^[100] data view
Bee Swarm	Opens the Bee Swarm ^[102] data view
Scan Path	Opens the Scan Path ^[108] data view

Focus Map	Opens the Focus Map ^[119] data view
Heat Map	Opens the Heat Map ^[125] data view
Key Performance Indicators	Opens the Key Performance Indicators ^[130] data view
Gridded AOIs	Opens the Gridded AOIs ^[137] data view
AOI Sequence Chart	Opens the AOI Sequence ^[146] data view
Binning Chart	Opens the Binning Chart ^[150] data view
Event Statistics	Opens the Event Statistics ^[152] data view
Reading Statistics	Opens the Reading Statistics ^[184] data view
Line Graph	Opens the Line Graph ^[203] data view

Export

Export [...] Video	Exports the currently displayed gaze replay, bee swarm, scan path, focus map, heat map, kpis or gridded aois to a video file. These Menu commands are available only if the corresponding data views are activated.
Save Image to File...	Saves the graph/chart from the currently selected view to an image file. The following file formats are supported: BMP, JPG, PNG.
Copy Image to Clipboard	Copies the graph/chart from the currently selected view to clipboard. Afterwards, it can be pasted into other third party applications.
Export Raw Data to File...	Opens the Raw Data Export dialog, which allows the creation of text files from the raw data of an experiment

Export Event Data to File... Opens the **Event Export** dialog, which allows the creation of text files from the computed event data of an experiment

Help

Help Topics Opens this manual

About BeGaze 2.4 Shows general information about BeGaze 2.4 (see [About Box](#) ²⁵⁶)

9.2 The Toolbar

The toolbar is at the top of the workspace. It gives you short-cuts to important features.



Here is an overview of the buttons and its meanings:

General buttons



Starts the [Create Experiment wizard](#) ²⁷ to create a new experiment for standard data



Opens a dialog to select an existing experiment



Saves the current experiment



Prints the current diagram.



Opens a dialog to remove existing experiment(s)

Diagram selection



[Gaze Replay](#)^[100]: displays a quick gaze data overlay over all the stimulus images in the experiment



[Bee Swarm](#)^[102]: displays raw gaze data overlay over the stimulus image



[Scan Path](#)^[108]: displays gaze data overlay over the stimulus image



[Focus Map](#)^[119]: shows gaze patterns over the stimulus image visualized as a transparent map



[Heat Map](#)^[125]: shows gaze patterns over the stimulus image visualized as a colored map



[Key Performance Indicators](#)^[130]: displays relevant statistical data for each defined AOI over the stimulus image



[Gridded AOIs](#)^[137]: displays relevant statistical data for an automatically defined AOI grid over the stimulus image



[AOI Sequence Chart](#)^[146]: displays AOI hit order over time



[Binning Chart](#)^[150]: gives a statistical overview of AOI hits per binning frame



[Event Statistics](#)^[152]: computes diverse statistics based on events and AOI hits



[Reading Statistics](#)^[184]: computes diverse statistics based on events and AOI hits on text for reading experiments



[Line Graph](#)^[203]: displays x and y directions of gaze data plotted as graphs over time and events displayed in a timeline

Export buttons



Opens a dialog that allows to export raw data to file



Opens a dialog that allows to export events to file

Other commands




Opens the [AOI Editor](#) 

9.3 Hotkeys Overview

Several functions of BeGaze 2.4 can be executed using keyboard commands. The following tables give you an overview.

General keyboard commands

Keys	Description
[CTRL] + [O]	opens a dialog box to select a saved Experiment from the Database 
[CTRL] + [W]	closes the view of the selected data view
[CTRL] + [SHIFT] + [W]	closes all views of opened plug-ins
[CTRL] + [B]	closes all views of opened data views but selected one
[CTRL] + [G]	saves current settings globally

Keys	Description
[CTRL] + [E]	saves current settings for the current experiment
[CTRL] + [C]	copies selected diagram to clipboard, so it can be pasted into other third-party applications
[CTRL] + [S]	saves selected diagram to an image file
[F1]	opens this help file
[CTRL] + [X]	opens and closes the stimuli selection
[CTRL] + [TAB]	steps forward through the data view tabs
[CTRL] + [SHIFT] + [TAB]	steps backwards through the data view tabs
[CTRL] + [MOUSEWHEEL]	only when zoom ^[69] is available: zooms in and out

[AOI Editor](#)^[75] keyboard commands

Keys	Description
[DEL]	deletes selected AOIs
[HOME]	jumps to first key frame
[END]	jumps to last key frame
[PG Up]	goes to next key frame
[PG Dn]	goes to previous key frame
[CTRL] + [Z]	undo action
[CTRL] + [Y]	redo action
[V]	toggles the visibility of the selected AOI
[D]	deletes current keyframe

[SHIFT] + [MOUSEWHEEL]	changes the size of a selected AOI
-------------------------------	------------------------------------

Video keyboard commands

The following keyboard commands are available to navigate in a video (see [Player Control](#)^[66]). They are available in the [AOI Editor](#)^[75], [Scan Path](#)^[108], [Attention Map](#)^[119] and [Key Performance Indicators](#)^[130] data views.


Keys	Description
[SPACE]	plays/pauses the presentation
Right arrow key	moves presentation one step forward according to the selected step size
Left arrow key	moves presentation one step backward according to the selected step size
Arrow up key	increases the step size
Arrow down key	decreases the step size
[CTRL] + [HOME]	jumps to the begin of the trial resp. the selected time window
[CTRL] + [END]	jumps to the end of the trial resp. the selected time window
[B]	set and resets a bookmark
[CTRL] + arrow left	jump to previous bookmark
[CTRL] + arrow right	jump to next bookmark
[ALT] + arrow right	Jumps to the next user event
[ALT] + arrow left	Jumps to the previous user event
[SHIFT] + arrow right	Jumps to the next annotation
[SHIFT] + arrow left	Jumps to the previous annotation
[CTRL] + [ENTER]	Add/Edit annotation

Line Graph  keyboard commands

Keys	Description
Left arrow key	moves selected time cursor to the left
Right arrow key	moves selected time cursor to the right

Appendix

Chapter



10 Appendix

10.1 About Box

To get general information about BeGaze 2.4 go to the Help menu of the [Menu Commands](#) [246] and select About BeGaze 2.4.




- **BeGaze 2.4 Version:** The line displays the current version number.
- **Copyright:** The line displays copyright information.
- **Home Page:** Here you can visit our home page.
- **Licensed data view packages:** BeGaze 2.4 is licensed to one computer only. Here you can see a list with all licensed data view packages.
- **Copy to Clipboard:** In a service case please click here to copy to clipboard detailed information about each licensed data view and

report this to the customer support and service team of your local distributor or [SMI](#) ^[275].

10.2 Dongle - Installation and Troubleshooting

BeGaze 2.4 is dongle-protected. You may have to place the USB-dongle in the appropriate PC before you can start the program. If BeGaze 2.4 displays a message box stating **HASP SRM Protection System: The software requires a hardware key (dongle)**, check the following:


1. The activity LED of the USB-dongle should show a red light if the dongle is plugged in.
2. If the activity LED does not show a red light, check the USB port status in the Windows hardware settings dialog. Open the Windows **Control Panel** and double click the **System** icon. Switch to the **Hardware** tab and click on the **Device Manager** button. Verify, that the **Universal Serial Bus controllers** tree does not show any yellow warning signs (⚠). The screen shot below shows a functional USB port with a correct Windows driver installation.



If the dialog displays a warning sign (⚠) for a driver, right click the entry and select the **Update Driver...** command from the context menu.

3. Verify, that the dongle driver is installed properly. Open the Windows **Control Panel** and double click the **Add or Remove Programs** icon.

Check if the list shows the **HASP SRM Run-time** entry.



Note, that the **HASP SRM Run-time** is installed during the installation of BeGaze 2.4. Do not deny the installation of this software during installation when prompted.



Type and status of your licenses are stored on the dongle device, not on the PC on which BeGaze 2.4 is installed. With the license update procedure, the dongle is updated. That means, that you can run BeGaze 2.4 on any PC when the dongle is plugged in.

10.3 Experiment Types

The eye tracking experiments fall into two major groups:

- experiments with eye tracking data (standard data)
- experiments with eye tracking and head tracking data

Dependent on the type of experiment the way data is collected differs slightly.

10.4 Database

All BeGaze 2.4 experiments will be collected in a database. Once you imported the data files, images and AOI files in BeGaze 2.4, you will no longer have to keep in mind the location of these files as they are stored bundled in the database.

The path where the database is located can be changed by going to the **File** menu and selecting **Change Data Storage Location**.

Initially, the database is located in the user's data folder. This corresponds to "Application Data" folder in Windows XP and "AppData \Roaming" folder in Windows Vista. For example, if your computer is running Windows XP and your user name is "BegazeUser", the complete path to the database will be: C:\Documents and Settings\BegazeUser\Application Data\SMI\BeGaze 2\BeGaze 2 Data Base.

If more users decide upon sharing the data base, they should change data storage location to a local folder where all have enough security rights.

Due to performance and concurrent access issues, a common network folder should not be used.



Note that the Change Data Storage Location menu command is available only if all experiments are closed.

10.5 System Requirements

Hardware requirements

BeGaze 2.4 should be installed on a personal computer or laptop with the following minimum requirements:

- OS: Windows XP Service Pack 2 / Windows Vista / Windows 7
- CPU: AMD or Intel Dual Core with 2.6 GHz

RAM: minimum 2 GB

VGA: 3D accelerated, 512 MB RAM, DirectX 9 Compatible, OpenGL V1.2 compatible

HDD: at least 10 GB of free hard disk space

For best views the monitor should be of size 19" or bigger with a minimum resolution of 1280x1024 pixels.

For database backups a DVD writer is recommended.

Some functions of BeGaze 2.4 need a printer connected.



Graphic card compatibility with OpenGL

BeGaze 2.4 is using OpenGL functionality in order to achieve best performance. The graphic card needs to be compliant with the OpenGL standard V1.2. Unfortunately not all graphic card drivers fully support this OpenGL standard, even though they are giving compliance statements to OpenGL. This might result in corrupted visualizations in the scan path and attention map views.

The OpenGL version can be verified with the Extension Viewer from RealTech VR:

<http://www.realtech-vr.com/glview/index.html>

Compliant and non-compliant graphic cards for Experiment Center and BeGaze

The following list contains the tested graphic card models that are compliant (recommended = yes) and non compliant (recommended=no) with Experiment Center and BeGaze.

(This list is not intended to be complete)

Recommended	Vendor	Model	Memory (MB)	Shared Memory	OpenGL Version
yes	NVIDIA	GeForce 7600 GS	256	No	2,1
yes	NVIDIA	GeForce 8500 GT	512	No	2,1
yes	NVIDIA	GeForce 9600 GT	512	No	3,0

yes	NVIDIA	GeForce 6200	128		2,1
yes	Intel	GMA 3100	384	Yes	1,4
yes	NVIDIA	GeForce 9800 GT	512	No	3,1
yes	Winfast	Geforce 8800 GTS	320	No	2,1
yes	ATI	Radeon X1050	256		2,1
yes	NVIDIA	GeForce 8600 GT	256	no	3,2
yes	NVIDIA	GeForce 9500 GT	512	no	3
yes	NVIDIA	GeForce 9400	512	no	3,2
yes	ATI	Mobility Radeon 9000 IGP	128		1,3
no	NVIDIA	GeForce 5200 FX	128	No	2,1
no	NVIDIA	GeForce 8800 GTS	320	No	2,1
no	ATI	FireGL V 3400			
no	NVIDIA	GeForce 8400			
no	NVIDIA	Quadro FX1700			
no	NVIDIA	Quadro FX570			
no	NVIDIA	Quadro FX5500			
no	Matrox	Orion	32 MB		
no	ATI	FireGL V 3100	128 MB		
no	Matrox	G550 DH	32 MB	no	

10.6 Limits

SMI guarantees BeGaze to work within the following limits:

Max. number of stimuli in one experiment	250
Max. number of trials per stimuli	250
Max. length of video / max. number of videos	2h / 5
Max. length of video / max. number of videos	1min / 200


Max. file size of video	1GB
Max. number of subjects per experiment	200
Max. length per trial / max. number of stimuli	2h / 5
Max. length per trial / max. number of stimuli	10min / 200
Max. number of AOs per stimulus	250
Max. stimulus size (excl. Web)	1680x1050
Max. stimulus size for Web	1680x10.00 0



Due to a limitation of Microsoft video handling, avi files are limited to 1GB of size. HED videos are split into multiple trials when the 1GB file size is exceeded.

10.7 Program Installation

The product installation media (CD-Rom) offers suitable software packages to install. Please run the auto-start application from the installation medium and click on the respective buttons to install necessary software.



The Experiment Suite 360° includes the BeGaze 2.4 as well as the Experiment Center 2.4 software. To install the Experiment Suite 360°, proceed as follows:

1. Insert the installation media (CD-Rom).
The auto-start application opens.
2. Click on the **Install from CD** button.
Follow the steps of the installation wizard.



While installing the Experiment Suite 360°, the USB dongle driver (HASP SRM Run-time) is installed or updated. You may need to

update the USB dongle license information. Refer to [Dongle Protection and License Update](#) ⁽¹²⁾ for details.

The Microsoft .NET Framework, the Microsoft Internet Explorer, and the Microsoft Media Player software components are available from the BeGaze 2.4 installation media. These software components are also available from the Microsoft web site where you can download them for installation to the desired PC workstation. Both software components will inspect your PC workstation during installation and may issue warning messages if the PC resources do not meet the necessary performance.



Please use always the latest versions that are available for download from the Microsoft web site.

Copyright and Trademarks

Chapter



XI

11 Copyright and Trademarks

Copyright

The SOFTWARE is owned by SensoMotoric Instruments GmbH or its suppliers and is protected by the Federal Republic of Germany copyright laws and international treaty provisions. Therefore, you must treat the SOFTWARE like any other copyrighted material (e.g. book or musical recording) except that you may either

- a) make one copy of the SOFTWARE solely for backup or archival purposes or
- b) transfer the software to a single hard disk provided you keep the original solely for backup or archival purposes.

You may not copy the written materials accompanying the SOFTWARE.

The user is not entitled to allow a third party to use the software simultaneously without written approval of SensoMotoric Instruments GmbH. Independent branch offices or subsidiary companies are also understood to be a third party in this sense. SensoMotoric Instruments GmbH and/or its supplying firm remain the owners of the delivered software, even if it is altered.

Trademarks

BeGaze is a trademark of SensoMotoric Instruments GmbH

iView X is a trademark of SensoMotoric Instruments GmbH

Experiment Center is a trademark of SensoMotoric Instruments GmbH

SensoMotoric Instruments is a trademark of SensoMotoric Instruments GmbH.

Microsoft, Windows are registered trademarks of Microsoft Corporation.

License Agreement and Warranty

Chapter



XII

12 License Agreement and Warranty

IMPORTANT – PLEASE READ CAREFULLY: This license agreement (“Agreement”) is an agreement between you (either an individual or a company, “Licensee”) and SensoMotoric Instruments GmbH (“SMI”). The “Licensed Materials” provided to Licensee subject to this Agreement include the software programs, which were installed on the equipment Licensee purchased from SMI (the “Designated Equipment”) at the time of delivery of the Designated Equipment and/or other software programs SMI has granted Licensee access to (the “Software”) as well as any “on-line” or electronic documentation associated with the Software, or any portion thereof (the “Documentation”), as well as any updates or upgrades to the Software and Documentation, if any, or any portion thereof, provided to Licensee at SMI’s sole discretion.

By installing, downloading, copying or otherwise using the Licensed Materials, you agree to abide by the following provisions. This Agreement is displayed for you to read prior to using the Licensed Materials.

If you do not agree with these provisions, do not download, install or use the Licensed Materials. If you have already paid for the Licensed Materials, you may return them for a full refund to: SensoMotoric Instruments, Warthestraße 21, 14513 Teltow, Germany. If the Licensed Materials were installed on Designated Equipment, you may also return such Designated Equipment for a full refund.

1. **License.** Subject to the terms of this Agreement, SMI hereby grants and Licensee accepts a non-transferable, non-exclusive, non-assignable license without the right to sublicense to use the Licensed Materials only (i) for Licensee’s business operations, (ii) on the Designated Equipment, if any, or on a back-up equipment while the Designated Equipment is inoperable, and (iii) in accordance with the Documentation. Licensee may make one copy of the Software in machine readable form for backup purposes only; every notice on the

original will be replicated on the copy. Installation of the Software, if any, is Licensee's sole responsibility.

2. **Rights in Licensed Materials.** Title to and ownership in the Licensed Materials and all proprietary rights with respect to the Licensed Materials and all copies and portions thereof, remain exclusively with SMI. The Agreement does not constitute a sale of the Licensed Materials or any portion or copy of it. Title to and ownership in Licensee's application software that makes calls to but does not contain all or any portion of the Software remains with Licensee, but such application software may not be licensed or otherwise transferred to third parties without SMI's prior written consent.
3. **Confidentiality.** Licensed Materials are proprietary to SMI and constitute SMI trade secrets. Licensee shall maintain Licensed Materials in confidence and prevent their disclosure using at least the same degree of care it uses for its own trade secrets, but in no event less than a reasonable degree of care. Licensee shall not disclose Licensed Materials or any part thereof to anyone for any purpose, other than to its employees and sub-contractors for the purpose of exercising the rights expressly granted under this Agreement, provided they have in writing agreed to confidentiality obligations at least equivalent to the obligations stated herein.
4. **No Reverse Engineering.** License shall not, and shall not allow any third party to, decompile, disassemble or otherwise reverse engineer or by any means whatsoever attempt to reconstruct or discover any source code or underlying ideas, algorithms, file formats or programming or interoperability interfaces of the Software or of any files contained or generated using the Software. The obligations under this Section shall survive any termination of the Agreement. Licensee may apply to SMI at the address above for further information on the interface between the Software and other applications; SMI will on reasonable terms and conditions provide such information as the Licensee may reasonably require to enable the Software to interoperate with its applications.
5. **Warranty.**
 - a) If Licensed Materials show a defect during the limitation period

- specified in Sec. 5 d), which already existed at the time of transfer of risk, SMI will provide Licensee with a corrected version of such defective Licensed Material (“Repeat Performance”).
- b) Licensed Materials are defective, if they do not conform to the written agreement between SMI and Licensee, if any, or to SMI’s specification.
 - c) Licensee’s claims for defects are excluded:
 - for insignificant divergences from the agreement or specification;
 - for damages caused after the transfer of risk by external influences such as fire, water, currency surge, etc.; improper installation, operation, use or maintenance; use in fields of application and environmental conditions other than those expressly specified by SMI; and use in combination with other products not approved by SMI for this purpose, excessive stress or normal wear and tear;
 - for defects of the Software, which cannot be reproduced;
 - for Licensee’s or a third party’s modifications of Licensed Materials and results therefrom;
 - in as far as Licensee has not notified SMI of apparent defects within 5 business days of delivery and of hidden defects within 5 business days of discovery.
 - d) The limitation period for warranty claims is 12 months. This does not apply in cases of fraud, intent, gross negligence and personal injury. Repeat performance does not restart the limitation period.
 - e) Licensee has to always first provide SMI the opportunity for Repeat Performance within a reasonable period of time. Repeat Performance does not constitute acceptance of a legal obligation.
 - f) If the Repeat Performance fails even within an additional period of time of reasonable length to be granted by Licensee, Licensee shall, notwithstanding its claims under Section 6, if any, be entitled to rescind the contract or to claim a price reduction.
 - g) If the analysis of an alleged defect shows that it is not covered by the above warranty, SMI may charge for the failure analysis and correction of the defect, if any, at SMI’s then applicable rates.
6. **Liability Limitations.**
- a) SMI is only liable for damages caused by slight negligence if such are due to a material breach of duty, which endangers the achievement of the objective of the Agreement, or the failure to

comply with duties, the very discharge of which is an essential prerequisite for the proper performance of the Agreement.

- b) In cases of Section 6 a), the liability is limited to the damage, which is typical for contracts like this Agreement and which could have been foreseen.
 - c) SMI's liability is also limited to the damage, which is typical for contracts like this Agreement and which could have been foreseen for damages caused by the gross negligence of an agent or an employee of SMI, who is not an officer or executive of SMI.
 - d) In cases of Sections 6 a) and c), SMI's liability is limited to a maximum amount of EURO 500.000, respectively EURO 100.000 for financial losses.
 - e) Licensee's claims for damages caused by slight negligence or caused by the gross negligence of an agent or an employee of SMI, who is not an officer or executive of SMI, which are not based on defects and thus not subject to time-barring in accordance with Sec. 5 d) are time-barred at the latest 2 years from the point in time Licensee became aware of the damage and regardless of Licensee's awareness at the latest 3 years after the damaging event.
 - f) With the exception of liability under the Product Liability Law, for defects after having given a guarantee, for fraudulently concealed defects and for personal injury, the above limitations of liability shall apply to all claims, irrespective of their legal basis, in particular to all claims based on breach of contract or tort.
 - g) The above limitations of liability also apply in case of Licensee's claims for damages against SMI's employees or agents.
7. **SMI Indemnity.** SMI will defend or settle any action brought against Licensee to the extent based on a claim that Licensed Materials, used within the scope of the license granted in this Agreement, infringe any copyright valid in the European Union or the European Economic Area and will pay the cost of any final settlement or judgment attributable to such claim, provided (i) Licensee has given prompt notice to SMI of such claim, (ii) Licensee has not recognized an infringement of the third party's copyright, and (iii) the entire defense and the settlement negotiations are reserved for SMI with Licensee's full cooperation and support. If Licensee discontinues the use of Licensed Materials for damage reduction or other important reasons, Licensee is obligated to point out to the third party that the

discontinuation of use does not constitute an acknowledgement of a copyright infringement. If SMI believes Licensed Materials are likely to be the subject of an infringement claim, it may elect to replace or modify such Software or Documentation to make it non-infringing or terminate the Agreement on written notice to Licensee. SMI shall have no obligation to defend (or any other liability) to the extent any claim involves a Software release other than the current, unaltered release, if such would have avoided infringement, or use of the Software in combination with non-SMI programs or data, unless the infringement would also incur without such combination. The foregoing states the entire obligation and liability of SMI with respect to any infringement by Licensed Materials of any intellectual property rights or other proprietary rights of Licensee or a third party without prejudice to any claims for damages in accordance with Section 6.

8. **Licensee Indemnity.** Licensee will defend and indemnify SMI, and hold it harmless from all costs, including attorney's fees, arising from any claim that may be made against SMI by any third party as a result of Licensee's use of Licensed Materials, excluding claims for which SMI is obligated to defend or indemnify Licensee under Section 7.
9. **Export Restriction.** Licensee will not remove or export from Germany or from the country Licensed Materials were originally shipped to by SMI or re-export from anywhere any part of the Licensed Materials or any direct product of the Software except in compliance with all applicable export laws and regulations, including without limitation, those of the U.S. Department of Commerce.
10. **Non-Waiver; Severability; Non-Assignment.** The delay or failure of either party to exercise any right provided in this Agreement shall not be deemed a waiver. If any provision of this Agreement is held invalid, all others shall remain in force. Licensee may not, in whole or in part, assign or otherwise transfer this Agreement or any of its rights or obligations hereunder.
11. **Termination.** This Agreement may be terminated without any fee reduction (i) by Licensee without cause on 30 days notice; (ii) by SMI, in addition to other remedies, if Licensee fails to cure any breach of its obligations hereunder within 30 days of notice thereof; (iii) on notice by either party if the other party ceases to do business in the

normal course, becomes insolvent, or becomes subject to any bankruptcy, insolvency, or equivalent proceedings. Upon termination by either party for any reason, Licensee shall at SMI's instructions immediately destroy or return the Licensed Materials and all copies thereof to SMI and delete the Software and all copies thereof from the Designated Equipment.

12. **Entire Agreement; Written Form Requirement.** There are no separate oral agreements; any supplementary agreements or modifications hereto must be made in writing. This also applies to any waiver of this requirement of written form.
13. **Notices.** All notices under the Agreement must be in writing and shall be delivered by hand or by overnight courier to the addresses of the parties set forth above.
14. **Applicable Law and Jurisdiction.** German law applies with the exception of its conflict of laws rules. The application of the United Nations Convention on Contracts for the International Sale of Goods (CISG) is expressly excluded. The courts of Berlin, Germany, shall have exclusive jurisdiction for any action brought under or in connection with this Agreement.

Teltow, Germany, 2004-2010
SensoMotoric Instruments GmbH

About SMI

Chapter



XIII

13 About SMI

SensoMotoric Instruments (SMI) is a world leader in dedicated computer vision applications, developing and marketing eye & gaze tracking systems and OEM solutions for a wide range of applications.

Founded in 1991 as a spin-off from academic research, SMI was the first company to offer a commercial, vision-based 3D eye tracking solution. We now have over 17 years of experience in developing application-specific solutions in close collaboration with our clients.

We serve our customers around the globe from our offices in Teltow, near Berlin, Germany and Boston, USA, backed by a network of trusted local partners in many countries.

Our products combine a maximum of performance and usability with the highest possible quality, resulting in high-value solutions for our customers. Our major fields of expertise are:

- Eye & gaze tracking systems in research and industry
- High speed image processing, and
- Eye tracking and registration solutions in ophthalmology.

More than 4,000 of our systems installed worldwide are testimony to our continuing success in providing innovative products and outstanding services to the market. While SMI has won several awards, the largest reward for us each year is our trusted business relationships with academia and industry.

Please contact us:

Europe, Asia, Africa, South America, Australia

SensoMotoric Instruments GmbH (SMI)
Warthestraße 21
D-14513 Teltow
Germany
Phone:+49 3328 3955 0
Fax:+49 3328 3955 99
email: info@smi.de

North American Headquarters

SensoMotoric Instruments, Inc.
75 Arlington Street
Boston, MA 02116
USA
Phone: +1 (857) 241 3865
Fax: +1 (857) 241 3601
Toll-Free: 888 SMI USA1
email: info@smiusa.com

Please also visit our home page: <http://www.smivision.com>

Copyright © 2010 SensoMotoric Instruments GmbH

Last updated: February 2010

Index

- A -

About

- BeGaze 2.4 256
- Manual 4
- SMI 275

Adjust Event Detection 216

Analysis Menu 246

AOI

- Change Position 80
- Change Shape/Size 80
- Create 78
- Delete 93
- Edit 80
- Edit Properties 86
- Format Description 94
- Load 94
- Priority 80
- Rename 80
- Save 94
- Visibility 90

AOI Editor

- Key Frames 92
- Open 78
- Overview 75
- Toolbar 77

AOI Sequence Chart

- Main Pane 147
- Plug-in 146

Attention Map 119, 125

- Video Export 239

AVI, Optimizing 241

- B -

Backup Experiment 43

Bee Swarm

- Appearance 106
- Main Pane 104
- Plug-in 102
- Settings 105

Bee Swarm Settings

- Bee Swarm Tab 106
- Cursor Tab 107

BeGaze 2.4 Version 256

Binning Chart

- Main Pane 151
- Plug-in 150

Built-In Detector 215

- C -

Chart Display Modes 74

Control Pane 66

Converting Video 242

Copyright 266

Cursor Tab, Bee Swarm 107

Cursor Tab, Scan Path 116

Cursors 208

Custom Map 121, 126

- D -

Data Import 30

Details

- Events 63
- Subject 60

Diagrams 51

Dongle

- Dongle
 - General Information 12
 - Installation 257
 - Troubleshooting 257
- E -**
- Event
 - Details 63
 - Detection, Experiment Wizard 37
- Event Statistics
 - Overview 152
 - Selection Trees 153
 - Templates 155
 - Time Interval 157
- Events
 - Export 225
 - Selection 61
- Experiment
 - Backup 43
 - Close 43
 - Delete 46
 - Load from Folder 28
 - Menu 246
 - Modify 41
 - Open 42
 - Restore 45
 - Save 42
- Experiment Wizard
 - Data Import 30
 - Event Detection 37
 - General Settings 29
 - Measurement Scenario 39
 - Signal 39
 - Stimuli 34
 - Stimulus Association 36
- Export 236
- Events 225
- File Format 228
- Menu 246
- Overview 225
- Raw Data 232
- Raw File Format 236
- Video 239
- F -**
- File Export 225, 228, 229, 236, 237
- File Format
 - Export 228
 - Header 229
 - Raw File Format 236
 - Trial Section 229
- File Menu 246
- Filters, Subjects 56
- Fixation
 - Detection 215, 219
 - Parameters 216
- Fixations, Scan Path 117
- Focus Map 121, 126
 - Main Pane 121
 - Plug-in 119
 - Settings 122
- G -**
- Gaze Replay
 - Plug-in 100
- General Settings, Experiment Wizard 29
- Getting Started 21

Gridded AOIs

 Main Pane 139

 Plug-in 137

Group Properties 56

- H -

Header, File Format 229

Heat Map 121, 122, 126

 Main Pane 126

 Plug-in 125

 Settings 128

Help Menu 246

High Speed Event Detection 221

Hotkeys 251

How BeGaze works 15

- I -

Installation 263

Introduction 2

- K -

Key Features 8

Key Frames 92

Key Performance Indicators

 Main Pane 131

 Plug-in 130

 Settings 133

 Video Export 239

- L -

License

 Agreement 268

 Update 12

Line Graph

Data 208

Diagram Cursors 208

Events 204

Main Pane 206

Miniview 209

Plug-in 203

Settings 210

Time 208

Low Speed Event Detection 219

- M -

Main Pane

 AOI Sequence Chart 147

 Binning Chart 151

 Focus Map 121

 Gridded AOIs 139

 Heat Map 126

 Key Performance Indicators 131

 Line Graph 206

 Scan Path 104, 110

Manual 4

Measurement Scenario 39

Menu Commands 246

Modify, Experiment 41

- O -

Online Help 4

Optimizing, AVI 241

- P -

Peak Velocity Window 216

Playback Control 67

Player Control 66

Plug-in

- Plug-in
 - AOI Sequence Chart 146
 - Bee Swarm 102
 - Binning Chart 150
 - Event Statistics 152
 - Focus Map 119
 - Gaze Replay 100
 - Gridded AOIs 137
 - Heat Map 125
 - Key Performance Indicators 130
 - Line Graph 203
 - Operating Panes 54
 - Overview 50
 - Panes Overview 53
 - Scan Path 108
 - Selection 51
- Primary Event
 - Fixation 219
 - Saccade 221
- Print Menu 246
- Product Variants 11
- Program
 - Installation 263
- R -**
- Raw Data Export 232
- Raw File Format 237
- Requirements 259
- Restore, Experiment 45
- S -**
- Saccade
 - Detection 215, 221
 - Parameters 216
- Scan Path
 - Appearance 113
 - Main Pane 110
 - Plug-in 108
 - Settings 113
 - Video Export 239
- Scan Path Settings
 - Cursor Tab 116
 - Fixations Tab 117
 - Scan Path Tab 113
- Selection
 - Events 61
 - Plug-in 51
 - Subjects 56
- Selection Trees, Statistics 153
- Settings
 - Bee Swarm 105
 - Focus Map 122
 - Heat Map 128
 - Key Performance Indicators 133
 - Scan Path 113
- Shortcuts 251
- Signal 39
- Stimulus 34
 - Association 36
 - Selection 55
- Subjects
 - Details 60
 - Filters 56
 - Selection 56
- System Requirements 259
- T -**
- Templates, Event Statistics 155

Thumbnail Control
 Context Menu 72
 Overview 70

Time Interval, Event Statistics 157

Toolbar
 AOI Editor 77
 Common 249

Trademarks 266

Trial Section
 File Format 229
 Raw File Format 237

- U -

Use Cases 18

- V -

Video
 Background Information 244
 Converting 242
 Export 239
 Optimizing 241

View Menu 246

VirtualDub 242

Visibility, AOI 90

- W -

Wizard 27

- Z -

Zoom Control 69

